
# **VivoHeart Intravascular Ultrasound Imaging System**

## **INSTRUCTION FOR USE (Indonesian) INSTRUKSI PENGGUNAAN**

*Insight Lifetech*

## Introduction of Instruction for Use

This manual includes instructions for the operator, including system introduction, system usage, exclusion of basic malfunction, and basic maintenance.

This system needs to be simultaneously connected to other devices during operation. Please refer to the instruction manuals of these devices for detailed information on operation and safety.

**Instructions of warning, caution, and attention are as follows:**

**Warning:** This may involve the personal safety of the patient or user. Ignoring this information may result in impairment to the personal safety of the patient or user.

**Caution:** To avoid damage to the equipment, you must follow special services or precautions indicated in these instructions.

**Attention:** This document provides information that may be particularly useful to you or help maintain or clarify important instructions.

Product is under review by Ministry of Health Indonesia

## Manufacturer Information

Name of Manufacturer	Insight Lifetech Co., Ltd.	
Address of Registrant/Manufacturer	2102, Tower A, Hongrongyuan North Station Center, North Station Community, Minzhi, Longhua District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA	
Manufacturing Site	Zone A, F, 2 <sup>nd</sup> Floor, Zone C, E, F, 3 <sup>rd</sup> Floor, Zone A, 4 <sup>th</sup> Floor, and Zone A, B, 10 <sup>th</sup> Floor, Building 3, Tinwe Business Park, No.6 Liufang Road., Baoan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA	
Address of after-sales service provider	Insight Lifetech Co., Ltd.	
Address of after-sales service provider	Zone E, 3 <sup>rd</sup> Floor, Building 3, Tinwe Business Park, No. 6 Liufang Rd., Baoan District, 518000, Shenzhen, PEOPLE'S REPUBLIC OF CHINA	
Number of Medical Device Registration Certificates	国械注准 20223060975	
Number of Product Technical Requirements	国械注准 20223060975	
Number of Manufacturing Licenses	粤食药监械生产许 20204044 号	
IFU prepared on	2023.06.26	
Version of the IFU	01	
Software	VivoHeart software (console)	VivoHeart software (pullback unit)
Software Version	V1.1	V1.1
Product Service Life	8 years	
Service Life of the Pullback Unit	2 years	

## Contents

1. Product Overview	1
1.1. Operating principles	1
1.2. Product description	1
1.3. Intended purpose	1
1.4. Contraindications	1
1.5. Adverse events	2
1.6. Intended users	2
1.7. Model	2
1.8. General warnings, cautions, precautions, and safety information	3
1.9. Symbols with description	7
2. Structural characteristics of the IVUS System	9
2.1. Appearance	9
2.2. Cart system setting	9
2.2.1 Displayer	11
2.2.2 Casters	11
2.3. Pullback unit	12
3. Preparation of the IVUS System	13
3.1. Movement of the IVUS System	13
3.1.1 Device transportation	13
3.1.2 Device movement	13
3.2. Installation of the IVUS System	13
3.1.3 Initial inspection	13
3.1.4 Connection of alternating current power cord	14
3.3. Power on and off	14
4. Use of the IVUS System	16
4.1. Overview	16
4.2. Preparation for usage of the pullback unit and IVUS Catheter	16
4.2.1 Preparation of the pullback unit	16
4.2.2 Connect the catheter to the pullback unit	17
4.2.3 Catheter replacement	18
4.2.4 Connect the pullback unit cable to the cart	18
4.3. Starting the measurement	19
4.3.1 Logging into the IVUS System	19
4.3.2 Creating patient data	21
4.3.3 Real-time imaging	22
4.3.4 Viewing data	27
4.3.5 Image measurement	34
4.3.6 Viewing archive data	40
4.3.7 Shutting down the console system	41
5. Settings	42
5.1. General settings	42
5.2. Account (user) settings	43

5.3.	Doctor settings . . . . .	44
5.4.	Comment settings. . . . .	45
5.5.	Measurement settings. . . . .	46
5.6.	DICOM settings. . . . .	46
5.7.	Log settings . . . . .	46
5.8.	Factory settings . . . . .	47
6.	Technical Specifications . . . . .	47
6.1.	Performance requirements . . . . .	47
6.2.	Device classification and regulatory compliance . . . . .	47
6.3.	Environmental conditions . . . . .	48
6.4.	Port Specifications . . . . .	48
6.5.	RFID emitter . . . . .	49
6.6.	Cable information . . . . .	49
7.	Basic troubleshooting . . . . .	50
8.	Cleaning and Maintenance . . . . .	57
8.1.	Cleaning . . . . .	57
8.2.	Maintenance. . . . .	57
8.2.1	Inspection. . . . .	57
8.2.2	Maintenance plan. . . . .	58
8.2.3	Fuse replacement . . . . .	58
8.2.4	Printer maintenance . . . . .	58
9.	Warranty and After-sales Services . . . . .	62
9.1.	Warranty. . . . .	62
9.2.	After-sales services . . . . .	62
10.	Acoustic Output Information . . . . .	63
11.	Electromagnetic Compatibility Information. . . . .	66

## List of Figures

Figure 1 Front view of cart . . . . .	9
Figure 2 Rear view of cart . . . . .	10
Figure 3 Side view of cart. . . . .	10
Figure 4. Front view of displayer . . . . .	11
Figure 5. Rear view of displayer. . . . .	11
Figure 6. Casters . . . . .	11
Figure 7. Structure of the pullback unit . . . . .	12
Figure 8. Illustration of inserting the pullback unit into its sterile bag. . . . .	16
Figure 9. Illustration of catheter connection to the pullback unit. . . . .	17
Figure 10. Illustration of successful catheter connection . . . . .	17
Figure 11. Connecting the pullback unit cable and the cart . . . . .	18
Figure 12. Illustration of successful connection . . . . .	19
Figure 13. Login interface. . . . .	20
Figure 14-1. Real-time interface. . . . .	20
Figure 14-2. Real-time interface. . . . .	21
Figure 15. Patient interface. . . . .	22
Figure 16. Imaging display interface of hidden longitudinal section view. . . . .	23
Figure 17. Imaging display interface with longitudinal section view enabled. . . . .	24
Figure 18. Gain adjustment interface . . . . .	26
Figure 19. Graymap adjustment interface. . . . .	26
Figure 20. Viewing the exam . . . . .	27
Figure 21. Viewing the list interface. . . . .	28
Figure 22. Editing exam information interface. . . . .	29
Figure 23. Image annotation interface . . . . .	30
Figure 24. Bookmark interface . . . . .	31
Figure 25. Image comparison interface . . . . .	33
Figure 26. Length measurement interface . . . . .	35
Figure 27. Longitudinal section length measurement interface . . . . .	36
Figure 28. Area measurement interface (area of a surface plotted by cursor trajectory) . . . . .	38
Figure 29. Area measurement interface (stipple area) . . . . .	38
Figure 30. Archive interface . . . . .	40
Figure 31. General settings interface . . . . .	42
Figure 32. Account settings interface . . . . .	43
Figure 33. Doctor settings interface . . . . .	44
Figure 34. Comment settings interface. . . . .	45
Figure 35. Measurement settings interface. . . . .	46
Figure 36. DICOM settings interface . . . . .	46

## **1. Product Overview**

### **1.1. Operating principles**

VivoHeart Intravascular Ultrasound Imaging System (hereinafter called IVUS System) is used in conjunction with the TrueVision Intravascular Ultrasound Imaging Catheter (hereinafter called IVUS Catheter) for imaging of vascular structures. The IVUS System emits electric pulses to stimulate the ultrasonic transducer at the front end of the IVUS Catheter to generate ultrasonic waves. The ultrasonic waves propagate in the ultrasonic medium in beams and will be scattered and reflected when they encounter the interface of two different media with different acoustic impedance. Then the reflected ultrasonic waves will be received by the ultrasonic transducer to convert into electrical signals, which will be delivered to the IVUS System to transform into graymap images. Different types of tissues respond differently in ultrasound absorption and reflection, and the IVUS System displays different gray scales based on the intensity of the received ultrasonic signals, to form graymap images reflecting the vascular structure.

### **1.2. Product description**

The VivoHeart Intravascular Ultrasound Imaging System consists of the imaging console (including the cart, displayer, keyboard, and mouse), and the pullback unit.

### **1.3. Intended purpose**

The VivoHeart Intravascular Ultrasound Imaging System is intended to be used in medical institutions in conjunction with the TrueVision Intravascular Ultrasound Imaging Catheters (model: 02-6F40-01 and 02-6F60-01) for ultrasonic examination of intravascular lesions in coronary arteries. Intravascular ultrasound imaging is indicated in patients who are candidates for transluminal coronary interventional procedures.

### **1.4. Contraindications**

This IVUS Imaging System must be used in conjunction with the IVUS Catheter. Therefore, the IVUS Imaging System should not be used if the use of any catheter constitutes a safety risk to the patient. The contraindications generally include the following symptoms:

- Bacteremia or sepsis
- Severe coagulation disorder
- Contraindications for coronary intervention surgery
- Severe hemodynamic instability or shock
- Patients diagnosed with severe coronary artery spasms
- Life-threatening arrhythmia

### 1.5. Adverse events

Vascular imaging has all the risks and discomfort related to catheterization. These risks or discomfort may occur at any time, with various frequencies and severity. In addition, these complications may result in the need for additional treatments such as surgical intervention and may lead to death in rare cases, possible adverse reactions are as follows:

- Anaphylaxis
- Angina
- Cardiac arrest
- Arrhythmia (including but not limited to ventricular tachycardia, atrial/ventricular fibrillation, and complete atrioventricular block)
- Cardiac tamponade or pericardial effusion
- Death
- Device jam (requiring surgical intervention)
- Embolism (air, foreign body, tissue, or thrombus)
- Bleeding or hematoma
- Hypotension
- Infection
- Myocardial infarction
- Myocardial ischemia
- Stroke and transient cerebral ischemia
- Thrombosis
- Vascular occlusion and sudden closure
- Vascular injury (including but not limited to dissection and perforation)

### 1.6. Intended users

**General operator:** The VivoHeart Intravascular Ultrasound Imaging System should be used by a physician trained for the SOP of the catheterization lab or under his/her guidance.

**Maintenance:** Sales representatives or technicians designated by Insight Lifetech Co., Ltd. (Shenzhen) (hereinafter referred to as Insight Lifetech).

**Installation:** To be carried out by hospital technicians or IT personnel, sales representatives, or technicians designated by Insight Lifetech.

### 1.7. Model

VH-60

## 1.8. General warnings, cautions, precautions, and safety information

### Warning

1. The IVUS System should not be used at home.
2. Staff using the IVUS System must pass professional training before they are allowed to use the device.
3. Read this manual and the attached instruction manual for the IVUS Catheter before using the IVUS System.
4. Please do not modify the IVUS System.
5. Please insert the plug for the IVUS System into a grounded electrical socket. Otherwise, there will be a risk of electrocution.
6. Except for the alternating current fuse, the IVUS System does not contain any parts that require servicing by the operator. Please do not open the protective lid of the IVUS System to avoid electrocution. If the device requires servicing, please seek help from the staff authorized by Insight Lifetech. Otherwise, the service agreement will be invalid.
7. Please only use the fuses listed in this manual to replace the fuse to prevent fires from occurring. Please do not use a fuse with a higher current range than that stipulated in this manual.
8. Explosions may occur if flammable anesthetics are present in the environment.
9. The pullback unit of the IVUS System is connected using the floating double insulation (type of CF) which limits the current leakage to circuits in the stipulated range of GB 9706.1-2007, which prevents the defibrillator from affecting the device during electrical discharge.
10. Please use an IVUS Catheter that is compatible with the IVUS System. If the connected catheter cannot be recognized by the device, please do not continue using that catheter.
11. Please ensure that the pullback unit is at a stable level before and during automatic withdrawal to prevent it from turning over. Otherwise, it will threaten the safety of the patients.
12. JPEG compression will decrease the spatial effects and graymap quality of stored images, even when the highest (100%) settings are used. This may affect the diagnostic results of the images. Please do not use this system setting without approval from your institution.
13. Please do not touch the blade at the front of the printer as it may result in physical injury.
14. Please do not incinerate discarded printing paper.
15. The IVUS System should form an equipotential body when connected with other devices (effectively connected by equipotential).
16. If there is a need for inspection or change of the fuse at the alternating current input end, the power supply should be disconnected from the alternating current mains to prevent electrocution.
17. Please adjust all settings carefully to avoid covering signals with weak intensity but having diagnostic value. Improper settings will severely affect image quality.
18. If the alternating current power supply of the medical institution is unstable (power dip or surge), unforeseen situations may occur in the device which may affect image reliability.
19. Please do not place the pullback unit in an autoclave or any solution, or carry out disinfection for it. Please store the pullback unit in a sterile bag. Meanwhile, as the device is extremely sensitive, please gently take and release it to avoid collisions or falling. Please do not use the

pullback unit that is possibly damaged. Except for the compatible IVUS Catheter connector, do not insert any other objects in the pullback unit.

20. Please use a medical display that conforms to the requirements of GB 9706.1-2007 and YY0505-2012 when connecting to an external display. Any connected non-medical devices should meet the corresponding product safety requirements. Devices that only rely on basic insulation should not be used in the system. Please contact Insight Lifetech or the local agent if you have any questions.
21. Please do not draw self-intersecting geometric shapes when measuring the area. Otherwise, the geometric area reported by the system may be inaccurate.
22. Please do not use the IVUS System if there is liquid infiltration into the outer casing or power supply unit as this may cause electrocution to the user or patient. Please contact Insight Lifetech or the local agent if this situation arises.
23. Please do not sterilize the IVUS System or any components. Please do not use the IVUS Imaging System and components if they have undergone sterilization.
24. After the device and reusable accessories have passed their service life, they should be disposed of according to local laws or returned to the manufacturer for recycling. Improper waste disposal may lead to environmental pollution. Please contact the local municipal government to obtain more information on waste recycling.
25. Before using the device, the device and the cables must be inspected for any damage that may affect the safety of the patient. If obvious damage or aging is found, the component must be replaced before use to ensure that the device can be used normally and safely.
26. An electromagnetic field may interfere with the normal performance of the device. Therefore, please ensure that all other devices near the VivoHeart Intravascular Ultrasound Imaging System comply with relevant electromagnetic compatibility requirements.

### **Caution**

1. If the seal of the hemostatic valve is too tight, the signal transmission component of the IVUS Catheter will be compressed, which may result in device damage or image artifacts.
2. Please be careful when inserting the IVUS Catheter back into the opening of the sterile bag of the pullback unit to ensure that the sterile bag is not caught by the IVUS Catheter and the pullback unit.
3. Please do not connect or disconnect the IVUS Catheter when the pullback unit is still operating. Otherwise, the IVUS Catheter or pullback unit will be damaged.
4. If an unapproved catheter is connected or the IVUS Catheter is not connected correctly, the identification data and depth of the corresponding IVUS Catheter will not be displayed. At this point, the device will disable the imaging function. Please continue operation after solving the problem.
5. Please do not allow any object to block the linear motion path of the pullback plate, otherwise, the pullback operation of the device will be affected.
6. Except for a dedicated circular label, please do not paste any adhesive label on the CD. These

- labels will detach or unbalance the CD, thereby damaging the optical disc drive or optical disc.
7. Please store the recordable DVD media only in the DVD box. An ordinary CD box may compress the DVD too tightly. When the DVD is removed from the CD box, the DVD may be excessively bent, thereby resulting in data loss.
  8. Please do not paste any adhesive label on the recordable DVD as it will damage the DVD. Adhesive labels for CDs cannot be used on DVDs.
  9. Written labels on the CD or DVD can only be made using a pen with ribbon-made point or quick-drying permanent marker. Please do not use a slow-drying marker as the ink may damage the optical disc, thereby resulting in data loss. Never use a ballpoint pen or other pens with sharp tips as it may scratch the protective layer of the optical disc and damage the optical disc, thereby resulting in data loss.
  10. Please do not turn off the printer or system power supply when the printer is still running as it may damage the print head.
  11. It is recommended that the cart be pushed rather than pulled. When moving the cart, the tilt angle of the cart monitor can be rotated to avoid blocking the visual field. Please do not move the cart over objects to prevent the cart from tipping over or possible damage to the device.
  12. The display screen of the monitor contains an anti-glare coating that may be scratched by sharp or coarse objects. Please avoid touching the surface of the monitor. Please do not spray any substance directly on the monitor surface or outer casing.
  13. Please handle the printing paper with care to avoid dust and fingerprint oil stains accumulating on the paper.
  14. Please use the cleaning paper to clean the print head only when spots or streaks are present in the output image. Over-frequent cleaning will damage the print head.
  15. Please avoid stepping or pressing the power cable for the pullback unit during use.
  16. Never turn off the power supply to the monitor when the IVUS System is connected to the power supply. Otherwise, the resolution of the monitor may be affected. Please contact an Insight Lifetech sales representative when this problem occurs.
  17. Never press the power supply switch for the image processor when the device is in use. The power supply is controlled by software and the main alternating current supply switch. If the power supply switch of the image processor is pressed during exam recording, recorded data will be damaged and unusable.
  18. When turning off the power supply of the device, turn off the console first before turning off the power supply switch of the main alternating current. Turning off the power supply switch first will result in system malfunction.
  19. Before turning off the switch for the main alternating current power supply, shutdown the

application software of the device first. After the main alternating current power supply switch is at the OFF position, please wait for at least 10 seconds before turning the switch to the ON position.

20. The sterile bag of the pullback unit should not obstruct the insertion of the catheter head into the pullback unit.
21. If it is necessary to install the IVUS System near any other equipment, observe and ensure that it can work normally. If the IVUS System is found to be affected, please increase its distance from other equipment.

#### **Attention**


1. Ensure that the IVUS Catheter moves following a straight line. Excessive bending of the proximal end of the catheter will affect pullback measurement accuracy, thereby decreasing imaging quality.
2. Please do not touch or damage the non-labeled surface of the DVD or CD. Please comply with the processing, labeling, and storage instructions provided by the optical disc manufacturer.
3. Please check carefully when removing the label to ensure that the correct label is removed. There will not be any warnings on the system and this operation cannot be reversed.
4. The alternating current power supply should be a typical power supply that is suitable for a hospital environment. When the device is not working normally, you may have to turn off the system and restart it so that it resumes normal operations. Please comply with the normal shutdown procedure in the device chapter.
5. Special EMC preventive measures are required for medical devices. The device should be installed and used according to the EMC information in the attached system document. Portable and mobile RF communications may affect the performance of the electrical medical device.
6. Please label, process, and store the corresponding storage media according to the CD, DVD, or USB mobile storage device manufacturer's instructions to prevent data loss or damage over time. Otherwise, this may result in complete loss or damage to the saved image data.
7. Exam deletion cannot be reversed. Please confirm before deletion.
8. If the following problems occur when attempting to write DICOM files into portable media:
  - a) Medium type error
  - b) No medium is present in the driver
  - c) Insufficient storage space in the inserted CD, DVD, or USB mobile storage device

You must solve the corresponding problem if a text prompt appears on the device. If the problem is unsolved, the text prompts will be continuously shown on the screen. In this case, please determine the type of the problem from the aforementioned 3 problems and carry out corresponding measures (insertion of the correct medium, inserting the medium in the driver,

and ensuring that a blank medium is inserted). After the medium problem is solved, the device will resume normal operation.

9. Please move the cart by using the handle only.

### 1.9. Symbols with description

Symbol	Description
	Caution, consult accompanying documents
	Follow Instructions for Use
	Serial number
	Date of manufacture
	Defibrillation-proof type CF applied part
	Alternating current
	Equipotential grounding
	Separate collection for waste of electrical and electronic equipment
	On
	Off
	Non-ionizing radiation (RFID radiofrequency emitter)
	Temperature limit
	Humidity limitation
	Atmospheric pressure limitation
	This side up

Symbol	Description
	Fragile
	Keep dry
	Keep away from sunlight
	Do not roll
	Do not stack
	Warning hands pinching

## 2. Structural characteristics of the IVUS System

### 2.1. Appearance

- External dimensions of the VivoHeart Intravascular Ultrasound Imaging System [length × width × height]: (638±10) \*(551±10)\*(1515±50) mm
- Device weight: 80±5 kg

### 2.2. Cart system setting

The IVUS System provides broad and important functions to assist the IVUS Catheter in interventional intravascular ultrasonic imaging.


Figure 1 Front view of cart


Figure 2 Rear view of cart


Figure 3 Side view of cart

**Note:**

Please move the cart by using the handle only.

Please use a medical display that conforms to the requirements of GB 9706.1-2007 and YY0505-2012 when connecting to an external display using the DVI port. It is recommended that the Beacon medical display be used, with model No. of CL1930\*JABXUPO01.

### 2.2.1 Displayer

The displayer is a component of the IVUS System and there are 5 buttons labels in the lower right corner of the displayer, locating at the bottom of the lower right corner, with 1 indicator light, and 3 ports on the back of the displayer, as shown:


Figure 4. Front view of displayer


Figure 5. Rear view of displayer

### 2.2.2 Casters

The IVUS System is equipped with 4 casters and all casters contain brakes:

- To lock the casters, step on the levers (ON) on the 4 casters to immobilize the system.
- To unlock the caster, step on the levers (OFF) on the 4 casters to allow the system to move.


Figure 6. Casters

### 2.3. Pullback unit

The pullback unit is a component of the IVUS System and is the main control device for intravascular positioning and movements of the IVUS catheter. After catheter positioning is completed based on the instructions of the catheter, initial imaging can be performed to explore the target region, adjust image quality, and examine the morphology carefully. The following figure shows the structure of the pullback unit.


Figure 7. Structure of the pullback unit

The electronic and mechanical components of the pullback unit include:

Region No.	Description	Function
1	Pullback distance display	This component displays the movement distance/speed of the catheter ultrasonic transducer during the exam (precision: 0.1 mm, maximum value: 150 mm).
2	Pullback control button	Pressing this button causes the pullback unit to be manually moved along the track in the pullback plate.
3	Catheter release button	Pressing this button enables the catheter to be removed from the catheter interface.
4	Catheter port	Connects to the catheter.
5	Catheter connection indicator light	This indicator light is lit up after the catheter connection is successful.
6	[Reset] Distance display reset button	The pullback distance reading is zeroed.
7	[Image] Imaging start/stop button	Start/stop ultrasound imaging by the IVUS Catheter.
8	[Pullback] Pullback start/stop button	Start/stop one automatic pullback exam. When the pullback control button is pressed for the manual pullback, press this button to start or stop data recording.
9	[Speed] Pullback speed adjustment button	Set the speed of the automatic pullback

### 3. Preparation of the IVUS System

#### 3.1. Movement of the IVUS System

Caution must be exercised when moving the IVUS System and rules in this instruction manual must be strictly adhered to.

##### 3.1.1 Device transportation

- Please comply with the symbol requirements on the packaging during transportation.
- Do not transport goods that may damage the device together with it.
- Do not place any heavy objects on the packaging.
- When using a vehicle for transportation, ensure that the vehicle can bear the weight of the device components and passengers. If lifting equipment is used, ensure that the equipment can bear the weight of the device.
- Do not load or unload this device under sunlight, rain, or snow without any shelter.

##### 3.1.2 Device movement

- Before moving the device, please ensure that the device is turned off and the power cord has been disconnected.
- Before moving the device, ensure that the casters are unlocked. (refer to 2.2.2. “casters” for unlocking).
- Before moving the device, ensure that the IVUS Catheter is disconnected to avoid pressing on the power cable of the pullback unit.
- Please move the cart using the handle only.
- Please be careful and slow when crossing uneven elevator entrances, steps, or other obstacles and avoid tilting.

#### **Note:**

*Transportation or movement in conditions other than the aforementioned conditions may result in a tilting hazard, which may cause injury to the patient or operator and damage the equipment.*

#### 3.2. Installation of the IVUS System

##### **Attention**

To ensure that the IVUS System can work normally, please read this instruction manual and the instruction manual of the IVUS Catheter first. Install the device following the requirements.

##### 3.1.3 Initial inspection

Before opening the packaging, inspect the package to see if there is any damaging. If the packaging is damaged, request damage compensation from the carrier before repackaging.

Carefully open the packaging and take out the IVUS System and accessories. Check if the device is complete based on the following packing list, check whether the optional parts and accessories are correct, and check if there are any damages on the device, the optional parts, and the accessories. Please contact the supplier immediately if there is any problem.

**Packing list:**

No.	Description	Quantity (EA)
1	VivoHeart Intravascular Ultrasound Imaging System	1
2	Pullback unit	1
3	Printer (optional)	1
4	Printing paper (roll $\phi 50$ mm*110mm), (Optional)	1
5	Three-pin power cord (3.0m, $3 \times 1.5$ mm <sup>2</sup> )	1
6	Console quality certificate	1
7	Instruction manual	1
8	Cross screwdriver	1
9	Allen key with ball end (4 mm)	1

**3.1.4 Connection of alternating current power cord**

Use the three-pin power cord that comes with the VivoHeart Intravascular Ultrasound Imaging System. Insert one end of the power cord into the main power supply port at the back of the IVUS System and insert the other end of the power cord into a grounded flat 3-pin socket. When connecting, the power cord should be easily pulled out from the power supply port and the power supply socket. Press the main power switch at the back of the device to connect to the power supply of the IVUS System, the application software of the IVUS System will prepare the device for clinical use.

The IVUS System can be connected to the power supply for a long period. This enables the system to be in standby mode, thereby shortening the total time for ultrasound imaging. The IVUS System uses a Linux-based operating system. We recommend that the system be restarted manually periodically (connecting to the power supply after disconnecting) to achieve the best performance.

**Note:**

*The rated working voltage and frequency of the VivoHeart Intravascular Ultrasound Imaging System are 100-240V, 50/60Hz, respectively.*

*Connect the power cord to a socket dedicated to the hospital to ensure the reliability of the power connection, and to prevent abnormal power interruption during measurement and affecting imaging.*

**3.3. Power on and off****Turning on the device**

Ensure that the power cord is connected and the main power supply on/off button has been turned on before pressing the console switch at the front of the cart to initiate the IVUS System. The screen will remain black for a few seconds before lighting up and the console system login interface will appear.

**Shutting down the device**

Click the mouse to shut down the system software of the IVUS System before pressing the power supply switch at the back of the IVUS System to turn off the device.

**Caution**

Electricity is still present in the power cord and the console of the IVUS System after the device is turned off. Remove the power cord plug from the socket to completely disconnect the power.

Product is under review by Ministry of Health Indonesia

## 4. Use of the IVUS System

### 4.1. Overview

Two operators are required when using the IVUS System for diagnosis: a sterile operator and a non-sterile operator. All steps requiring contact with the IVUS Catheter or exterior of the pullback unit must be carried out by the sterile operator. All steps involving direct contact with the keyboard or mouse must be carried out by the non-sterile operator.

### 4.2. Preparation for usage of the pullback unit and IVUS Catheter

Before using, ensure that the pullback unit, pullback unit sterile bag, and TrueVision Intravascular Ultrasound Imaging Catheter (IVUS Catheter) are prepared.

#### 4.2.1 Preparation of the pullback unit

- 1) Insert the pullback unit into the sterile bag of the pullback unit. Refer to the instruction manual for the sterile bag of the pullback unit attached to the IVUS Catheter.

**Note:**

**Two people are required to insert the pullback unit into its sterile bag. One person should be inside the sterile zone and one person outside of the sterile zone.**

- 2) Insert the pullback unit sterile bag buckle into the pullback unit.
- 3) Slide the pullback unit back and forth on the pullback plate to ensure that the pullback unit sterile bag does not obstruct the gear mechanism.

① The sterile staff opens the sterile bag

② The non-sterile staff holds the pullback unit and inserts it into the sterile bag

③ The sterile bag buckle is inserted into the pullback unit.

④ Ensure that the sterile bag does not obstruct the pullback unit

Figure 8. Illustration of inserting the pullback unit into its sterile bag

#### 4.2.2 Connect the catheter to the pullback unit

**Note:**

Connection of the catheter to the pullback unit must be performed by the operator in the sterile zone. For detailed information, refer to the instructions attached to the sterile bag of the pullback unit.

- 1) Hold the purple proximal strain-relief tube of the catheter and ensure that the hub and the white clasp of the pullback unit sterile bag are in the same direction;
- 2) Insert the hub into the port of the pullback unit with the pullback unit sterile bag until the case is locked into place with a clear audible locking sound;
- 3) When the indicator light on the pullback unit lights up (green) and a message displays on the screen showing “Catheter connected”, the catheter is connected successfully;
- 4) Insert the catheter clamp at the front groove of the pullback unit. A clear locking sound will be heard once insertion is completed.



Figure 9. Illustration of catheter connection to the pullback unit

① The catheter indicator light on the pullback unit is light up

② The display interface pops up "Catheter connected"

Figure 10. Illustration of successful catheter connection

#### 4.2.3 Catheter replacement

- 1) Press the catheter release button in front of the pullback unit before removing the catheter with some force.
- 2) Insert a new catheter into the catheter port of the pullback unit until it is locked into place and the message “Catheter connected”, as well as the type of catheter is displayed on the screen.

#### 4.2.4 Connect the pullback unit cable to the cart

- 1) Find the pullback unit port at the front of the cart system and pay attention to the red dot on top of the jack. At the same time, pay attention to the corresponding red dot at the pullback unit cable connector.
- 2) Match the red dot on the pullback unit cable connector to the red dot at the port of the pullback unit at the front of the cart, and gently insert the cable connector into the jack of the port of the pullback unit connector at the front of the cart until the knurled tube of the pullback unit connector is locked into place. At this time, the pullback unit cable is connected to the IVUS System.
- 3) To disconnect the pullback unit cable from the cart system, gently pull the knurled tube outwards to unlock the tube and then directly pull out the connector from the jack.

Figure 11. Connecting the pullback unit cable and the cart


Figure 12. Illustration of successful connection

### 4.3. Starting the measurement

Before starting measurement on a patient, you must get to know the entire process of measurement: from logging into the IVUS System to creating patient data, imaging, viewing of data, and filing to shutting down the console. The entire procedure can be summarized as:

- 1) Logging into the IVUS System
- 2) Creating the patient data
- 3) Imaging
- 4) Viewing data
- 5) Viewing archive
- 6) Shutting down the console system

#### 4.3.1 Logging into the IVUS System

After the login page is displayed on the screen, please get prepared to log into the IVUS System and start exam.

#### Attention

When logging into the system for the first time, the user name is displayed as “**Admin**” and the password is blank. Please contact the manufacturer for the password. The system administrator will designate a user name, which cannot be changed by the user. The system administrator can add or delete users at any time on the settings menu and change the password to log into the IVUS Imaging System.


Figure 13. Login interface

- 1) Input your user name in the **User** input box.
- 2) Input your password in the **Password** input box.
- 3) Click the **Login** button to enter the real-time interface.

#### Attention

It is recommended that a unique user ID (user name and password) be used to ensure that your settings can always be used and not accidentally changed by a third-party. If you do not have a unique user ID, contact your system administrator.

The display after entering the real-time interface is as follows:


Figure 14-1. Real-time interface


Figure 14-2. Real-time interface

Region No.	Description
1	Live - Real-time display of images during the exam
2	Review - Edit and view current exam information
3	Pullback Speed - Set the pullback speed for the pullback unit
4	Image depth - Set the depth of the image display, the imaging depth can be adjusted
5	Pullback method- Set the pullback method for the pullback unit. Options include automatic and manual pullback (shutting down automatically).
6	Image - Start image acquisition
7	Record - Pull back the pullback unit
8	Bookmark - Add a bookmark for the current image
9	Graymap - Adjust the graymap of the image display
10	Gain - Adjust the gain of the image display

#### 4.3.2 Creating patient data

Click the mouse to select the real-time interface and use the keyboard for manual entry or selection at the right patient interface:

- [Name]
- [ID]
- [Date of birth]
- [Gender]
- [Anonymity]
- [Physician]
- [Exam]
- [Vessel]

-- [Phase]

Information can be filled in before, during, or after a pullback (View page) at any time. After the information is filled in, click at any place outside of **Create new patient** to save the information. Click on the newly created information in the archive for viewing. Click **View** and the information of the patient will be shown on the right side of both the Live interface and View interface. The **Anonymity** function can be enabled or disabled by single-clicking to show or hide the patient name, patient ID, date of birth, or gender.


Figure 15. Patient interface

#### Attention

The patient ID is unique and cannot be re-entered or blank.

#### 4.3.3 Real-time imaging

Real-time imaging can be conducted by operating the software interface or the pullback unit. Before imaging, ensure that the catheter is connected normally to the pullback unit. Imaging is classified as automatic imaging and manual imaging. By default, it is automatic imaging.

#### Regarding interface imaging display

The ultrasound image will be shown on the display of the cart system, and there are two imaging display layouts for your selection: Imaging display with hidden longitudinal section view and imaging display with active longitudinal section view. By default, it is an imaging display with a hidden longitudinal section view.

#### Imaging display with hidden longitudinal section view

When the longitudinal section view is hidden, the transverse section ultrasound view will be displayed on the screen, as shown in the following figure.


Figure 16. Imaging display interface of hidden longitudinal section view

Region No.	Description
1	Full-length transverse section ultrasound image display
2	Patient name, patient ID, exam name, blood vessel, and stage information
3	Frame slide
4	Record exam time and catheter model number
5	Pullback distance and pullback frame number information
6	Activate the longitudinal section view button
7	Activate the double-transverse sections view button
8	Print - Print the image in the current frame
9	Export - Export current exam image
10	System menu - Includes settings, archive, logout, and shutdown functions

### Imaging display with active longitudinal section view

In the basic imaging procedure, the transverse section ultrasound image will be displayed on the display screen. Clicking the Activate longitudinal section view button during imaging will cause the longitudinal section view to be simultaneously displayed on the screen. This is shown in the following figure:


Figure 17. Imaging display interface with longitudinal section view enabled

Region No.	Description
1	Transverse section ultrasound image display under the longitudinal section view
2	Patient name, patient ID, exam name, blood vessel to be examined, and the stage information
3	Longitudinal section rotary guide ruler
4	Longitudinal section ultrasound image
5	Frame slide
6	Pullback frame information and pullback distance
7	Exam recording time and catheter model number
8	Exam name

#### Attention

The **Activate longitudinal section view** button only changes the image view. The longitudinal section information of the longitudinal section view function can only be obtained after the pullback unit is started.

#### 4.3.3.1 Manual imaging

- 1) After the catheter tip is inserted into the distal end of the target region, press the Image button on the pullback unit. The display screen will show the transverse section image and Image depth (4, 5, 6, 7, 8 mm) can be selected for depth adjustment.
- 2) Hold the handle of the pullback unit and press the pullback control button on the pullback unit. Slowly move the pullback unit to the right to pull back the catheter imaging assembly. To save an image during manual imaging, click on the Pullback button on the software interface or the

Pullback icon on the pullback unit.

- 3) Press the Image button on the pullback unit and the motor will stop rotating when manual imaging is completed.

#### 4.3.3.2 Automatic imaging

##### Automatic imaging using the pullback unit

- 1) After the catheter tip is inserted into the distal end of the target region, in the real-time interface, adjust **The pullback speed** (0.5, 1, 5, 10 mm/s), **Image depth** (4, 5, 6, 7, 8 mm), **Graymap**, and **Gain** on the real-time interface as needed.
- 2) Click the **Acquire** button with the mouse or press the **Image** button on the pullback unit, and the system will start automatic imaging based on the parameters just set. After imaging is started, the background of the **Acquire** button on the real-time interface changes to black.
- 3) Click the **Record** button with the mouse or press the **Pullback** button on the pullback unit, and the pullback unit will start pullback based on the parameters just set. After the pullback is started, the background of the **Record** button on the real-time interface changes to black.
- 4) Click the **Record** button with the mouse or press the **Pullback** button on the pullback unit, and the pullback unit will stop pulling back.
- 5) Click the **Acquire** button with the mouse or press the **Image** button on the pullback unit to complete the entire imaging process.

The acquisition can also be stopped by skipping step (4) and directly performing step (5).

#### 4.3.3.3 Gain adjustment

##### Regarding gain adjustment

Click the **Acquire** button with the mouse or press the **Image** button on the pullback unit, click the blue cursor of **Gain** and shift left or right with the mouse, or use the left and right arrow keys at the lower right corner of the keyboard to move, the image brightness will change as the adjusted value changes. Gain can only be adjusted before and during acquisition. Gain cannot be adjusted during and after the pullback.


Figure 18. Gain adjustment interface

#### 4.3.3.4 Graymap (brightness adjustment)

##### Regarding graymap adjustment

Clicking the **Acquire** button with the mouse or pressing the **Image** button on the pullback unit, clicking the blue cursor of **Graymap** and shifting left or right with the mouse, or using the left and right arrow keys at the lower right corner of the keyboard to move, the image contrast will change as the adjusted value changes. Graymap can be adjusted before and during acquisition and pullback.


Figure 19. Graymap adjustment interface

#### 4.3.3.5 Quick pullback bookmarking

During pullback, click **Bookmark** with the mouse to quickly label frames of interest for rapid localization during playback.

#### 4.3.4 Viewing data

After imaging is completed, the system automatically changes to the view page to view the latest patient data. Click **Archive** for earlier patient data.

##### 4.3.4.1 Viewing the exam

##### Regarding viewing the exam

The function of viewing previous exams provides greater flexibility for understanding acquired data.

When viewing the exam, you can:

- Quickly navigate to any bookmark added during recording.
- Create a new bookmark or navigate to any single frame in all exams to view static images.
- Perform length or area measurement.
- Add annotation.
- Select any frame as the start/stop point or center point for replay or dynamic play functions.


Figure 20. Viewing the exam

Region No.	Description
1	Review - Edit current exam information
2	Switch exam - Switch to another exam of the current patient
3	Exam information - View and edit exam information of the current patient

Region No.	Description
4	Print - Print the image in the current frame
5	Export - Export current exam image
6	Frame slide
7	Graymap adjustment - Adjust the graymap of the image

### Selecting the exam to be viewed

When one or more exams are performed for one patient on the same day, the upper left corner of the display screen will show the exam list. Move the mouse pointer to the exam to be viewed and click to select that exam. The screen will display the ultrasound image of the currently selected exam.


Figure 21. Viewing the list interface

### Deleting the currently selected exam

The upper right corner of every currently selected exam will show a delete icon. Click the **Delete** button with the mouse and select the **OK** button in the popup dialog box to delete the currently selected exam.

#### 4.3.4.2 Editing exam information

### Regarding editing exam information

Exam names are used to help in identifying specific exams for specific patients. During viewing the exam, you can access the exam list before selecting the exam(s) for viewing. Under the default settings, the system will name exams as exam 1, exam 2, exam 3, etc. If required, you can define individual exam labels, as shown in the following figure:


Figure 22. Editing exam information interface

### Editing exam information

Click the mouse to select exam information in the view interface. Use the keyboard to input the name you want in the exam name field, which can be used to modify the name of the current exam. Following that, select the corresponding physician, surgical procedure, and blood vessel type. Click on the software interface outside the input field to store the just edited exam information.

#### 4.3.4.3 Exporting the viewed exam

Click the mouse to select the  button on the **Archive** interface, or right-click patient information in Archive, to click the **Export** button. Following that, select the export method (DVD/CD, network, portable storage media) and data type (DCM, AVI, BMP, JPG) in the popup dialog box, then click the **OK** button to export the currently viewed exam.

#### 4.3.4.4 Printing data

Select the specific transverse or longitudinal section ultrasound image frame for printing. Click the mouse to select the  button and the printer will promptly print out the selected ultrasound image frame.

#### 4.3.4.5 Image annotation

##### Regarding text annotation

You can directly add text to static transverse or longitudinal section images. If images are not clear

at the first viewing, adding text annotation will aid in subsequent image viewing. You can select pre-defined text from **Comment**  under the **Review** menu. Otherwise, you can self-define text to provide meaningful descriptions of specific images and patients' conditions or record information other than normal clinical observation.


Figure 23. Image annotation interface

### Adding annotation

Click the **Comment**  button under the **Review** menu to add annotation on the screen.

- 1) Use the mouse to manipulate the navigation pin to navigate to the specific frame to which you want to add an annotation. This frame can be either a full-length cross-section ultrasound image or a combination of transverse and longitudinal section ultrasound images.
- 2) Click the mouse to select the **Comment**  button under the **Review** menu.
- 3) Move the mouse pointer to the spot on the image where you want to annotate. Key in the annotation description that you want to add using the keyboard or select pre-defined text at the drop-down list beside the annotation button.

### Moving the annotation

- 1) Move the mouse pointer to the annotation text box that you wish to move and the text box will be highlighted.
- 2) Hold down the left mouse button to drag the annotation text box to the new location you want.

### Amending the annotation

- 1) Move the mouse pointer to the annotation text box that you wish to amend and the text box will be highlighted.

- 2) Click the mouse to allow the annotation text box to enter the editing mode and key in the amended annotation using the keyboard.

#### Deleting the annotation

- 1) Move the mouse pointer to the annotation text box that you wish to delete and the text box will be highlighted.
- 2) Click the mouse to allow the annotation text box to enter the editing mode and right-click and press the **Delete** button to delete the annotation you want to delete.

#### 4.3.4.6 Bookmark

##### Regarding the bookmark

A bookmark refers to capturing a snapshot of all information in the cross-section ultrasound image region for future viewing.


Figure 24. Bookmark interface

##### Creating a bookmark

- 1) Use the mouse to manipulate the navigation pin to navigate to the frame where you want to create a bookmark.
- 2) Click the mouse to select the **Bookmark** button under the **Review** menu. The display screen will show “Bookmark successfully created”.

##### Deleting a bookmark

- 1) Move the mouse pointer to the upper right corner of the thumbnail and a delete button will appear. Click the **Delete** button to delete the bookmark.

#### 4.3.4.7 Dynamic play

##### Regarding the dynamic play function

The dynamic play function enables you to select any frame in the exam and then select a series of frames in the range before and after that frame, so that all frames within that range are shown in a scanning loop (from the distal end to the proximal end and back again).

The dynamic play function can highlight the difference between fast and slow moving acoustic reflectors on the image. Vascular wall motion (pulsation and rotation) and vascular plaque reflection are considered "slow" motion (frame by frame) while vascular plaque reflection significantly changes between frames. By playing a short sequence of continuous frames back and forth again and again, an observer can easily detect these differences. A smaller range of frame loops can minimize the interference of vascular wall motion artifacts. Highlighting these differences will enable the observer to easily identify boundaries between blood flowing and not flowing.

##### Starting and stopping dynamic play function

- 1) Click the mouse to select the Dynamic play  button under the Review menu, and a total of 15 frames before and after the frame (as the center frame), where the cursor locates, will be played in a loop.
- 2) When clicking the mouse to select the  button again, the dynamic play will stop.

##### Attention

You can set the number of frames shown in the Settings interface if you want the loop to show 3, 5, 7, 9, 11, 13, and 15 frames.

#### 4.3.4.8 Image comparison

##### Regarding image comparison

Image comparison refers to the comparison of transverse section images from 2 different frames.


Figure 25. Image comparison interface

### Activate image comparison

- 1) Click the mouse to select the **Activate view of double transverse sections** button, and the image region will show transverse section images from the left and right frames, and two frame sliders will show on the longitudinal section image.
- 2) Click the mouse to select one transverse section image from the left frame and move the mouse pointer to the corresponding frame slider. Click to activate the frame slider, move it to the desired location, and click the mouse again to fix the frame slider.
- 3) Click the mouse to select one transverse section image from the right frame and move the mouse pointer to the corresponding frame slider. Click to activate the frame slider, move it to the desired location, and click the mouse again to fix the frame slider.
- 4) Compare the transverse section images selected from two frames.
- 5) The length between the two frames will be displayed on the longitudinal section image.

### Disabling image comparison

When other statuses are active, the image comparison state will automatically be disabled.

#### 4.3.4.9 Longitudinal section display

Display or hide the longitudinal section image.


The rotating ruler can be rotated at the transverse section to show different longitudinal sections.

#### 4.3.4.10 Replay, forward, and backward functions

### Regarding replay

Different from the dynamic play function, the replay function allows you to specify a frame at any position of the current exam and continuously play from this frame to the last frame. Please follow the following procedure to set up and operate the replay function.

### Replay the exam

- 1) Select the starting frame that you want to view in a continuous loop.
- 2) In the view interface, click the **Replay**  button with the mouse to start the continuous loop from the beginning frame to the last frame.
- 3) Click the **Replay**  button with the mouse to stop playing the continuous loop.
- 4) Click the **Forward**  or **Backward**  button with the mouse to view the image of the previous or the next frame.

The image displayed on the transverse section can be selected by clicking the mouse on the longitudinal section.

### 4.3.5 Image measurement

Select the specific ultrasound image frames for transverse or longitudinal sections that you want to measure on the viewing exam page. The measurement performed will be related to the selected frame and will be shown on the display screen when viewing this frame.

The VivoHeart Intravascular Ultrasound Imaging System uses length and area as measurement parameters.

- Length - The length measurement error in the same transverse section view is within +/- 10% or +/- 300  $\mu\text{m}$  (whichever is larger). The length measurement error in the longitudinal section view is within +/- 10% or +/- 1.5 mm (whichever is larger).
- Area - Area measurement error  $\leq 1.5\text{mm}^2 + \text{reading} \times 10\%$ .

#### 4.3.5.1 Length measurement

Important information on various morphology appearing in the blood vessel, lumen, or outside of the vessel can be obtained from length measurement. In most situations, the length measurement is used to measure the dimensions of vascular structures in the transverse section view to help confirm plaques or other restrictive occlusions.


Figure 26. Length measurement interface

Region No.	Description
1	Comment - Text description of the current frame image
2	Length - Measurement of diameter length, longitudinal section length in the image
3	Draw Area - Area measurement (curve)
4	Area – Area measurement (stipple curve)
5	Bookmark – Bookmark the current image
6	Replay – Play continuously from the current frame image to the last frame image
7	Backward – View the previous frame image
8	Forward – View the next frame image
9	Dynamic play – The dynamic play range is 3 to 15 frames. Each selection increases or decreases by 2 frames.
10	Bookmark Length– can disable or enable displaying the length to the current frame/bookmark

Starting the length measurement

- 1) Use the mouse to manipulate the navigation pin to navigate to the specific frame of morphology which you want to measure.
- 2) Click the mouse to select the **Length** button under **Review**.
- 3) Move the mouse pointer to the starting point on the image where you want to start the measurement and click to confirm the starting point.
- 4) Continue to move the mouse pointer, and when the pointer reaches your desired end point, click again to confirm the end point. The image will show a straight line, and at the same time, the measured length will also be shown beside the straight line.

- Reset the starting and ending points as needed. The straight line and length of each measurement will be displayed in different colors.

### Changing the end point

If you wish to change the start or end point of a measurement line, please move the pointer to any point. After the pointer changes into a crossed arrow, click to obtain control of that point, and move it to a new position. Click to release that point and fix it to the image displayed.

### Completing the length measurement

- Continue to measure the image diameter length for the specific frame until you have obtained the information you need.
- If required, use the **Forward** and **Backward** buttons under the **Review** menu to measure any other frame.
- After the measurement is completed, you can return to the imaging screen, complete the exam, edit the exam list, or bookmark it according to your needs.
- All measurement information will be saved as part of the exam. You can re-visit these measurements during the exam at any time.

### Deleting the length measurement

- Move the mouse pointer to the straight line on the screen that you want to delete, and the selected line will be highlighted. Right-click the **Delete** icon to delete the selected length.


Figure 27. Longitudinal section length measurement interface

### Starting the longitudinal section length measurement

- Use the mouse to manipulate the navigation pin to navigate to the specific frame of morphology

which you want to measure.

- 2) Click the mouse to select **Activate longitudinal section view button** on top of the screen.
- 3) Click the mouse to select the **Length** button under the **Review** menu.
- 4) Move the mouse pointer to the starting point on the longitudinal section image where you want to start the measurement and click to confirm the starting line.
- 5) Continue to move the mouse pointer, and when the pointer reaches your desired end point, click again to confirm the end line. Part of the length will be shown on the image and the measured length will also be shown in the middle of the measured length of the longitudinal section.
- 6) Reset the starting and ending lines according to needs. The longitudinal section length of each measurement will be displayed in different colors.

#### **Changing the end point**

If you wish to change the starting or ending line of a longitudinal section length measurement, move the pointer to any endpoint. After the pointer changes into a left and right arrow, click to obtain control of that endpoint line, and move it to a new position. Click to release that endpoint line and fix it to the longitudinal section display image.

#### **Completing the length measurement of the longitudinal section**

- 1) Continue to measure the length of the longitudinal section until you have obtained the information you need.
- 2) After the measurement is completed, you can return to the imaging screen, complete the exam, edit the exam list, or bookmark it according to your needs.
- 3) All the measurement information will be saved as part of the exam. You can re-visit these measurements during the exam.

#### **Deleting the length measurement of the longitudinal section**

- 1) Move the mouse pointer to the measurement length on the longitudinal section that you want to delete, and the selected length will be highlighted. Right-click the **Delete** icon to delete the selected length.

#### **4.3.5.2 Area measurement**

A closed curve needs to be drawn using the mouse pointer for manual area measurement.


Figure 28. Area measurement interface (area of a surface plotted by cursor trajectory)


Figure 29. Area measurement interface (stipple area)

### Moving the mouse pointer to plot a curve for area measurement

- 1) Use the mouse to manipulate the navigation pin to navigate to the specific frame of morphology which you want to measure.
- 2) Click the mouse to select the **Area** button under the **Review** menu.
- 3) Move the mouse pointer to the point on the image where you want to start plotting the curve and click to confirm the area starting point.
- 4) Move the pointer to plot the curve until the starting point is reached, and then click to form a

closed curve. At this point, the area, longest diameter, and shortest diameter of the closed curve will be automatically shown in the lower right corner of the image.

- 5) Repeat these steps if required. Peripheral trajectories with different colors will be plotted each time. The measurement for each color will be displayed at the lower right corner of the image and every measurement will show the total area, longest diameter, and shortest diameter.

### **Connecting points to generate a curve for area measurement**

- 1) Use the mouse to manipulate the navigation pin to navigate to the specific frame of morphology which you want to measure.
- 2) Click the mouse to select the **Area** button under the **Review** menu.
- 3) Move the mouse pointer to the starting point on the image where you want to plot the peripheral trajectory and click to confirm the area starting point. Then the screen will show the starting point for the curve.
- 4) Move the mouse pointer along the periphery of the region you want to measure and click at each point. The screen will show each point gradually. Then right-click to close the region for measurement or click the **Complete** icon in the middle. At this point, the lower right corner of the image will automatically show the area, longest length, and shortest diameter measured in the periphery.
- 5) Repeat these steps if required. Peripheral trajectories with different colors will be plotted each time. The measurement for each color will be displayed at the lower right corner of the image and every measurement will show the total area, longest diameter, and shortest diameter.

### **Performing other area measurements**

If there is a need to perform other area measurements in the same exam, use the **Forward** and **Backward** buttons under the **Review** menu to measure any other frame.

### **Editing area measurement**

If you wish to change a part of the area trajectory, use the mouse pointer to drag and redraw the trajectory for this part.

- 1) Move the mouse pointer on the screen to select the trajectory that you want to edit.
- 2) After the mouse pointer changes from an arrow to a cross, click to redraw the part of the curve based on your needs.
- 3) Click again after obtaining your desired trajectory to stop plotting this part of the curve.

### **Deleting area measurement**

- 1) Select the circle drawn and right-click for the **Delete** icon to appear.
- 2) Click on **Delete** and select **OK** in the popup dialog box and that area measurement will be deleted.

## 4.3.6 Viewing archive data

### 4.3.6.1 Viewing records in the archive

Click the mouse to select the **Archive** menu. The archive interface that pops up will show patient information. Click the mouse to select the patient record and click the **Review** button to view all exams of that patient; right-click on the patient for viewing and click on the view option to view all exams of that patient; or double-left-click on the patient to enter the view page. The available space of the hard disk drive is displayed on the left bottom of the screen. The capacity is 1T (1024GB).


Figure 30. Archive interface

### 4.3.6.2 Creating new exams for existing patients

Click the mouse to select the patient that you want to create a new exam, then double left-click to enter the view page and click real-time in this state to create a new exam for the current patient.

### 4.3.6.3 Editing and deleting records in the archive

Click the mouse to select the record in the archive for editing. After that, move the mouse pointer and click on the **Edit** button or right-click to select the edit option. The patient's name, patient ID, date of birth, and gender can be edited in the popup dialog box.

Click the mouse to select the record in the archive for deletion. After that, move the mouse pointer and click on the **Delete** button or right-click to select the delete option. Select the **OK** button in the popup dialog box to delete this exam record.

### **Attention**

Records in the archive are listed by exam date and patient information (there may be multiple exams in one day). The exams on the Review page are listed by each exam.

#### **4.3.6.4 Exporting records in the archive**

Click the mouse to select the record in the archive for exporting. Then move the mouse pointer and click on the **Export** button or right-click to select the export option. Select the export method (DVD/CD, network, portable storage media) and data type (DCM, AVI, BMP, JPG) in the popup dialog box. After that, click the **OK** button to export all exam records of the patient for that date.

#### **4.3.6.5 Patient search**

Click the mouse to select the search input box in the archive and enter the patient information to be searched. The archive page will only show information about the searched patient.

#### **4.3.7 Shutting down the console system**

Click the mouse to select the **Shutdown** button under the **main menu** option at the right upper corner of the real-time or view interface. Select **OK** in the popup dialog box to shut down the console.

## 5. Settings

Click the mouse to select the **System Menu** icon on the right top of the screen in the real-time interface to enter the Settings interface. There are 8 settings, namely General, Account, Doctor, Comment, Measurement, DICOM, and Log, and factory settings (by default it is hidden).

### 5.1. General settings


Figure 31. General settings interface

In this interface, you can:

- Set the language used for the display screen: Languages available are Chinese and English. System default settings: Chinese. Click the mouse to select English and the system will pop up a prompt stating that the language switch will be effective after restarting the system.
- Set date and time: Input the current date and time. Use the mouse to click the date or the up and down arrows on the right side of the date to change the date or time as you desired. Click to exit the settings page and the modified date or time takes effect.
- Set the standby time: Options available are 15 minutes, 30 minutes, 1 hour, and never. Default system setting: Never.
- Set the hospital name: Input the hospital name
- Set display mode: turn on and off
- Set the dynamic play speed: The dynamic play speed range is 3 to 15, and each selection increases or decreases by 2.
- Image mode: Available image modes: standard, high-resolution, and sharp-lumen image modes.
- Set hotkeys F1 and F2.
- Display the software version.
- Set reset to factory settings: By clicking the **Restore Factory Settings**, the system will reset parameters to factory settings on the premise of preserving patient data.

## 5.2. Account (user) settings


Figure 32. Account settings interface

### Creating a new user

Click the mouse to select the **Account** menu in the General settings interface, then select the **New user** button in the account settings interface. After inputting the user name and password in the popup dialog box and confirming the password, select the **OK** button to create new user information.

### Editing the existing user information

Click the mouse to select the existing user for information editing. After selecting the **Edit** button and re-entering the name and password in the popup dialog box and confirming the password, the existing information can be re-edited.

### Deleting the existing user information

Click the mouse to select the existing user for deletion and select the **Delete** button, then select the **OK** button in the popup dialog box to delete that user information.

### Account permission

There are two types of users including the administrator and ordinary users, with different account permissions.

- Administrator account: possesses all functional permissions, such as creating and deleting ordinary users, modifying all user passwords, which requires a password to enter the system operation interface;
- Ordinary account: possesses all functional permissions except creating and deleting ordinary

users, and also can modify current user password, which requires a password to enter the system operation interface;

User names and passwords are needed for all users to log in.

### 5.3. Doctor settings


Figure 33. Doctor settings interface

#### Creating a new doctor

Click the mouse to select the **Doctor** menu in the settings interface, then select the **Create new doctor** button in the doctor setting interface. Input the doctor's name and remarks in the popup dialog box, and select the **OK** button to create new doctor information.

#### Editing existing doctor information

Click the mouse to select the existing doctor for information editing. Then select the **Edit** button and re-enter the doctor's name and remarks in the popup dialog box and select the **OK** button to re-edit the existing doctor information.

#### Deleting the existing doctor's information

Click the mouse to select the existing doctor for deletion and select the **Delete** button, then select the **OK** button in the popup dialog box to delete the existing doctor information.

## 5.4. Comment settings


Figure 34. Comment settings interface

### Creating new annotation

Click the mouse to select the **Comment** menu in the settings interface, then select the **Create Annotation** button in the doctor annotation interface. Input the annotation name in the popup dialog box, and select the **OK** button to create a new comment. New annotation information will be shown within the list of optional annotation frames on the left.

### Editing the existing annotation information

Click the mouse to select the existing annotation for information editing. Select the **Edit** button and re-enter the annotation information in the popup dialog box and select the **OK** button to re-edit the existing annotation information.

### Deleting the existing annotation information

Click the mouse to select the existing annotation for deletion and select the **Delete** button, and select the **OK** button in the popup dialog box to delete the existing annotation information.

### Moving the annotation

Annotations in the optional annotation frame on the left can be moved to the current frame on the right by using the **Push into Current** button. All annotations shown in the current annotation frame on the right will be displayed in the view interface for selection. Annotations in the current annotation frame on the right can also be moved to the left optional annotation frame using the **Pull out of Current** button for editing and deletion.

## 5.5. Measurement settings


Figure 35. Measurement settings interface

## 5.6. DICOM settings


Figure 36. DICOM settings interface

In this interface, you can carry out local settings and also add, edit, delete, connect, and disconnect from the server.

## 5.7. Log settings

This log module is mainly used by the maintainer to locate problems.

## 5.8. Factory settings

The factory settings interface is available only to maintenance personnel.

# 6. Technical Specifications

## 6.1. Performance requirements

No.	Performance indicator	40 MHz system	60 MHz system
1	Acoustic working frequency	40 MHz, error $< \pm 15\%$	60 MHz, error $< \pm 15\%$
2	Imaging radius	$\geq 6$ mm	$\geq 4$ mm
3	Axial resolution	$\leq 50$ $\mu\text{m}$ (depth: 3 mm) $\leq 50$ $\mu\text{m}$ (depth: 1.5 mm)	$\leq 40$ $\mu\text{m}$ (depth: 2 mm) $\leq 40$ $\mu\text{m}$ (depth: 1.5 mm)
4	Lateral resolution	$\leq 600$ $\mu\text{m}$ (depth: 3 mm) $\leq 300$ $\mu\text{m}$ (depth: 1.5 mm)	$\leq 300$ $\mu\text{m}$ (depth: 2 mm) $\leq 250$ $\mu\text{m}$ (depth: 1.5 mm)
5	Pullback direction resolution	$\leq 600$ $\mu\text{m}$ (depth: 3 mm)	$\leq 500$ $\mu\text{m}$ (depth: 2 mm)
6	Longitudinal geometric position accuracy	$\leq 10\%$ (within the scope of imaging radius)	
7	Transverse geometric position accuracy	$\leq 10\%$ (within the scope of imaging radius)	
8	Pullback direction geometric position accuracy	$\leq 10\%$	
9	Image geometric distortion	$\leq 10\%$ (within the scope of imaging radius)	

## 6.2. Device classification and regulatory compliance

Type of Protection Against Electric Shock	Class I equipment
Degree of protection against electric shock	Defibrillation-proof CF-type application parts 
Operating voltage	100-240V, 50/60Hz
Rated Input Power	300VA

Degree of Protection Against Liquid Entry	IPX0
Operation Mode	Intermittent loading and continuous operation (The pullback unit continuously operates for 10 minutes at a pullback speed of 0.5 mm/s, 1.0 mm/s, and 5.0 mm/s and rests for 10 minutes, and operates for 3 minutes at a pullback speed of 10.0 m/s and rests for 5 minutes).
Classified by the degree of Safety in the Environment Containing Flammable Anesthetics	Equipment that cannot be used in flammable anesthetic gas mixed with air or flammable anesthetic gas mixed with oxygen or nitrous oxide
Standards Complied	GB 9706.1-2007, GB9706.9-2008, GB/T 14710-2009, YY1420-2016, and YY0505-2012

### 6.3. Environmental conditions

Operating Environment Temperature	+10 °C to +35 °C
Operating Environment Humidity	30 % RH to 75 % RH (no condensation)
Atmospheric Pressure of Operating Environment	525 mmHg to 795 mmHg (700 hPa to 1060 hPa)
Transportation and Storage Temperature	-20 °C to +55 °C
Relative Humidity for Transportation and Storage	10 % RH to 85 % RH (no condensation)
Pressure of Transportation and Storage Environment	525 Hg to 795 mmHg (700 hPa to 1060 hPa)

### 6.4. Port Specifications

DVI	The device uses a standard DVID-D video output port (imaging display) and transmits digital video signal, and its maximum resolution supported by the port is 1920×1200@60 HZ
USB Port	2.0/3.0 protocol, type A port, female connector
Pullback unit Port	Can only connect to the pullback unit, GB 9706.1-2007 isolation port, signal input/output port
Main Power Port	C14 three-pin plug, alternating current power input
Network Port	RJ-45, TCP/IP protocol, transmission speed: 10/100/1000 Mbps, local area network connection (private)

### 6.5. RFID emitter

The catheter connector of the pullback unit has an RFID identification function which acts by the label on the catheter and the coil emission on the pullback unit. The RFID emission frequency is 13.56 MHz, and modulated by ASK/OOK(100%) modulation, with a wavelength of approximately 22m. Except for metallic materials, this wavelength can penetrate most materials but will usually decrease the reading distance. Although the magnet field area of this frequency decreases rapidly, it can generate a relatively uniform reading zone and its effective radiation power is 50 nW.

Warning: Even if other devices comply with pertinent national emission standards, the RFID emitter may still experience interference from other devices.

### 6.6. Cable information

No.	Name of Cable	Length	Shield Cable
1	Interthree-pin power cord	3 m	No
2	pullback unit cable	4 m	Yes

## 7. Basic troubleshooting

### Introduction of the troubleshooting

This section describes basic troubleshooting techniques to confirm, solve, and correct device issues. Please contact an Insight Lifetech sales representative if technical support is required.

For hospital local area network issues related to the DICOM network archive, please contact the IT department of the hospital first for troubleshooting. If the problem cannot be resolved, contact an Insight Lifetech sales representative.

### Troubleshooting method

The user must understand the normal operation of the device to confirm and solve operation problems. The user must adopt rational methods to solve the problem step by step.

### Problem resolution

- 1) Identify and confirm whether the problem is in the imaging operation, software, or hardware.
- 2) If the problem can be solved, you can quickly restore the normal operation of the device.

Please attempt to diagnose and solve the problem using this manual before contacting the Insight Lifetech sales representative. Please always inspect the power source and imaging controls first, and then read this manual for troubleshooting techniques for specific components.

### Messages displayed on the screen

When product malfunction occurs, the malfunction information in the following table will be shown on the display screen. This information is only displayed for a few seconds and will repeatedly appear until the user successfully performs corrective actions.

Malfunction Code	Malfunction Description	Possible Cause	Corrective Action
101	Patient ID cannot be blank	The user performs the pullback operation (by clicking on the pullback button on the pullback unit or the record button on the console) when the patient ID is blank	Re-enter patient ID before repeating the operation
102	Patient ID cannot be repeated	The user performs a pullback operation (by clicking on the pullback button on the pullback unit or the record button on the console) when the newly inputted patient ID is duplicated with the patient ID stored in the archive.	Verify the patient ID. Input the new ID if the failure is due to an error in inputting the new ID; If the newly created exam operation is for the same patient, enter "Archive" and create a new exam under the patient's existing record, and then perform the pullback operation.
104	Data exporting	Data was not successfully	Retry the data export function after

Malfunction Code	Malfunction Description	Possible Cause	Corrective Action
	failure	saved or data export failed since the system is busy	a while
105	Bookmark already exists	The new bookmark overlaps with the old bookmark at the longitudinal axis	Adjust the bookmark location and re-create the bookmark
107	Bookmark creation failure	Bookmark was created at the wrong position or the bookmark creation failed since the system is busy	Re-perform bookmark creation operation. Contact the manufacturer if multiple attempts fail.
109	Printing failure	Printing failure caused by no paper inside the printer, incorrect printing paper placing, data transmission error, etc.	Re-perform the printing operation after checking the printer status. If multiple attempts fail, please contact the manufacturer.
110	Viewing exam failure	Data saving failure, or viewing data failure caused by data saving failure	Re-perform the exam or view other exams in the archive to see if the function is normal; re-start the device and repeat exam recording and other operations; please contact the manufacturer for assistance if all these actions fail.
111	The recording frames have reached their maximum number	The number of image frames recorded has reached the preset maximum value of the system	Stop pullback and check the storage status of recorded data. Re-run pullback and recording if necessary.
112	Abnormal data cache[112]	Aging system, decreased performance, inability to quickly save data under high-speed pullback.	Stop pullback, check whether the previous exam data is normal, and contact the manufacturer for assistance.
114	Data storage failure [114]	The inability to save data is normally caused by abnormal system database or insufficient hard disk space	Check the data storage status and if there is sufficient space; if there is sufficient space, re-perform the pullback operation; please contact the manufacturer for assistance if multiple attempts fail.
116	Factory default reset failure	Database configuration parameter restoration failure	Re-perform reset to factory settings. Please contact the manufacturer for assistance if multiple attempts fail.
117	Abnormal image data acquisition [117]	The pullback unit rotary encoder abnormality or abnormal data transmission	Please contact the manufacturer for assistance.
125	Failed to connect to the DICOM server	DICOM settings error or failure in the connection	Please contact the manufacturer for assistance.
127	Failed to connect to the Worklist server	DICOM settings error or failure in the connection	Please contact the manufacturer for assistance.
140	Insufficient disk space, please clear the disk in time!	Insufficient disk space	Clear the disk

Malfunction Code	Malfunction Description	Possible Cause	Corrective Action
141	Failed to restore to the factory settings	The default database is missing, and the document is damaged.	Please contact the manufacturer for assistance.
143	Abnormal in the printer system	Abnormal in the printer system	Please contact the manufacturer for assistance.
144	The printer paper is not installed properly	The printer paper is not installed properly	Input the print paper
145	The printer hatch is not closed	The printer hatch is not closed	Close the printer hatch
301	Device abnormality [301]	PCIE disconnection caused by the console PCIE device missing or damaged.	Please contact the manufacturer for assistance.
302	Device abnormality [302]	Console PCIE device abnormality, DMA flow control credit runs out	
304	Device abnormality [304]	pullback unit power supply error caused by damaged pullback unit, or PCIE device malfunction in the console.	
305	Device abnormality [305]	Send timing error	
306	Device abnormality [306]	Emission configuration timing exceeds the scanning cycle	
307	Device abnormality [307]	Receiving configuration timing exceeds the scanning cycle	
308	Device abnormality [308]	Malfunctioning FPGA board in the console, an abnormal system temperature alarm	
309	Device abnormality [309]	Abnormal Vccint power supply alarm	
310	Device abnormality [310]	Abnormal Vccaux power supply alarm	
311	Device abnormality [311]	Abnormal Vccbram power supply alarm	
402	The pullback unit is not connected	pullback unit disconnected from the main imaging unit	Check if the pullback unit is connected to the main imaging unit
404	Catheter not connected	The catheter is disconnected from the pullback unit	Pull out the catheter and reinsert it
405	Unable to identify the current catheter	The batch number burned on the catheter chip does not match protocol requirements or the catheter of this product. The aging pullback unit is unable to identify the catheter.	Change the catheter. If the malfunction is still present, please contact the manufacturer for assistance

Malfunction Code	Malfunction Description	Possible Cause	Corrective Action
406	pullback unit connection timed out [406]	Abnormal serial port or system is busy.	1. Remove and re-insert the pullback unit; 2. Restart the device; 3. Contact the manufacturer for assistance
501	Even the clutch is pressed, the automatic pullback cannot be carried out	Click the pullback unit button of the pullback device or record button in the interface when the clutch is pressed.	Inspect the clutch and carry out the pullback operation again. Please contact the manufacturer if multiple attempts fail.
502	Pullback has reached the furthest position	The pullback unit has reached the maximum pullback position;	Stop pullback and check the storage status of recorded data. Re-run pullback and recording if necessary.
503	Catheter not connected	The Image or Pullback button on the pullback device was pressed when the catheter was not connected.	Check if the catheter is connected
506	Abnormal rotating motor temperature [506]	High-speed rotation was carried out for a long duration. High motor temperature alarm	Turn off the pullback unit and rest for more than 5 minutes before use. If the warning still appears after rest, increase the rest duration. Please contact the manufacturer if the warning still appears after multiple rests.
507	Rotating motor initialization is not completed	Rotating motor initialization is not completed	Check whether the display screen on the pullback unit motor is lit up. If not, remove and re-insert the pullback unit
508	Pullback has reached the maximum distance	The pullback device has reached the maximum distance limit during the recording	Push the pullback unit to the front end
600	Pullback unit malfunction [600]	Abnormal pullback unit voltage 26V	Remove and reinsert the pullback unit Restart the device. Please contact the manufacturer if the same problem occurs after multiple attempts
601	Pullback unit malfunction [601]	Abnormal pullback unit voltage A5V	
602	Pullback unit malfunction [602]	Abnormal pullback unit voltage 3V3	
603	Pullback unit malfunction [603]	Abnormal pullback unit voltage 5V	
604	Pullback unit malfunction [604]	Abnormal pullback unit voltage 12V	
605	pullback unit malfunction [605]	Abnormal pullback unit voltage 5V5	
606	Abnormal pullback motor [606]	Maximum current exceeded in pullback motor	Check if the pullback unit is stuck by the sterile bag or other objects. If so, remove the obstruction.
607	Abnormal	Maximum temperature	Otherwise, please contact the

Malfunction Code	Malfunction Description	Possible Cause	Corrective Action
	pullback motor [607]	exceeded in pullback motor	manufacturer for assistance.
608	Abnormal rotating motor [608]	Maximum current or voltage exceeded in rotating motor	Check if the catheter is jammed. If so, please replace the catheter; Remove the catheter and check whether a foreign object has jammed the catheter port on the pullback unit. If so, please remove the foreign object. Please contact the manufacturer for assistance if the malfunction is not due to the above causes.
609	Abnormal rotating motor [609]	Maximum temperature exceeded in rotating motor	
610	Abnormal rotation speed [610]	Unstable rotation speed	
611	Abnormal pullback speed [611]	Damaged pullback motor resulting in unstable pullback speed	Please contact the manufacturer for assistance
612	Abnormal pullback speed [612]	Large pullback speed fluctuation or the external encoder of the pullback unit is damaged	Inspect the pullback process of the pullback unit to see if the sterile bag is jammed. Move the pullback unit to the initial position and perform pullback again; Please contact the manufacturer if this prompt appears after multiple attempts.
613	Abnormal rotating motor [613]	Abnormal rotating motor driver board	Remove and re-insert the pullback unit and restart the device. Please contact the manufacturer for assistance if this prompt still appears.
614	Communication abnormality [614]	The pullback unit was not normally used by the system.	Remove and re-insert the pullback unit and restart the device. Please contact the manufacturer for assistance if this prompt still appears.
615	Failure to open the pullback unit power supply [615]	Poor contact between the pullback unit and the console connection port, cable malfunction, or the pullback unit malfunction.	Remove and re-insert the pullback unit and restart the device. Please contact the manufacturer for assistance if this prompt still appears.
616	Abnormal rotating motor temperature [616]	The operation duration of the pullback unit is too long or, the rotating motor load is too high, which causes the rotating motor to exceed the protective temperature.	1. Stop the motor for 2-5 minutes before restarting; 2. Check if the catheter is jammed. If so, please replace the catheter; 3. Remove the catheter and check whether a foreign object has jammed the catheter port on the pullback unit. If so, please remove the foreign object. 4. If none of the above, contact the manufacturer for assistance.
617	Abnormal rotating motor [617]	The rotary encoder has detected abnormal rotation speed	1. Examine if the catheter is jammed. If so, replace the catheter; 2. Remove the catheter and check whether a foreign object has

Malfunction Code	Malfunction Description	Possible Cause	Corrective Action
			jammed the catheter port on the pullback unit. If so, please remove the foreign object. 3. If none of the above, contact the manufacturer for assistance.

Product is under review by Ministry of Health Indonesia

In addition, prompt messages on the executing status on part of the functions are also displayed, and there is no need to take any actions for such messages, as follows:

Status prompt code	Description of the prompt messages
103	Data exported successfully
106	Bookmark created successfully
108	Printing completed
121	The subject user cannot be deleted.
122	The default annotation cannot be edited.
123	The default annotation cannot be deleted.
124	The DICOM server is connected successfully.
126	The Worklist server is connected successfully.
136	The measurement item already exists.
137	Boundary not recognized
134	Succeeded in switching the Chinese input method
135	Succeeded in switching the English input method
141	The printer is disconnected.
401	Pullback unit connected
403	Catheter connected

## 8. Cleaning and Maintenance

### 8.1. Cleaning

#### Warning

- Before cleaning the device, the power must be turned off and the power cord must be disconnected from the socket.
- Make sure that all parts are dry before reconnecting the power cord.
- Please do not continue using the device if you suspect that liquid has penetrated the IVUS System and contact an Insight Lifetech sales representative or local distributor for technical support.

#### Attention

Always use a clean and non-corrosive soft cloth or paper towel when cleaning the device with detergent or when wiping off residual detergent after cleaning. It is recommended to clean the device once a week.

The following detergents have been verified for the cleaning of the IVUS System:

- Ethanol (75%)
- Isopropyl alcohol (70%)

Follow the following steps to clean the surface of the IVUS System;

1. Turn off the power supply and unplug the power cord.
2. Use a soft cloth to absorb an appropriate amount of cleaning solution, and wipe all external surfaces of the IVUS System, including the display screen.
3. After cleaning, take a new cloth or paper towel to absorb some tap water to wipe off the residual cleaning solution until no obvious residual cleaning solution is left.
4. Place the IVUS System in a cool and well-ventilated place to dry.
5. Reconnect the power cord and press the power button to turn on the IVUS System.

### 8.2. Maintenance

#### Warning

A maintenance plan shall be established in hospitals or medical institutions that use the device to prevent device ineffective and unforeseeable results that may threaten personal safety.

#### 8.2.1 Inspection

A comprehensive inspection shall be performed by the maintenance personnel to ensure the normal operation of the IVUS System before using, after continuously using it for 12 months, or when the device is repaired or upgraded.

The following items are inspected:

- The environment and power supply meet the requirements
- The power cord has no abrasion and its insulation performance is good
- No mechanical damage to the equipment
- The specified accessories are used
- All functions of the IVUS System are good
- Leakage current meets the requirements

### 8.2.2 Maintenance plan

The following tasks can only be completed by the maintenance staff. Please contact the maintenance staff if the following maintenance is required. Please clean the device before testing or maintenance.

Inspection of Maintenance Items	Frequency
Perform leakage current, dielectric strength, and ground resistance tests according to GB 9706.1-2007	The tests are performed once a year. Testing is required if the IVUS System is dropped or electrical repairs were performed.

### 8.2.3 Fuse replacement

The alternating current power supply of the cart system contains a fuse, which is located in the combined filter of the power grid input end. Please inspect if the fuse has blown or if the IVUS System is malfunctioning.

#### Cart system fuse replacement

- 1) Ensure that the main power supply of the cart is at the OFF position (the switch is located at the lower part behind the cart controller) and remove the power cord.
- 2) Use a straight screwdriver to push out the fuse holder.
- 3) Replace the fuse. Please use a high-breaking capacity ceramic fuse, with model No. of T3.15AH250V, dimensions:  $\varnothing 5 \times 20$ mm.
- 4) Reconnect the power cord, turn on the power supply switch, and test if the device is operating normally. If the fuse is blown again, contact the Insight Lifetech sales representative.

### 8.2.4 Printer maintenance

#### Attention

**Please do not attempt to pull out the paper from the printer when the printer panel is closed.**

#### Warning

**Please do not touch the blade at the front of the printer as it may result in physical injury.**

### **Printer paper storage and precautions during usage**

#### **Storage of unused printing paper**

- Please store unused printing paper in a cool and dry environment with a temperature below 30 °C.
- Please do not store unused printing paper near places with volatile liquids. Please do not let the printing paper contact any organic volatile liquids, cellophane paper, or vinyl chloride.

#### **Loading printing paper**


- Be careful when loading the printing paper and ensure that your fingers do not contact the print surface. Medical exam gel, fingerprints, or sweat stains will blur the images.
- When loading the printing paper into the carton, ensure that the front page is up (the lateral surface), otherwise, the image cannot be printed.
- After tearing out the label at the leading edge of the printing paper, pull 15-20 cm of printing paper before printing. Residual label adhesive on the printing paper will decrease the quality of the printed image.
- If previously used printing paper is reused, please pull out 15-20 cm of printing paper before printing. Dust may decrease the quality of the printed image.

#### **Storage of printouts**


- To prevent fading or discoloration in printouts, please store them in a cool and dry environment with a temperature below 30 °C.
- Please store printouts in a polypropylene packaging bag, or between plastic-free papers.
- Please do not store printouts in places that are under direct sunlight or damp.
- Please do not store printouts near places with volatile liquids. Please do not let printouts contact any organic volatile liquids, cellophane paper, or vinyl chloride.
- To prevent fading in printouts, please do not place them above or below diazo paper.
- If printouts are to be placed between other sheets of paper, use double-sided tape or water-based adhesive.
- Please do not burn discarded printing paper.

#### **Loading the printing paper**



- 1) Press the power switch to turn on/off the power supply. If no printing paper has been loaded, the LCD backlight is orange and a message is displayed as “EMPTY”.


- 2) Press the OPEN button to open the printer cover.

- 3) Tear out the label at the leading edge of the printing paper, pull 15-20cm of the printing paper, and load the printing paper on the paper tray in the direction as shown below. The printer does not work normally if the loading direction is wrong.


- 4) Push the printer cover to close it.
- 5) After printing, hold the edges of the printer paper and cut the printer paper upward in an arc. If you pull the printer paper horizontally when cutting the printer paper, the printer paper inside the device will be skewed, which may cause a printing error. If the printer paper is skewed, repeat steps 1) to 4) to reinstall the printer paper correctly.


#### Attention

- Please do not let the printing paper drop or place it face down on the table or other objects when loading or removing printing paper from the tray.

- Ensure that the printing paper is loaded straightly. If the printing paper is crooked, it may result in a paper jam.

- If the paper is not completely pulled out in step (3), press the FEED button so that the left part of the printing paper comes out of the device. Do not manually pull out the printing paper.
- Press OPEN or FEED to cancel the ongoing printing operation.

## 9. Warranty and After-sales Services

### 9.1. Warranty

The warranty for this device covers all device malfunctions due to component failure or manufacturing processes. During the 1-year warranty period, all the malfunctioning parts can be repaired or replaced free of charge.

#### **Manufacturing process and raw materials**

Insight Lifetech guarantees that we use materials and processes that conformed to requirements in the manufacturing of this device. Under normal usage and maintenance, if there are confirmed malfunction reports during the warranty period starting from the date of shipment, the hardware will be repaired or replaced.

#### **Software or firmware**

For software or firmware installed in Insight Lifetech hardware products, if a confirmed malfunction report is received during the warranty period starting from the date of shipment, Insight Lifetech will repair or replace the software or firmware. However, Insight Lifetech cannot guarantee that the use of the hardware, software, or firmware products are not be interrupted or error-free.

**Description: Insight Lifetech hereby declares that the obligations of this warranty do not include transportation and other costs.**

Except for the fuse, there are not any other components that the user can repair by himself/herself. All repairs must be performed by technical staff approved by Insight Lifetech. Insight Lifetech is not responsible for the direct, indirect, or final damage, or delay caused by the following circumstances:

- The components are disassembled, stretched, or re-commissioned
- The device is repaired or modified by staff that was not authorized by Insight Lifetech
- Damage arising from non-normal usage that exceeded the stipulated usage conditions
- Modification or removal of the original serial number or manufacturing label
- Improper operation

### 9.2. After-sales services

If you have any questions during the use of this device, please contact Insight Lifetech or the local distributor immediately.

**After-sales service provider (manufacturer):** Insight Lifetech Co., Ltd.

**Contact address:** Zone E, 3rd Floor, Building 3, Tinwe Business Park, No.6 Liufang Road, Baoan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

**Postal Code:** 518000

**Email:** service@insight-med.com

**Tel:** +86 (755) 23074845

## 10. Acoustic Output Information

The working mode of acoustic output: B automatic scanning mode

Application: Coronary imaging

Note: In the document "Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers", the intravascular ultrasound of coronary angiography is classified to determine the maximum allowable sound output energy. This catheter is not applicable for fetal imaging.

40 MHz system acoustic output table

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum index value		0.11	N/A	0.032	N/A	N/A	N/A
Index component value			0.032		N/A	N/A	N/A
Acoustic Parameters	$p_{ra}$ at $Z_{MI}$ (MPa)	0.68					
	P (mW)		0.16				N/A
	$P_{1 \times 1}$ (mW)		0.16				
	$Z_s$ (cm)			N/A			
	$Z_{bp}$ (cm)						
	$Z_b$ (cm)					N/A	
	$Z_{MI}$ (cm)	0.1					
	$Z_{pii, \alpha}$ (cm)	0.1					
Other Information	$f_{awf}$ (MHz)	39.35	N/A	39.35	N/A	N/A	N/A
	$p_{rr}$ (Hz)	29790					
	$S_{rr}$ (Hz)	29.93					
	$n_{pss}$	1					
	$I_{pi, \alpha}$ at $Z_{pii, \alpha}$ (W/cm <sup>2</sup> )	10.95				N/A	
	$I_{spta, \alpha}$ at $Z_{pii, \alpha}$ or $Z_{sii, \alpha}$ (mW/cm <sup>2</sup> )	19.99					
	$I_{spta}$ at $Z_{pii}$ or $Z_{sii}$ (mW/cm <sup>2</sup> )	27.63	N/A	Fixed	N/A	N/A	N/A
Operating control conditions	$p_r$ at $Z_{pii}$ (MPa)	0.79	N/A	Fixed	N/A	N/A	N/A
	Focus (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
	Depth (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
	Frequency (MHz)	40	N/A	40	N/A	N/A	N/A
Note: N/A means that there is no corresponding intended use for the probe or there is no data to be filled in.							
Uncertainty: sound intensity: $\pm 25.12\%$ , sound pressure: $\pm 12.56\%$ , sound power: $\pm 25.12\%$ , central frequency: $\pm 0.13\%$							

60 MHz system acoustic output table

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum index value		0.18	N/A	0.011	N/A	N/A	N/A
Index component value			0.011		N/A	N/A	N/A
Acoustic Parameters	$p_{ra}$ at $Z_{MI}$ (MPa)	1.17					
	P (mW)		0.055				N/A
	$P_{1 \times 1}$ (mW)		0.055				
	$Z_s$ (cm)			N/A			
	$Z_{bp}$ (cm)						
	$Z_b$ (cm)					N/A	
	$Z_{MI}$ (cm)	0.1					
	$Z_{pii, \alpha}$ (cm)	0.1					
Other Information	$f_{awf}$ (MHz)	59.18	N/A	59.18	N/A	N/A	N/A
	$p_{rr}$ (Hz)	29790					
	$S_{rr}$ (Hz)	29.23					
	$n_{pss}$	1					
	$I_{pi, \alpha}$ at $Z_{pii, \alpha}$ (W/cm <sup>2</sup> )	12.73				N/A	
	$I_{spta, \alpha}$ at $Z_{pii, \alpha}$ or $Z_{sii, \alpha}$ (mW/cm <sup>2</sup> )	12.57					
	$I_{spta}$ at $Z_{pii}$ or $Z_{sii}$ (mW/cm <sup>2</sup> )	17.47	N/A	Fixed	N/A	N/A	N/A
Operating control conditions	$p_r$ at $Z_{pii}$ (MPa)	0.88	N/A	Fixed	N/A	N/A	N/A
	Focus (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
	Depth (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
Operating control conditions	Frequency (MHz)	60	N/A	60	N/A	N/A	N/A
Note: N/A means that there is no corresponding intended use for the probe or there is no data to be filled in.							
Uncertainty: sound intensity: $\pm 25.12\%$ , sound pressure: $\pm 12.56\%$ , sound power: $\pm 25.12\%$ , central frequency: $\pm 0.13\%$							

Terms	Definition	Unit
MI	Mechanical index	n/a

TIS	Soft tissue thermal index	n/a
$p_{ra}$	Peak negative pressure after attenuation	MPa
$f_{awf}$	Acoustic working frequency	MHz
$P$	Output power	Mw
$z$	Point at specific length from the acoustic source on the acoustic axis	cm
$t_d$	Pulse duration	us
$prr$	Pulse repetition frequency	Hz
$p_r$	Peak negative pressure	MPa
$I_{pa,\alpha}$	Mean sound intensity of pulse after attenuation	W/cm <sup>2</sup>
$I_{pi}$	Pulse sound intensity integral	J/m <sup>2</sup>
$I_{pi,\alpha}$	Pulse sound intensity integral after attenuation	J/m <sup>2</sup>

## 11. Electromagnetic Compatibility Information

### Warnings

- The IVUS System shall not be used close to or stacked on other equipment. If it must be used close or stacked on other equipment, observation shall be carried out to verify that it can operate normally under the used configurations;
- Except for the cables sold by the manufacturer of the IVUS System as spare internal parts, the use of accessories and cables other than those specified may result in increased emission of the IVUS System or decreased immunity.
- Combined use of the IVUS System and accessories, transducer, or cable outside of specifications may result in increased emission of the IVUS System or decreased immunity.

### Attention

- The user shall install and use the equipment according to the electromagnetic compatibility information provided in the accompanying documents;
- Portable and mobile RF communication equipment may affect the performance of the IVUS System. When using it, avoid strong electromagnetic interference, such as that generated close to mobile phones and microwave ovens.

### Basic performance

- When real-time imaging acquisition is selected in the IVUS System, the IVUS Catheter will continuously output sound under the default acoustic working frequency with maximum gain and maximum imaging depth set.

Guidelines and Manufacturer's Declaration - Electromagnetic Emission		
The IVUS System (model No.: VH-60) is expected to be used in the following electromagnetic environment. The purchaser or user of the IVUS System (model No.: VH-60) should ensure that the device is used in this electromagnetic environment:		
Emission test	Conformance	Electromagnetic environment - guidelines
Radiofrequency emission GB 4824	1 set	Radiofrequency energy is only used for the internal functioning of the IVUS System (model No.: VH-60). Therefore, its RF emissions are very low and has a low possibility of interfering with nearby electronic equipment.

Radio-frequency emission GB 4824	Class A	The IVUS System (model No.: VH-60) is suitable for use in all establishments other than domestic and those directly connected with the public low-voltage power supply network that supplies buildings for domestic purposes.	
Harmonic emission GB 17625.1	Not applicable		
Voltage fluctuation/ flicker emission GB 17625.2	Not applicable		
Guidelines and Manufacturer's Declaration - Electromagnetic Immunity			
The IVUS System (model No.: VH-60) is intended to be used in the following electromagnetic environment. The purchaser or user of the IVUS System should ensure that it is used in the following electromagnetic environment:			
Immunity test	IEC 60601 test level	Conformance level	Electromagnetic environment - guidelines
Electrostatic discharge GB/T 17626.2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient burst GB/T 17626.4	±2 kV for the power supply lines ±1 kV for input/output supply lines	±2 kV power supply lines ±1 kV input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge GB/T 17626.5	±1 kV line to line ±2 kV line to ground	±1 kV line to line ±2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips and interruption GB/T 17626.11	< 5 % $U_T$ , for 0.5 cycle (>95% dip at $U_T$ ) 40 % $U_T$ , for 5 cycles (>60% dip at $U_T$ ) 70 % $U_T$ , for 25 cycles (>30% dip at $U_T$ ) < 5 % $U_T$ , for 5 s (>95% dip at $U_T$ )	< 5 % $U_T$ , for 0.5 cycle (>95% dip at $U_T$ ) 40 % $U_T$ , for 5 cycles (>60% dip at $U_T$ ) 70 % $U_T$ , for 25 cycles (>30% dip at $U_T$ ) < 5 % $U_T$ , for 5 s (>95% dip at $U_T$ )	Mains power quality should be that of a typical commercial or hospital environment. It is recommended that an uninterruptible power supply or batteries be used if the user requires continuous operation of the IVUS System during

			a power outage.
Power frequency magnetic field (50 Hz) GB/T 17626.8	3 A/m	3 A/m, 50 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Note: U <sub>T</sub> refers to alternating current voltage before test voltage is applied.			


Guidelines and Manufacturer's Declaration - Electromagnetic Immunity			
The IVUS System (model No.: VH-60) is intended to be used in the following electromagnetic environment. The purchaser or user of the IVUS System (model No.: VH-60) should ensure that it is used in the following electromagnetic environment:			
Immunity test	IEC 60601 test level	Conformance level	Electromagnetic environment - guidelines
Conducted RF GB/T 17626.6 Radiated RF GB/T 17626.3	3 V (effective value) 150 kHz–80 MHz  3 V/m 80 MHz–2.5 GHz	3 V (effective value)  3 V/m	<p>Portable and mobile RF communications equipment should not be used near any component of the IVUS System (model No.: VH-60), including cables, at a distance less than the recommended distance. This distance should be calculated based on the formula for the corresponding emission frequency.</p> <p><b>Recommended separation distance</b></p> $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P} \quad 80 \text{ MHz} - 800 \text{ MHz}$ $d = 2.3\sqrt{P} \quad 800 \text{ MHz} - 2.5 \text{ GHz}$ <p>In the equation,  <i>P</i> —Maximum rated output power of the emitter provided by the emitter manufacturer, in watts (W);  <i>d</i>—recommended isolation distance, in meters (m)<sup>b</sup>.</p> <p>The field strength of the fixed radiofrequency emitter is determined through a survey (°) of the electromagnetic field. Every frequency range (°) should be lower than the conformance level.</p> <p>Interference may occur near devices with the following symbols.</p> 
<p>Note 1: The high-frequency band formula is used for frequencies of 80 MHz and 800 MHz.</p> <p>Note 2: These guidelines may not apply to all situations as electromagnetic transmission is affected by the absorption and reflection of buildings, objects, and humans.</p>			

Recommended isolation distance between portable and mobile radiofrequency communications equipment and the IVUS System (model No.: VH-60)			
The IVUS System (model No.: VH-60) is intended to be used in an electromagnetic environment in which radio-frequency radiation interference is controlled. Based on the maximum rated output power of the communications equipment, the purchaser or user can prevent electromagnetic interference based on the following recommended minimum distance between portable and mobile radio-frequency communications equipment (emitters) and the IVUS System (model No.: VH-60).			
Maximum rated output power of emitter (W)	Isolation distance from emitters with different frequencies/m		
	150 kHz–80 MHz $d = 1.2\sqrt{P}$	80 MHz–800MHz $d = 1.2\sqrt{P}$	800 MHz–2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00
<p>For emitters whose maximum rated output power are not listed in the above table, it is recommended that the separation distance is d and meters (m) as the unit. The formula in the corresponding emitter frequency field can be used to determine the isolation distance. In this equation, P is the maximum rated output power of the emitter provided by the manufacturer of the emitter, with watts (W) as the unit.</p> <p>Note 1: The high-frequency range formula is used for frequencies of 80 MHz and 800 MHz.</p> <p>Note 2: The above guidelines may not apply to all situations as electromagnetic transmission is affected by absorption and reflection of buildings, objects, and humans.</p>			

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