# **TRUEVISION Intravascular Ultrasound Imaging Catheter**

## **INSTRUCTION FOR USE (Indonesian) INSTRUKSI PENGGUNAAN**

*Insight Lifetech*

## Contents

<b>Instruction for Use of the TrueVision Intravascular Ultrasound Imaging Catheter .....</b>	<b>2</b>
<b>1. Product name.....</b>	<b>2</b>
<b>2. Catalogue number .....</b>	<b>2</b>
<b>3. Users .....</b>	<b>2</b>
<b>4. Diagram .....</b>	<b>3</b>
<b>5. Packaging contents.....</b>	<b>3</b>
<b>6. Product description.....</b>	<b>4</b>
<b>7. Intended purpose.....</b>	<b>4</b>
<b>8. Contraindications.....</b>	<b>5</b>
<b>9. Adverse reactions .....</b>	<b>5</b>
<b>10. Warnings.....</b>	<b>6</b>
<b>11. Precautions .....</b>	<b>7</b>
<b>12. Directions for Use .....</b>	<b>8</b>
<b>13. Product performance .....</b>	<b>11</b>
<b>14. Environmental requirements .....</b>	<b>12</b>
<b>15. Packaging and labels.....</b>	<b>12</b>
<b>16. Manufacturer information .....</b>	<b>14</b>
<b>17. Electromagnetic compatibility information .....</b>	<b>14</b>
<b>18. Acoustic Output--Conformance to GB9706.9-2008 .....</b>	<b>17</b>
<b>19. After-sales services .....</b>	<b>20</b>
<b>Instruction Manual for the Pullback Unit Sterile Bag.....</b>	<b>21</b>

## **Introduction of Instruction for Use**

This Instruction for Use contains usage instructions for the operator, including the use of the catheter and pullback unit sterile bag. This catheter must be simultaneously connected to the VivoHeart Intravascular Ultrasound Imaging System (VH-60) (hereafter referred to as the IVUS System) for operation. Please refer to the instruction manual of the IVUS System for details on the operation and safety information of the IVUS System.

### **Warning and attention are defined as follows:**

- 1) **Warning:** This may involve the personal safety of the patient or user. Ignoring this information may result in impairment to the personal safety of the patient or user.
- 2) **Attention:** This document provides information that may be particularly useful to you or help maintain or clarify important instructions.

# **Instruction for Use of the TrueVision Intravascular Ultrasound Imaging Catheter**

## **Warning**

- Please read this instruction for use carefully before using the catheter. Please comply with all contraindications, warnings, attention, and adverse events mentioned in the instruction manual. Otherwise, this may lead to complications in the patient. The physicians confirm and assess foreseeable surgery risks and inform the patient.
- The product provided has been sterilized using ethylene oxide. Please do not use it if the sterile barrier system is damaged. Call Insight Lifetech or contact your local sales representative if the product is found damaged.
- The TrueVision Intravascular Ultrasound Imaging Catheter (hereafter referred to as the IVUS Catheter) is a disposable medical device. Please do not reuse, rework, or re-sterilize it. Reuse, rework, or re-sterilization may damage the structural integrity of this product and/or cause device malfunction, thereby resulting in injury, disease, or even death to the patient. Reuse, rework, or re-sterilization may also cause the device to become contaminated and/or cause infection or cross-infection in patients. This includes but is not limited to the transmission of infectious disease from one patient to another patient. Device contamination may result in injury, disease, or even death to the patient.
- After usage, used products and packaging should be disposed of according to pertinent regulations of the hospital, competent authority, and/or local government.

## **1. Product name**

TrueVision Intravascular Ultrasound Imaging Catheter (IVUS Catheter)

## **2. Catalogue number**

02-6F40-01

02-6F60-01

## **3. Users**

This product should be performed in a specific manner by a physician skilled in interventional cardiology or interventional radiology and intravascular ultrasound examination techniques

#### 4. Diagram


Figure 1. TrueVision Intravascular Ultrasound Imaging Catheter

Table 1 Name description

No.	Name	No.	Name
1	Distal assembly	9	Ultrasonic transducer
2	Imaging window shaft	10	Rotating drive cable
3	Proximal shaft	11	Connector assembly
4	Distal stress relief tube	12	Distal tip
5	Telescoping shaft assembly	13	Guidewire exit port
6	Proximal stress relief tube	14	Distal fluid exit port
7	One-way stopcock	15	Transducer housing
8	Hub	16	Marker band

#### 5. Packaging contents

- One TrueVision Intravascular Ultrasound Imaging Catheter (IVUS Catheter)
- One pullback unit sterile bag
- One extension tube
- One 3 ml syringe and one 10 ml syringe respectively
- One three-way stopcock
- One IFU
- One certificate of qualification

## **6. Product description**

The TrueVision Intravascular Ultrasound Imaging Catheter consists of the catheter and the accessories, of which the catheter consists of an outer sheath and an imaging core assembly. The accessories consist of the extension tube, syringe (3ml and 10ml), and three-way stopcock. The catheter and accessories are single-use and EO sterilized, and the shelf life is two years.

The proximal end of the catheter has a connection port that is linked to the pullback unit so that the catheter and pullback unit can be connected. The pullback unit catheter connection port is composed of a mechanical drive socket and motor adapter. In addition, the proximal connector is embedded with a radiofrequency identification (RFID) tag to facilitate matching with the pullback unit.

The distal shaft lumen and the proximal shaft lumen constitute the effective working length of the catheter and the telescopic assembly partially extends outside the guiding catheter. A marker is embedded in the catheter at 0.5cm from the distal end. A guidewire exit port is present in the catheter at 1.6cm from the distal tip. A 28-cm long imaging window shaft is present at the distal end of the catheter, which is connected to the proximal shaft. In addition, there are two markers on the outer wall of the proximal shaft at 90 cm and 100 cm from the distal end, which are used to help estimate the position of the catheter relative to the distal end of the guide catheter. The telescopic assembly enables the imaging core assembly to move 15 cm linearly forwards and backwards while the transducer can move accordingly within the imaging window shaft. The outer wall of the inner telescopic tube contains a series of proximal markers at 1cm intervals, which can be used to assess the length of the lesion. The one-way stopcock contains a flush port that can be used to flush the interior of the catheter to maintain its cleanliness and expel air from the catheter as much as possible. Before use, the catheter must be flushed with heparinized saline. The purpose of this is to provide an acoustic coupling medium required for ultrasound imaging. During use, the one-way stopcock also helps retain saline in the catheter.

## **7. Intended purpose**

The IVUS Catheter is intended to be used in medical institutions for the ultrasonic examinations of intravascular lesions in coronary arteries. Intravascular ultrasound imaging is indicated in patients who are candidates for transluminal coronary interventional procedures. The IVUS Catheter is intended to be used in conjunction with the VivoHeart Intravascular Ultrasound Imaging System VH-60 (hereinafter called the IVUS System).

## 8. Contraindications

This IVUS Catheter should not be used if the use of any catheter constitutes a threat to patient safety.

The contraindications generally include the following symptoms:

- Bacteremia or sepsis
- Severe coagulation disorder
- Contraindications for coronary intervention surgery
- Severe hemodynamic instability or shock
- Patients diagnosed with severe coronary artery spasms
- Life-threatening arrhythmia

## 9. Adverse reactions

Vascular imaging has all the risks and discomfort related to catheterization. These risks or discomfort may occur at any time, with various frequencies and severity. In addition, these complications may result in the need for additional treatments such as surgical intervention, and may lead to death in rare cases. The possible adverse reactions are as follows:

- Anaphylaxis
- Angina
- Cardiac arrest
- Arrhythmia (including but not limited to ventricular tachycardia, atrial/ventricular fibrillation, and complete atrioventricular block)
- Cardiac tamponade or pericardial effusion
- Death
- Device jam (requiring surgical intervention)
- Embolism (air, foreign body, tissue, or thrombus)
- Bleeding or hematoma
- Hypotension
- Infection
- Myocardial infarction
- Myocardial ischemia
- Stroke and transient cerebral ischemia
- Thrombosis
- Vascular occlusion and sudden closure
- Vascular injury (including but not limited to dissection and perforation)

## 10. Warnings

- This device should not be used after the expiry date on the label. The use of expired devices will lead to patient injury due to device degradation.
- Intravascular ultrasound imaging examinations should only be performed in a well-equipped cardiovascular intervention surgical room by physicians who are skilled in interventional cardiology, interventional radiology, or intravascular ultrasound examinations.
- The IVUS Catheter does not contain any component that the user can service. Never try to repair or modify any component of the catheter assembly. The use of modified catheters will result in poor imaging results or lead to complications in the patient.
- Please do not refit this device.
- Air infiltration of the IVUS Catheter and washing accessories may cause patient injury or death. Before inserting the IVUS Catheter into the IVUS System, it is important to confirm that the air in the catheter and washing accessories has been exhausted.
- Please do not squeeze, press, twist, or violently bend the catheter. Otherwise, the catheter performance may decrease, leading to vascular damage or complications in the patient.
- Please do not use coarse surface gauze to repeatedly rub the distal end of the catheter which has a hydrophilic coating to avoid disrupting the integrity of the coating. Otherwise, the hydrophilic coating will peel off and enter the blood vessels. This may lead to thrombosis or other complications in the patient.
- Please do not push or pull back the IVUS Catheter in the absence of an X-ray fluoroscopy. This may lead to vascular damage or complications in the patient.
- Please do not push the catheter if resistance is encountered. Please do not forcefully insert the catheter in a lumen that is smaller than the catheter or through a stenosed area. Forcefully pushing the catheter may lead to catheter damage, thereby leading to vascular damage or complications in the patient.
- If the catheter does not completely cover the guidewire when it is inserted into a blood vessel containing a stent, this may cause the stent to be stuck between the joint of the catheter and guidewire, resulting in catheter/guidewire jam, separation of the catheter tip, and/or stent migration.
- If resistance is encountered during catheter withdrawal, X-ray fluoroscopy should be used for examination, then the entire catheter and guidewire should be withdrawn together. Forceful withdrawal of the catheter may lead to vascular damage or complications in the patient.
- After the stent is deployed, the catheter must not be pushed forward when re-introduced to cross over the guidewire that may pass between one or more stent trabeculae. When it re-passes through the stent, the guidewire may pass through one or more stent struts. Subsequently,

pushing the catheter may cause the catheter and stent to be intertwined, leading to catheter/guidewire jam, separation of the catheter tip, and/or stent migration. Therefore, care should be taken when removing the catheter from a blood vessel where a stent was implanted.

- Stent malapposition or overlap and/or bending at the end of the blood vessel containing an implanted stent may lead to catheter and stent jamming. The separation between the guidewire and catheter or guidewire bending when the catheter is withdrawn may lead to guidewire kinking, causing damage to the distal tip and/or vascular damage. Guidewire bending or tip damage may lead to the catheter being stuck at the stent strut, leading to jamming.
- If the catheter needs to be inserted multiple times, the connection between the catheter and the pullback unit should not be severed to prevent disruption of catheter sterility.

## **11. Precautions**

- Please do not connect the catheter to electronic equipment other than the stipulated device as it may cause the catheter to be unable to work normally.
- Please do not attempt to connect or disconnect the catheter when the pullback unit is operating as this may damage the connector.
- If the guidewire cannot be successfully connected to the distal end of the catheter, please inspect the guidewire outlet for damage before inserting the catheter into the blood vessel. Guidewire outlet damage may increase resistance when the catheter is pushed forward or withdrawn.
- Please do not push the catheter forward in the absence of the guidewire as this may result in difficulty reaching the intended target region or cause the distal tip of the catheter to bend.
- Please do not push the catheter tip near the extremely soft end of the guidewire. This part of the guidewire cannot fully support the catheter. When withdrawing the catheter, the catheter that is pushed to this position may not follow the guidewire, causing the guidewire to bend into a ring. The catheter may drag the guidewire ring along the blood vessel and jam the tip of the guide catheter. If this situation occurs, the catheter assembly, guidewire, and guide catheter should be removed as a whole. If the catheter is pushed too close to the end of the guidewire, please ensure that the IVUS Imaging Catheter is stable when pushing the guidewire. If this operation does not work, withdraw the catheter and guidewire together.
- Please do not push or withdraw the IVUS Imaging Catheter when the imaging core assembly is not located at the furthest end of the imaging window lumen as this may cause the catheter to kink.
- Please inspect the catheter carefully during and after surgery. The catheter should be inspected for any possible damage during use. Multiple insertions may cause deformation of/changes to

the dimensions of the catheter outlet, thereby increasing the possibility of entanglement between the catheter and the stent. Please be extra careful when reinserting and/or withdrawing the catheter to avoid damaging the outlet.

- The pullback unit must be turned off before pushing and withdrawing the IVUS Catheter to prevent overloading of the pullback unit.

## **12. Directions for Use**

### **A. Inspection before use**

Before usage, please inspect the sterile barrier of the packaging for any damage and inspect the catheter and accessories for any damage or defects. Do not use any device that may be contaminated or defective. If the integrity of the sterile barrier system is disrupted or its contents are damaged, please contact an Insight Lifetech sales representative. Before imaging, all equipment used for the intervention should be carefully checked to ensure normal performance. Please contact an Insight Lifetech sales representative if the device is damaged.

Please ensure that the product has not expired before use. Please do not use products that have exceeded the expiration date. Please carry out processing and packaging according to the policies of the hospital, competent authority, and/or local government.

### **B. Preparation for catheter use**

- 1) For usage of the IVUS System and the pullback unit, read the section “Use of the IVUS System” in the IFU for the IVUS System.

***Note: The following steps must be carried out using an aseptic technique.***

- 2) Remove the catheter and accessories from the tray using an aseptic technique. Fully retract the active imaging core assembly to the proximal position through the telescoping shaft. Do not exert excessive force in retracting the imaging core assembly.
- 3) Connect a 3 ml (cc) and 10 ml (cc) syringe to the three-way stopcock before connecting the assembly to the extension tube. Fill these two syringes with heparinized saline. Connect the extension tube to the one-way stopcock on the connector. Use the 10 ml (cc) syringe as a reservoir to refill the 3 ml (cc) syringe.
- 4) Continuously flush the IVUS Catheter three times using 3 ml (cc) saline each time on the preparation table. Please do not use excessive pressure. Move the IVUS Catheter to the operating table. Ensure that all air in the catheter and accessories is expelled.
- 5) Connect the catheter to the pullback unit by aligning the proximal hub and pullback unit. Push the connector into the pullback unit until the hub clicks into place. Gently pull the hub to ensure that it is completely in the pullback unit.

- 6) Remove the catheter carefully from the carrier tubing. Verify that the imaging core assembly is in the fully retracted position and that the catheter is not tightly intertwined with the assembly. Open the pullback unit and observe the bright concentric ring images on the display to confirm that the catheter is working normally. Please contact an Insight Lifetech sales representative if the catheter is not recognized or wrongly identified.
- 7) Use the telescoping shaft to completely push the imaging core assembly to the distal end when using the pullback unit for imaging.

***Note: Please turn on the pullback unit before using the telescoping shaft to push the imaging core assembly forward into the catheter.***

- 8) Turn off the pullback unit. At this point, the pullback unit should remain turned off until the catheter is in the correct position.
- 9) According to requirements, fill the 10 ml (cc) syringe and reconnect to the 3-way stopcock, taking care not to introduce air into the catheter.
- 10) To prevent air from being introduced into the catheter, please do not retract the imaging core assembly before the catheter is inserted. Additional flushing will be required if the imaging core assembly is retracted before the catheter is inserted.

***Note:***

- a) ***If the pullback unit needs to be used, flush the catheter again when the imaging core assembly is located at the most distal end and the catheter is installed in the pullback unit.***
- b) ***If flushing cannot be carried out when the imaging core assembly is located at the most distal end, manually retract the imaging core assembly by 3-5mm and repeat flushing. Following that, manually push the imaging core assembly to the most distal end.***
- c) ***Pay extra attention during the operation to prevent the catheter telescopic tube from bending.***

### **C. Preparation for the pullback unit sterile bag use**

Please refer to the section “Instructions for use of the pullback unit sterile bag”.

### **D. Place guide catheter**

- 1) Use the catheter sheath guide to prepare the entry site according to standard procedure.
- 2) Ensure that the patient is prepared according to standard imaging procedure for intervention treatment before inserting the IVUS Imaging Catheter.
- 3) Insert the guide catheter and Y-connector. Insert the guidewire and push it to the target region.

### **E. Insertion of IVUS catheter into the guide catheter**

- 1) Use heparinized saline to moisten the distal end (around 28 cm) of the IVUS Catheter

sheath to activate the hydrophilic coating. Before placing the catheter on the guidewire, please wipe the guidewire with heparinized saline.

**Note: Please do not use medical gauze to repeatedly wipe the distal end of the IVUS Catheter.**

- 2) Replace the guidewire backwards into the distal end of the IVUS Catheter. Pass the guidewire through the IVUS Catheter until the guidewire passes through the guidewire outlet.

**Note: It is recommended to use a guidewire with high hardness at the distal tip.**

- 3) Continue pushing the catheter forward to the exit point of the guide catheter. Use the proximal marker as a reference if required. Tighten the hemostatic valve at the Y-connector of the guide catheter until it is tight enough to prevent liquid/blood leakage. As the hemostatic valve and the rotating drive cable are bound together, excessive tightening may distort the image.
- 4) Open the pullback unit to inspect if the catheter generates images. If the image flickers, this means there may still be air in the catheter. Flush the catheter again when imaging is started on the pullback unit. Please do not use excessive pressure. The image should show a single bright concentric circle. After confirming that the image is stable, stop imaging by pressing the Imaging button on the pullback unit.

## **F. Catheter placement and imaging**

- 1) Turn off the pullback unit and push the IVUS Catheter forward along the guidewire under X-ray fluoroscopy until the distal marker is at least 3 cm past the blood vessel/lesion.
- 2) Fix the catheter body and guidewire, turn on the pullback unit, and slowly retract the imaging core assembly along the length manually or through the selected automatic pullback speed (maximum retraction distance: 15 cm) to image the target region. Retract or push forward according to requirements.

**Note:**

- a) **The pullback unit must be turned on before pushing or retracting the imaging core assembly in the catheter.**
- b) **If the image is weakened please refer to sections G(3) and G(4).**
- 3) After imaging is completed, push the imaging core assembly to the most distal end and turn off the pullback unit. Maintain the position of the guidewire and remove the IVUS Catheter.
- 4) Please use a 3 ml (cc) syringe to flush and coil the catheter and place the pullback unit at one side before reinserting the IVUS Catheter.

**Note: The connection between the catheter and the pullback unit should not be interrupted if multiple insertions are performed to prevent disrupting catheter sterility.**

- 5) When ready to insert the IVUS Catheter, flush one more time with a 3 mL (cc) syringe.

- 6) Before reinserting the catheter, please inspect the guidewire outlet to ensure that no damage occurred during pullback.
- 7) After use, dispose of the product and packaging following the patient policy of hospital, competent authority, and/or local government.

## G Troubleshooting

- 1) Please stop imaging if vibrations occur in the telescoping part of the catheter when the imaging core assembly is pushed forward. Adjust the position of the outer sheath and ensure that the imaging window shaft is straight. Restart the pullback unit again and push the imaging core assembly forward again. Please pay attention if you want to continue using the catheter as the imaging core assembly may become entangled.
- 2) If the IVUS Catheter cannot be recognized by the IVUS System, contact an Insight Lifetech sales representative before operation.
- 3) Air bubbles may be present at the distal imaging window if the image is weakened during use or shadows still exist in the image after flushing the assembly at the original position. Remove the catheter and repeat the flushing procedure in B section "Preparatory work for catheter usage".
- 4) If the image cannot be restored after flushing, the imaging core assembly may have failed or the IVUS Catheter may have disconnected from the pullback unit. Stop imaging and check whether the hub is located inside the pullback unit. Please withdraw the catheter if the hub is inside the pullback unit but the malfunction is still present. Restart the pullback unit and visually inspect if the imaging core assembly is running. If the assembly is not running, please send the catheter back to an Insight Lifetech sales representative for analysis.

## 13. Product performance

No.	Performance Indicator	40 MHz System	60 MHz System
1	Acoustic working frequency	40 MHz, error < ±15%	60 MHz, error < ±15%
2	Imaging radius	≥ 6 mm	≥ 4 mm
3	Axial resolution	≤ 50 μm (depth: 3 mm) ≤ 50 μm (depth: 1.5 mm)	≤ 40 μm (depth: 2 mm) ≤ 40 μm (depth: 1.5 mm)
4	Lateral resolution	≤ 600 μm (depth: 3 mm) ≤ 300 μm (depth: 1.5 mm)	≤ 300 μm (depth: 2 mm) ≤ 250 μm (depth: 1.5 mm)
5	Pullback direction resolution	≤ 600 μm (depth: 3 mm)	≤ 500 μm (depth: 2 mm)
6	Longitudinal geometric position accuracy	≤ 10% (within the imaging radius range)	
7	Transverse		

	geometric position accuracy	$\leq 10\%$ (within the imaging radius range)
8	Pullback direction geometric position accuracy	$\leq 10\%$
9	Image geometric distortion	$\leq 10\%$ (within the imaging radius range)

#### 14. Environmental requirements

<b>Storage Environment Temperature</b>	Room temperature
<b>Storage Environment Humidity</b>	Uncontrolled
<b>Storage Environment Pressure</b>	Uncontrolled
<b>Operating Environment Temperature</b>	10 °C to 35 °C
<b>Operating Environment Humidity</b>	30% RH-75% RH (non-condensing)
<b>Operating Environment Pressure</b>	525 mmHg to 795 mmHg (700 hPa to 1060 hPa)
<b>Transportation Environment Temperature</b>	-29 °C to + 60 °C
<b>Transportation Environment Humidity</b>	30 % RH to 85 % RH
<b>Transportation Environment Pressure</b>	Uncontrolled

#### 15. Packaging and labels

Hospitals or surgeons should only use single-use intravascular ultrasound diagnostic catheters with factory packaging and complete identification. Please do not use the IVUS Catheter if the packaging is opened or damaged or the label is not complete or difficult to recognize. Symbols in the following table will appear on the packaging or label:

No.	Symbol	Explanation	No.	Symbol	Explanation
1		Packaging quantity: 1	10		Date of manufacture
2		Catalogue number	11		Use by date
3		Batch code	12		Do not use if package is damaged
4		Warning	13		Keep dry
5		Follow Instructions for Use	14		Keep away from light

6		Do not reuse	15		Defibrillation-proof type CF applied part
7		Do not resterilize	16		This side up
8		Sterilized using ethylene oxide	17		Fragile, handle with care
9		Non-pyrogenic	18	IPX7	In normal use, the degree of protection for the intake of fluid into the distal assembly and proximal shaft into the body part is IPX7

## 16. Manufacturer information

<b>Name of Registrant/Manufacturer</b>	Insight Lifetech Co., Ltd.
<b>Address of Registrant/Manufacturer</b>	2102, Tower A, Hongrongyuan North Station Center, North Station Community, Minzhi, Longhua District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
<b>Manufacturing Site</b>	Zone A, F, 2 <sup>nd</sup> Floor, Zone C, E, F, 3 <sup>rd</sup> Floor, Zone A, 4 <sup>th</sup> Floor, Zone A, B, 10 <sup>th</sup> Floor, Building 3, Tinwe Business Park, No.6 Liufang Road, Baoan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA 402, Building 1, Shenzhen Biomedicine Innovations Industrial Park, No. 14, Jinhui Road., Jinsha Community, Kengzi, Pingshan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
<b>Number of Product Registration Certificates</b>	国械注准 20223060974
<b>Number of Product Technical Requirements</b>	国械注准 20223060974
<b>Number of Manufacturing License</b>	粤食药监械生产许 20204044 号
<b>Version of IFU</b>	01
<b>Release Date of IFU</b>	2023-06-26
<b>Product Shelf Life</b>	2 years

## 17. Electromagnetic compatibility information

Special preventive measures for electromagnetic compatibility (EMC) are required for electrical devices for medical use. The IVUS Catheter should be installed and used according to the EMC information in the instruction manual of the IVUS System that it is used in conjunction with. Please refer to the requirements, guidelines, and manufacturer's declaration – electromagnetic radiation, as well as guidelines and manufacturer's declaration - electromagnetic immunity in the instruction manual for the IVUS System.

### Warnings

- The IVUS Catheter connected to the IVUS System shall not be used close to or stacked on other equipment. If it must be used close to or stacked on other equipment, observation shall be carried out to verify that it can operate normally under the used configurations.
- Except for the cables sold by the manufacturer of the IVUS Catheter as spare internal parts, the use of accessories and cables other than those specified may result in increased emission or decreased immunity when the IVUS Catheter is connected to the IVUS System.
- Combined use of the IVUS Catheter and accessories, transducer, or cable outside of specifications may result in increased emission of the IVUS System or decreased immunity.

**Attention**

- The user shall install and use the equipment according to the electromagnetic compatibility information provided in the accompanying documents;
- Portable and mobile RF communication equipment may affect the performance of the IVUS Catheter connected to the IVUS System during use, avoid strong electromagnetic interference, such as that generated close to mobile phones and microwave ovens, etc. when using it.

**Basic performance**

- When real-time imaging acquisition is used in the IVUS Catheter connected to the IVUS System, the IVUS Catheter will continuously output sound under the default acoustic working frequency with maximum gain and maximum imaging depth.

Guidelines and Manufacturer's Declaration - Electromagnetic Emission			
When the IVUS Catheter is connected to the IVUS System for use in the following electromagnetic environment as intended, the purchaser or user of the IVUS Catheter should ensure that it is used in this electromagnetic environment:			
Emission test	Conformance	Electromagnetic environment - guidelines	
Radiofrequency emission GB 4824	1 set	Radiofrequency energy is used for the IVUS Catheter connected to the IVUS System only for its internal function. Therefore, its radiofrequency emission is low and has a low possibility of interfering with nearby electronic equipment.	
Radiofrequency emission GB 4824	Class A	The IVUS Catheter connected to the IVUS System is suitable for non-domestic use and all facilities that are not directly connected to public low-voltage power supply networks for residences.	
Harmonic emission GB 17625.1	Not applicable		
Voltage fluctuation/flicker emission GB 17625.2	Not applicable		
Guidelines and Manufacturer's Declaration - Electromagnetic Immunity			
When the IVUS Catheter is connected to the IVUS System for use in the following electromagnetic environment as intended, the purchaser or user of the IVUS Catheter should ensure that it is used in this electromagnetic environment:			
Immunity Test	IEC 60601 test level	Conformance level	Electromagnetic environment - guidelines
Electrostatic discharge GB/T 17626.2	±6 kV contact discharge ±8 kV air discharge	±6 kV contact discharge ±8 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast	±2 kV for the power cord	±2 kV power cord	Mains power quality should be

transient/ burst GB/T 17626.4	±1 kV for input/output cable	±1 kV input/output cable	that of a typical commercial or hospital environment.
Surge GB/T 17626.5	±1 kV line to line ±2 kV line to ground	±1 kV line to line ±2 kV line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips and Interruptions GB/T 17626.11	< 5 % $U_T$ , for 0.5 cycle (>95% sag at $U_T$ ) 40 % $U_T$ , for 5 cycles (>60% sag at $U_T$ ) 70 % $U_T$ , for 25 cycles (>30% sag at $U_T$ ) < 5 % $U_T$ , for 5 s (>95% sag at $U_T$ )	< 5 % $U_T$ , for 0.5 cycle (>95% sag at $U_T$ ) 40 % $U_T$ , for 5 cycles (>60% sag at $U_T$ ) 70 % $U_T$ , for 25 cycles (>30% sag at $U_T$ ) < 5 % $U_T$ , for 5 s (>95% sag at $U_T$ )	Mains power quality should be that of a typical commercial or hospital environment. It is recommended that an uninterruptible power supply or batteries be used if the user requires continuous operation of the IVUS System during a power outage.
Power frequency magnetic field (50 Hz) GB/T 17626.8	3 A/m	3 A/m, 50 Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

*Note:  $U_T$  refers to alternating current voltage before test voltage is applied*

### **Guidelines and Manufacturer's Declaration - Electromagnetic Immunity**

When the IVUS Catheter is connected to the IVUS System for use in the following electromagnetic environment as intended, the purchaser or user of the IVUS Catheter should ensure that it is used in this electromagnetic environment:

<b>Immunity Test</b>	<b>IEC 60601 Test level</b>	<b>Conformance level</b>	<b>Electromagnetic environment - guidelines</b>
Conducted RF GB/T 17626.6	3 V (effective value) 150 kHz–80 MHz	3 V (effective value)	Portable and mobile RF communications equipment should not be used near any component of the VivoHeart Intravascular Ultrasound Imaging System, including cables, at a distance less than the recommended distance. This distance should be calculated based on the formula for the corresponding emission frequency. <b>Recommended separation distance:</b> $d = 1.2\sqrt{P}$ , 150 kHz–80 MHz; $d = 1.2\sqrt{P}$ , 80 MHz–800 MHz; $d = 2.3\sqrt{P}$ , 800 MHz–2.5 GHz; In the equation, $P$ -Maximum rated output power of the emitter provided by the emitter manufacturer, in watts (W); $d$ -recommended separation distance, in meters (m). The field strength of the fixed radiofrequency emitter is determined through a survey of the electromagnetic field. Every frequency range should be lower than the conformance level. Interference may occur near devices with the
Radiated RF GB/T 17626.3	3 V/m 80 MHz–2.5 GHz	3 V/m	

			following symbols. 
<p>Note 1: The high-frequency band formula is used for frequencies of 80 MHz and 800 MHz.</p> <p>Note 2: The guidelines are not suitable for all situations as electromagnetic transmission is affected by the absorption and reflection of buildings, objects, and humans.</p>			
<p><b>Recommended separation distance between portable and mobile radiofrequency communications equipment and the TrueVision Intravascular Ultrasound Imaging Catheter</b></p>			
<p>The IVUS Catheter is expected to be used in an electromagnetic environment in which radiofrequency radiation interference is controlled. Based on the maximum rated output power of the communications equipment, the purchaser or user can prevent electromagnetic disturbance based on the recommended minimum distances between portable and mobile radiofrequency communications equipment (transmitter) as follows and the IVUS Catheter.</p>			
Maximum rated output power of emitter (W)	The separation distance from emitters with different frequencies/m		
	150 kHz–80 MHz	80 MHz–800MHz	800 MHz–2.5 GHz
	$d = 1.2\sqrt{P}$	$d = 1.2\sqrt{P}$	$d = 2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.20	1.20	2.30
10	3.80	3.80	7.30
100	12.00	12.00	23.00
<p>Description: For emitters whose maximum rated output power is not listed in the above table, it is recommended that the formula in the corresponding emitter frequency field can be used to determine the separation distance, with recommended separation distance being d, and meters (m) as the unit. In this equation, P is the maximum rated output power of the emitter provided by the manufacturer of the emitter, with watts (W) as the unit.</p>			
<p>Note 1: The high-frequency range formula is used for frequencies of 80 MHz and 800 MHz.</p> <p>Note 2: The above guidelines may not apply to all situations as electromagnetic transmission is affected by the absorption and reflection of buildings, objects, and humans.</p>			

## 18. Acoustic Output--Conformance to GB9706.9-2008

- Acoustic output working mode: B automatic scanning mode

Application: Coronary angiography

Note: Intravascular ultrasound is classified under the document “Information for Manufacturers Seeking Marketing Clearance of Diagnostic Ultrasound Systems and Transducers” guided by the FDA to determine the maximum allowable acoustic output energy. And the catheter is not suitable for fetal imaging.

- Statement of cautious**

The operator of the VivoHeart Intravascular Ultrasound Imaging System is responsible for

understanding the risks of acoustic output generated by the IVUS Imaging Catheter and its related IVUS Imaging Catheter and adopting appropriate measures to decrease such risks.

- **Information on the IVUS Catheter**

IVUS Imaging Catheter with different model numbers has different acoustic outputs. Each Insight Lifetech IVUS Imaging Catheter contains an “Instruction Manual”, including declaration and forms on acoustic output. International test and terminology definition requirements can be found in GB9706.9-2008 – *Particular requirements for the basic safety and essential performance of ultrasonic medical diagnostic and monitoring equipment*.

40 MHz system acoustic output table

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum index value		0.11	N/A	0.032	N/A	N/A	N/A
Index component value			0.032		N/A	N/A	N/A
Acoustic Parameters	$p_{ra}$ at $Z_{MI}$ (MPa)	0.68					
	P (mW)		0.16				N/A
	$P_{1 \times 1}$ (mW)		0.16				
	$Z_s$ (cm)			N/A			
	$Z_{bp}$ (cm)						
	$Z_b$ (cm)					N/A	
	$Z_{MI}$ (cm)	0.1					
	$Z_{pii, \alpha}$ (cm)	0.1					
	$f_{awf}$ (MHz)	39.35	N/A	39.35	N/A	N/A	N/A
Other Information	$p_{rr}$ (Hz)	29790					
	$S_{rr}$ (Hz)	29.93					
	$n_{pss}$	1					
	$I_{pi, \alpha}$ at $Z_{pii, \alpha}$ (W/cm <sup>2</sup> )	10.95				N/A	
	$I_{spta, \alpha}$ at $Z_{pii, \alpha}$ or $Z_{sii, \alpha}$ (mW/cm <sup>2</sup> )	19.99					
	$I_{spta}$ at $Z_{pii}$ or $Z_{sii}$ (mW/cm <sup>2</sup> )	27.63	N/A	Fixed	N/A	N/A	N/A
	$p_r$ at $Z_{pii}$ (MPa)	0.79	N/A	Fixed	N/A	N/A	N/A
Operating control conditions	Focus (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
	Depth (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
	Frequency (MHz)	40	N/A	40	N/A	N/A	N/A
Note: N/A means that there is no corresponding intended use for the probe or there is no data to be							

filled in.  
Uncertainty: sound intensity:  $\pm 25.12\%$ , sound pressure:  $\pm 12.56\%$ , sound power:  $\pm 25.12\%$ , central frequency:  $\pm 0.13\%$

60 MHz system acoustic output table

Index label		MI	TIS		TIB		TIC
			At surface	Below surface	At surface	Below surface	
Maximum index value		0.18	N/A	0.011	N/A	N/A	N/A
Index component value			0.011		N/A	N/A	N/A
Acoustic Parameters	$p_{ra}$ at $z_{MI}$ (MPa)	1.17					
	P (mW)		0.055				N/A
	$P_{1 \times 1}$ (mW)		0.055				
	$z_s$ (cm)			N/A			
	$z_{bp}$ (cm)						
	$z_b$ (cm)					N/A	
	$z_{MI}$ (cm)	0.1					
	$z_{pii, \alpha}$ (cm)	0.1					
Other Information	$f_{awf}$ (MHz)	59.18	N/A	59.18	N/A	N/A	N/A
	$p_{rr}$ (Hz)	29790					
	$S_{rr}$ (Hz)	29.23					
	$n_{pss}$	1					
	$I_{pi, \alpha}$ at $z_{pii, \alpha}$ (W/cm <sup>2</sup> )	12.73				N/A	
	$I_{spta, \alpha}$ at $z_{pii, \alpha}$ or $z_{sii, \alpha}$ (mW/cm <sup>2</sup> )	12.57					
	$I_{spta}$ at $z_{pii}$ or $z_{sii}$ (mW/cm <sup>2</sup> )	17.47	N/A	Fixed	N/A	N/A	N/A
Operating control conditions	$p_r$ at $z_{pii}$ (MPa)	0.88	N/A	Fixed	N/A	N/A	N/A
	Focus (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
	Depth (mm)	Fixed	N/A	Fixed	N/A	N/A	N/A
Operating control conditions	Frequency (MHz)	60	N/A	60	N/A	N/A	N/A

Note: N/A means that there is no corresponding intended use for the probe or there is no data to be filled in.

Uncertainty: sound intensity:  $\pm 25.12\%$ , sound pressure:  $\pm 12.56\%$ , sound power:  $\pm 25.12\%$ , central frequency:  $\pm 0.13\%$

Terms	Definition	Unit
-------	------------	------

MI	Mechanical index	n/a
TIS	Soft tissue thermal index	n/a
$p_{ra}$	Peak negative pressure after attenuation	MPa
$f_{awf}$	Acoustic working frequency	MHz
$P$	Output power	Mw
$z$	Point at specific length from the acoustic source on the acoustic axis	cm
$t_d$	Pulse duration	us
$prr$	Pulse repetition frequency	Hz
$p_r$	Peak negative pressure	MPa
$I_{pa,\alpha}$	Mean sound intensity of pulse after attenuation	W/cm <sup>2</sup>
$I_{pi}$	Pulse sound intensity integral	J/m <sup>2</sup>
$I_{pi,\alpha}$	Pulse sound intensity integral after attenuation	J/m <sup>2</sup>

## 19. After-sales services

Contact Insight Lifetech or local distributor immediately if you have any questions during the use of this device. Insight Lifetech does not assume any liability for devices that are reused, reprocessed, or repeatedly sterilized.

**After-sales service provider (manufacturer):** Insight Lifetech Co., Ltd.

**Contact address:** Zone E, 3rd Floor, Building 3, Tinwe Business Park, No.6 Liufang Road, Baoan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

**Tel:** +86 (755) 23074845

**Email:** service@insight-med.com

**Website:** <http://www.insight-med.com>

## **Instruction Manual for the Pullback Unit Sterile Bag**

### **Warnings**

- Please read this instruction manual carefully before using the bag. Please comply with all contraindications, warnings, attention, and adverse events mentioned in the instruction manual. Otherwise, this may lead to complications in the patient. Insight Lifetech relies on physicians to confirm and assess all foreseeable surgery risks and informs each patient.
- The product provided has been sterilized using radiation. Please do not use it if the sterile barrier system is damaged. Please contact your sales representative if the product is damaged.
- This product is for single use only. Please do not reuse, rework, or resterilize. Reuse, rework, or resterilization may damage the structural integrity of this bag and/or cause device malfunction, thereby resulting in injury, disease, or even death to the patient. Reuse, rework, or resterilization may also cause the device to become contaminated and/or cause infection or cross-infection in patients. This includes but is not limited to the transmission of infectious disease from one patient to another patient. Device contamination may result in injury, disease, or even death to the patient.
- After usage, used products and packaging should be disposed of according to pertinent regulations of the hospital, competent authority, and/or local government.

### **1. Device description**

This product can provide effective and conformed coverage that is suitable for the pullback unit.

### **2. Packaging content**

Pullback unit sterile bag.

### **3. Intended Use/Indications**

The pullback unit sterile bag is used to cover the pullback unit during an intravascular ultrasound to maintain the sterility of the pullback unit and prevent microorganisms, liquids, and particulate matter from being transferred to patients and healthcare workers.

### **4. Contraindications**

Unknown.

### **5. Adverse events**

None.

## **6. Precautions**

None.

## **7. Packaging method**

Please do not use the bag if the packaging is opened or damaged. Please do not use the bag if the label is not complete or difficult to recognize.

## **8. Usage and storage**

Please place it in a cool and dry place away from light.

## **9. Instructions**

### **A. Inspection before use**

Before usage, please check the sterile barrier system of the packaging for any damage and check the bag for any damage or defect. Please do not use any potentially contaminated or defective device. If the integrity of the sterile barrier system is disrupted or its contents are damaged, please contact an Insight Lifetech sales representative.

### **B. Covering the pullback unit**

Two people are required to cover the pullback unit: one in the sterile zone and one outside the sterile zone. An appropriate aseptic technique should be used.

- 1) The non-sterile operator opens a sealed bag containing the bag.
- 2) The sterile operator then takes the bag out of the sealed bag.
- 3) The sterile operator opens the bag and makes the "Insert here" label located on one side while the panel is located on the opposite side.
- 4) The sterile operator puts his/her hands into the innermost layer of the bag and opens the bag so that the non-sterile operator can place the pullback unit. The opening of the bag can be easily located based on the "Insert here" label.
- 5) Then the non-sterile operator places the pullback unit into the bag opening and locates the pullback unit suitably so that the top of the device faces the "This side up" label.
- 6) Following that, the sterile operator grabs the covered pullback unit. The non-sterile operator grabs and pulls the tab until the bag is fully expanded.
- 7) After alignment, immediately stop pushing the pullback unit into the bag and fix the bag at a suitable position. There is no need for the sterile operator to push the pullback unit until the panel of the bag. It is normal for a space to be present between the front end of the pullback unit and the panel of the bag.

- 8) To fix the bag on the pullback unit, the sterile operator must align the panel of the bag with the front end of the pullback unit and gently push it to secure the connection. Bag material cannot be stuck between the panel and the device. It is normal to have excess bag material around and outside the connection.
- 9) Prepare a catheter that is compatible with the pullback unit for insertion, tear down the label at the panel, and dispose of it properly.
- 10) The pullback unit that is inserted into the bag is now prepared and ready for use.

Product is under review by Ministry of Health Indonesia

# **Kateter Pencitraan Ultrasonografi Intravaskuler TRUEVISION**

## **INSTRUKSI PENGGUNAAN**

Product is under review by Ministry of Health Indonesia

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## 16. Informasi Produsen

<b>Nama Pendaftar/Produsen</b>	Insight Lifetech Co., Ltd.
<b>Alamat Pendaftar/Produsen</b>	2102, Tower A, Hongrongyuan North Station Center, North Station Community, Minzhi, Longhua District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
<b>Lokasi Produksi</b>	Zone A, F, 2 <sup>nd</sup> Floor, Zone C, E, F, 3 <sup>rd</sup> Floor, Zone A, 4 <sup>th</sup> Floor, Zone A, B, 10 <sup>th</sup> Floor, Building 3, Tinwe Business Park, No.6 Liufang Road, Baoan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA 402, Building 1, Shenzhen Biomedicine Innovations Industrial Park, No. 14, Jinhui Road., Jinsha Community, Kengzi, Pingshan District, 518000 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
<b>Nomor Sertifikat Pendaftaran Produk</b>	国械注准 20223060974
<b>Nomor Persyaratan Teknis Produk</b>	国械注准 20223060974
<b>Nomor Lisensi Manufaktur</b>	粤食药监械生产许 20204044 号
<b>Versi IFU</b>	01
<b>Tanggal Rilis IFU</b>	2023-06-26
<b>Masa Penggunaan Produk</b>	2 tahun

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