
C €0123

**VivoCardio Cardiovascular Pressure Measurement (CPM)  
System**

**INSTRUCTION FOR USE  
(Indonesian)  
INSTRUKSI PENGGUNAAN**

*Insight Lifetech*

## CONTENTS

1. Product Summary.....	1
1.1 Product Description.....	1
1.2 Intended Use .....	1
1.3 Indication for Use .....	1
1.4 Contraindications .....	1
1.5 Adverse Events .....	2
1.6 Users.....	2
1.7 Model .....	2
1.8 Warnings, Precautions, Cautions and Safety Information.....	2
1.9 Symbols used on product, package and IFU.....	7
1.10 CPM System Introduction and Symbols .....	8
2. Installation.....	11
2.1 Initial Inspection.....	11
2.2 Installation.....	11
2.2.1 Install CPM System on a trolley .....	11
2.2.2 Install CPM System on a bedside holder .....	11
2.3 Connect CPM System with the power cable .....	11
2.4 Powering On and Shutting Off.....	12
3. FFR/cRR Measurement.....	13
3.1 FFR/cRR Measurement Screen.....	13
3.2 FFR/cRR Measurement Procedure .....	15
4. Create Patient Data .....	20
4.1 Patient Window.....	20
4.2 Create, Edit, Delete Patient Information .....	21
4.3 Create a new exam for existing patient.....	21
4.4 Rank Function.....	21
4.5 Button Display .....	21
5. Review Recordings .....	22
5.1 Review Window.....	22
5.2 Review Recordings .....	23
5.3 Export Recordings .....	24
6. Review Archive Studies .....	25
6.1 Archive Window .....	25
6.2 Review Files in Archive .....	26
6.3 Edit and Delete Files in Archive .....	26
6.4 Export Files from Archive .....	26
6.5 Rank Function.....	27
6.6 Single Selection, Multi Selection, Select All Functions.....	27
7. Settings.....	28
7.1 User settings.....	28
7.2 Institution Setup.....	30
8. Troubleshooting .....	31
8.1 On screen messages.....	31


8.2 Tracing EM interference .....	31
9. Cleaning and Maintenance .....	32
9.1 Cleaning.....	32
9.2 Maintenance .....	32
9.2.1 Review .....	32
9.2.2 Maintenance Schedule.....	33
10. Warranty and After-sales Service .....	34
10.1 Warranty .....	34
10.2 After-sales Service.....	34
11. Technical Specifications.....	35
12. EMC Information .....	37

# 1. Product Summary

## 1.1 Product Description

VivoCardio® Cardiovascular Pressure Measurement (CPM) System (hereinafter called CPM System) is a diagnostic computer designed to record, compute, display, and store data from TruePhysio® Rapid Exchange Pressure Microcatheter (hereinafter called Pressure Microcatheter) and other external transducers. The information is displayed as graphs and numerical values on the screen.

CPM System includes a Main Unit with a touch screen user interface and several input/output ports, a remote control, and a power cable.

**Figure 1. Cardiovascular Pressure Measurement (CPM) System**

## 1.2 Intended Use

CPM System is intended for use in catheterization and related cardiovascular specialty laboratories to compute and display systolic pressure, diastolic pressure, mean pressure, FFR and cRR based on the output from one or more electrodes, transducers or measuring devices.

## 1.3 Indication for Use

CPM System is indicated to provide hemodynamic information for the diagnosis and treatment of patients that undergo measurement of blood pressure with Pressure Microcatheter.

## 1.4 Contraindications

### CPM System

- CPM System has no patient alarm functions. Do not use it for cardiac monitoring.
- CPM System is contraindicated for use with the high frequency surgical equipment.

### Pressure Microcatheter

- Pressure Microcatheter is contraindicated for use in the cerebral vasculature.
- Pressure Microcatheter is contraindicated for use in the ventricles if the patient has a prosthetic mechanical valve.
- Pressure Microcatheter is contraindicated for neonatology/pediatric use.
- Pressure Microcatheter is contraindicated for use with the high frequency surgical equipment.

## **1.5 Adverse Events**

### **Pressure Microcatheter**

Potential complications which may be encountered during all catheterization procedure include: vessel dissection or occlusion, perforation, embolus, spasm, local and/or systemic infection, pneumothorax, congestive heart failure, myocardial infarction, hypotension, chest pain, renal insufficiency, serious arrhythmias or death.

## **1.6 Users**

CPM System should be handled by or supervised by a physician with training in catheterization laboratory procedures.

## **1.7 Model**

U152301

## **1.8 Warnings, Precautions, Cautions and Safety Information**

### **Warnings:**

#### **CPM System**

- Do not use CPM System in home healthcare environment.
- CPM System should be handled by or supervised by a physician with training in catheterization laboratory procedures.
- No modification of this CPM System is allowed.
- Do not open or remove access covers on CPM System unless specifically instructed by Insight Lifetech technical support to do so.
- Do not use CPM System if it has been dropped or in any other way exposed to mechanical or electrical damage. This can cause the user or patient to be exposed to electric shock or cause faulty readings of data. Contact Insight Lifetech Co., Ltd. or local distributor for further actions in above situation.
- Do not use CPM System if liquid is suspected to have penetrated the enclosure or the power supply unit. This can cause the user or patient to be exposed to electric shock. Contact Insight Lifetech Co., Ltd. or local distributor for further actions in above situation.
- Do not sterilize CPM System or any of its parts. Do not use CPM System or any of its parts if it has been sterilized.
- CPM System is only intended to be connected to Pressure Microcatheter (Model: 1-14-1).
- CPM System contains a lithium battery for the system real-time clock and is in danger of explosion. The battery is not intended to be replaced by users.
- When CPM System is connected with other medical equipment, an equipotential electric field should be formed among them.
- External equipment intended for connection to signal input, signal output or other connectors, shall comply with relevant IEC standard (e.g., IEC 60601 series for medical electrical equipment). In addition, all such combinations of systems shall comply with the standard IEC 60601-1-1, Safety requirements for medical electrical systems. Any external equipment intended for connection to signal output, or other connectors, shall comply all the related safety standards. If in doubt, contact Insight Lifetech Co., Ltd. or the local distributor.

- If a nasty smell is found from the battery of the remote control, please do not use the battery.
- CPM System and its recyclable accessories must be handled in accordance with local regulation or returned to manufacturer for recycling after expiry date. The wasted battery is dangerous. Please do not put the wasted battery together with domestic waste. When the battery is expired, please transfer it to a designated place for recycling. Improper disposal methods may cause the environmental pollution. Please contact local government to get more information for recycling of CPM System and the battery.
- Ensure that the total leakage current is within their allowable value when the patient is connected with multiple devices. The patient may suffer from electric shock.
- Before using equipment, please confirm if some damages exist in the CPM System or its accessories (e.g., P<sub>a</sub> output switch cable and AO input switch cable), because the damages may compromise the safety of patients. If the obvious signs of aging or damage is found, please replace this component with a new one to ensure the normal and safe operation of the CPM System.
- The performance of CPM System may be disturbed by electromagnetic environment. Please ensure that all other medical equipment close to the CPM System comply with the electromagnetic compatibility requirements. Because X-Ray or MRI equipment can release large amounts of electromagnetic radiation, they may be the source of electromagnetic interference.
- CPM System is a Type CF equipment and protected against the effects of a discharge of a defibrillator. CPM system will be restored after 10 seconds if defibrillation happens. Pressure Microcatheter readings may be affected by defibrillation. After defibrillation restart procedure, i.e. re-zero and re-equalize Pressure Microcatheter.
- To avoid the risk of electric shock, the CPM System must only be connected to a supply mains with protective earth.
- The plug (or appliance separable) is used as disconnect to the mains supply.
- Do not touch the SIP (signal input ports)/SOP (signal output ports) and the patient simultaneously.
- Do not touch the switching power adapter and the patient simultaneously.
- Do not touch the PressureCath IN or AO IN interface if the Pressure Microcatheter or AO transducer is disconnected from CPM System, otherwise the user may suffer from electric shock during defibrillation.
- Please restart the measuring process according to the Measurement procedure if abnormal waves are observed during measurement.

### **Pressure Microcatheter**

- Pressure Microcatheter is only intended to be connected to CPM System (Model: U152301) manufactured by Insight Lifetech. This CPM System complies with relevant IEC 60601 series standards.
- Pressure Microcatheter is intended for single use only. Do not resterilize or reuse. Reuse After cleaning attempts, resterilization and repackaging may result in patient/user infections and inaccurate pressure signals, thereby exposing the patients to increased risks for complications/adverse events.
- Positioning of guiding catheters, guidewires and Pressure Microcatheter in the ventricles is

potentially arrhythmogenic. It should never be done without ECG monitoring and presence of a functioning defibrillator.

- Pressure Microcatheter should only be used by physicians trained and experienced in interventional procedures.
- Do not use Pressure Microcatheter if the original sterile package is not complete. Inspect the device carefully prior to use to verify that all parts are present and undamaged.
- Do not proceed with the procedure when abnormal pressure value and waveform of  $P_a$  are exhibited on the CPM screen after connection of AO transducer with the AO IN signal cable. Any existence of bubbles in the tubing system of AO transducer may cause the abnormal waveform and errors in pressure readings. Please flush the tubing with saline solution if bubbles exist. If abnormal pressure value and waveform of  $P_a$  still exist after flushing, please replace the AO transducer with a new one.
- Do avoid the incarceration of coronary artery when engaging the guiding catheter to the ostium of targeted coronary artery. Please adjust the position of guiding catheter carefully and avoid deep engagement.
- Pressure Microcatheter should be handled carefully. Bending or excessive force during removal from carrier tubing may damage Pressure Microcatheter.
- Prior to use and when possible during the procedure, inspect Pressure Microcatheter carefully for bends, kinks or other damage. Do not readjust any bend or kink.
- Do not use damaged Pressure Microcatheter in any way, otherwise vessel/ventricle damage and/or inaccurate measurement may occur.
- Do not start the equalization procedure when the Pressure Microcatheter is inside of the guiding catheter. Please advance the Pressure Microcatheter and locate the sensor completely out of the tip of the guiding catheter for equalization.
- Observe all Pressure Microcatheter movements. Whenever Pressure Microcatheter is moved or torqued, the marker band should be examined under fluoroscopy.
- Excessive manipulation when sensor element or tip of Pressure Microcatheter is located in sharp bend, may cause damage or tip fracture. If the guiding catheter is in an anatomically severe or sharp bend, for example, a tortuous subclavian artery or adjacent vessel position, Pressure Microcatheter may be vulnerable to kinking or fracture.
- Always advance or withdraw Pressure Microcatheter slowly. Never advance or withdraw with excessive force if it meets resistance.
- Do not proceed with the procedure when abnormal pressure value and waveform of  $P_d$  are exhibited on the CPM System screen after connection of Pressure Microcatheter with the PressureCath IN interface. Please adjust the position of Sensor of Pressure Microcatheter carefully and slightly when abnormal waveform is exhibited. Avoid bending or kinking of the Pressure Microcatheter during adjustment.
- Please use Pressure Microcatheter before the expiry date. Otherwise it may cause patient infection.
- Please stop the measurement immediately if the intolerance to adenosine or ATP is observed.
- Do not record the FFR/cRR value when abnormal pressure value and waveform of  $P_d$  are exhibited on the CPM System screen. Adjust the position of the sensor by withdrawing or advancing the Pressure Microcatheter slightly to get normal readings.
- Do not record the FFR/cRR value when drift occurs. Please check if drift occurs by pulling the

Pressure Microcatheter back to the equalizing position and evaluate if  $P_d/P_a$  value is out of the range of 0.97~1.03. If the value is out of the range, please equalize the system again and restart the measuring process.

### **Precautions:**

#### **Pressure Microcatheter**

- Before using the Insight FFR/cRR System which consists of Pressure Microcatheter and CPM System, please read Pressure Microcatheter's IFU and CPM System's IFU and be sure that you understand them fully.
- Refer to instructions supplied with any interventional devices to be used in conjunction with Pressure Microcatheter for their intended uses, contraindications and potential complications.
- Confirm the compatibility of Pressure Microcatheter with other interventional devices (e.g. guiding catheters and guide wires). Free movement of Pressure Microcatheter along the guide wire within the guiding catheter is important. Make sure that Pressure Microcatheter does not meet excessive resistance during movement when used along with other interventional devices.
- Adjust or replace the hemostatic valve with an adjustable valve if it is found to inhibit the Pressure Microcatheter movement.
- Failure to achieve maximum coronary hyperemia may result in inaccurate FFR readings. The improper position of additional interventional devices, such as balloon catheters, may also result in inaccurate readings.
- Do not measure the blood pressure when the sensor of the Pressure Microcatheter is bent, since false pressure readings may occur.
- The Pressure Microcatheter readings may be affected by defibrillation. Recalibrate Pressure Microcatheter after defibrillation use.
- Make sure that the sensor does not touch the atrial or ventricular walls during movement, to avoid measurement deviation.
- Used Pressure Microcatheter may be a potential biological hazard. Handle and dispose the Pressure Microcatheter in accordance with medical practice, applicable local laws and regulations.
- Do not unpack the package too early, since it may cause patient infection.
- If the Pressure Microcatheter gets stuck or damaged during procedure, immediately disconnect the CPM System, and replace the damaged Pressure Microcatheter. The damage may include but not limited to: kinks, bends, sheared polymer coating, no pressure signal, or inaccurate pressure signal. Incorrect use may result in vessel/ventricle damage, arrhythmia induction, electric shock, or inaccurate pressure signals.

### **Cautions:**

#### **CPM System**

- If it is necessary to install CPM System adjacent to other equipment, the CPM System should be monitored to verify normal operation. If abnormal performance is observed, it is necessary

to increase the distance between the equipment.


- Please keep cleaning of the using environment and avoid vibration. The CPM System must be kept away from corrosives
- The applied part is not allowed to be contacted with other conductors (e.g., metal tap), because it will compromise the safety of CPM System.
- When the CPM system and its accessories are close to the expiry date, they must be handled according to local regulations or hospital regulations.
- When the battery is close to its expiry date, please replace it with a new one.

### **Safety Information:**

#### **CPM System**


- Please install the CPM system in a place where it can be observed, operated, and maintained easily.
- Please do routine preventive maintenance every two years. Besides, all special requirements in local regulations should be met.
- Please operate the CPM system under declared operation conditions.
- Please handle the expired CPM System and recyclable accessories in line with local regulations, or return to manufacturer for recycling usage.
- The wasted battery of remote control should be handled according to local regulations.
- Please use the AO output cable, P<sub>d</sub> output cable and other accessories from the manufacturer.
- The battery type of the remote control is AAA (DC, 1.5V). Make sure that the two batteries are replaced simultaneously if needed. Take out of batteries if they are not used for a long time. Please do not use the remote control if battery leakage happens. Please do not look straight at the infrared ray of the remote control to avoid potential visual impairment.
- Please use a special ground cable with yellow and green color during use, and connect the equipotential terminal to the equipotential network used in the catheter lab.
- Once the power system is involved in accident, the zeroed data can be recovered after restarting within 5 minutes.


## 1.9 Symbols used on product, package and IFU

Symbols	Explanation
	Keep dry
	This side up
	Fragile, handle with care
	Maximum stacking layer
	Caution
	Serial number
	Date of manufacture
	Defibrillation-proof type CF applied part
	“OFF” for a part of equipment
	“ON” for a part of equipment
	Alternating current
	Electronic waste bin. Disposal according to EU’s WEEE directive 2002/96/EC
	Manufacturer
	Unique device identifier
	Medical Device
	Importer
	Authorized Representative in the European Community
	Swiss Authorized Representative

Symbols	Explanation
	Refer to instruction for use

### 1.10 CPM System Introduction and Symbols

**Figure 2. CPM System and Remote Control**

No.	Symbol or Marking	Explanation
1	N/A	Main Unit
2	N/A	Display-Touch Screen
3	⊖/⊕	“OFF” for a part of equipment/“ON” for a part of equipment
4	N/A	Infrared Receiver
5	PressureCath™ IN	PressureCath IN on Pressure Microcatheter Interface-P <sub>d</sub> signal from Pressure Microcatheter. Connection to Applied Part.
6	N/A	Slot for storing Remote Control
7	N/A	Slot for storing Power Adapter
8	AO IN	AO IN on AO transducer Interface-P <sub>a</sub> signal from AO transducer. Connection to Applied Part.
9	PressureCath™ OUT	PressureCath OUT on Pressure Microcatheter Interface-P <sub>d</sub> signal to external monitor.
10	AO OUT	AO OUT on AO transducer Interface-P <sub>a</sub> signal to external monitor.
11	USB	USB interface

No.	Symbol or Marking	Explanation
12	VGA	VGA interface
13	N/A	Power Supply Cable
14		Equipotentiality
15	N/A	Remote Control
16		Move left button: Press this button to move to the previous record
17		Move right button: Press this button to move to the next record
18	OK	Enter button: Confirm the selection of menu button or function, or confirm entered value
19		Marker Button: Add marker during recording
20		Screenshot button: Intercept current screen and save it as a picture
21		Starting Recording, Stopping Recording function button
22	EQUALIZE	Equalize button: Equalize $P_a$ and $P_d$ value
23		Move to left menu button
24		Move to right menu button

## 2. Installation

### **Note:**

Please read this IFU and install the CPM System per this IFU to ensure that the CPM System can run normally.

### 2.1 Initial Inspection

Make sure the package is free from any damage before opening it. If there are any damage, please ask for claims against the carrier.

Open the package, then take out CPM System and its accessories carefully. Check all components have been supplied as the packing list. Finally, check whether all parts are undamaged. Contact Insight Lifetech Co., Ltd. or your local distributor for any question.

### Packing List

No.	Components Description	Qty (EA)
1	Main Unit	1
2	Remote Control	1
3	Instruction for use	1
4	AO input switch cable	1
5	P <sub>a</sub> output switch cable	1
6	P <sub>d</sub> output switch cable	1
7	Power Supply Cable	1
8	FFR&cRR cutoff value sticker	1

### 2.2 Installation

CPM System can be installed on a trolley or bedside holder for use.

#### 2.2.1 Install CPM System on a trolley

Please install CPM System on a trolley according to the trolley's installation instruction.

### **Warning:**

Make sure CPM System is securely mounted. If CPM System is not securely attached, it may fall over causing possible patient or operator injury and damage to the system.

### **Caution:**

Make sure that screws are well tightened to prevent CPM System from falling.

#### 2.2.2 Install CPM System on a bedside holder

Please install CPM System on a bedside holder according to the bedside holder's installation instruction.

### 2.3 Connect CPM System with the power cable

Plug the attached power cable into a power outlet, and connect the power adapter with the power cable. Connect in such way that the power cable can be easily removed from power outlet.

### **Note:**

*The rated voltage and frequency are 100-240 V and 50/60 Hz.*

*Plug the attached power cable into the specified power outlet. Make sure that power connection is stable to avoid that the power is interrupted abnormally during measurement.*

## **2.4 Powering On and Shutting Off**

### **Powering On**

Check that the power cable is connected with the power adapter and to the mains power outlet.

Press the On/Off switch to turn on CPM System. The screen will remain dark for a few seconds before a startup beep indicates that the system is booting up.

### **Shutting Off**

Long press the On/Off switch for 3 seconds to turn off CPM System.

### **Caution:**

The power still exists in the power cable, power adapter and some parts of main unit after turning off the CPM System. Pull out the power plug from the power outlet to get a complete power disconnection.

### **Note:**

*If the cable of power adapter is separated from power adapter, please contact Insight Lifetech Co., Ltd. or your local distributor for further instructions.*

*If the power cable does not work, please contact Insight Lifetech Co., Ltd. or your local distributor to replace for a new one.*

### 3. FFR/cRR Measurement

#### 3.1 FFR/cRR Measurement Screen


Figure 3. Measurement Screen (Live Window)

No.	Explanation
1	Information Bar (Patient Name, Patient ID, Catheter Room No., Date, Time)
2	P <sub>a</sub> , P <sub>d</sub> phasic and mean pressure
3	P <sub>a</sub> -P <sub>d</sub> trace or FFR trace
4	Mode selection button FFR- Fractional Flow Reserve (Default) cRR- constant Resistance Ratio
5	P <sub>a</sub> Value: mean/systolic/diastolic values
6	P <sub>d</sub> Value: mean/systolic/diastolic values
7	Equalize Value: Difference value of P <sub>a</sub> and P <sub>d</sub> before Equalize
8	P <sub>d</sub> /P <sub>a</sub> Value
9	Exam Summary: Summary list of all the FFR/cRR measurements in Exam for the day.
10	Starting Recording, Stopping Recording function button
11	Marker Button
12	Screenshot Button: Intercept the current screen and save it as a picture
13	Setup Bar: User, Institution, Factory Settings
14	Patient Bar: Create New patient, edit and delete patient information, New Exam
15	Live Bar: Live information during measurement
16	Review Bar: Review measurement information
17	Archive Bar: Review patient's archive
18	Click on the triangle module-zero P <sub>a</sub>
19	Click on the triangle module-zero P <sub>d</sub>
20	Click on the triangle module-Equalize P <sub>a</sub> and P <sub>d</sub>
21	P <sub>d</sub> /P <sub>a</sub> Value
22	cRR Value
23	cRR trace

**Color coding:**

P<sub>a</sub> values and graph are red

P<sub>d</sub> values and graph are cyan

P<sub>a</sub>-P<sub>d</sub> trace is green  
P<sub>d</sub>/P<sub>a</sub> Values are yellow  
FFR values and graph are yellow  
cRR values and graph are blue

### 3.2 FFR/cRR Measurement Procedure

**Note:**

*Please also refer to the instruction for use of Pressure Microcatheter during FFR/cRR Measurement.*

Overview of a measurement procedure:

- 1) Prepare CPM System
- 2) Prepare AO transducer
- 3) Prepare Pressure Microcatheter
- 4) Equalize P<sub>a</sub> and P<sub>d</sub>
- 5) Record FFR or cRR
- 6) Withdraw Pressure Microcatheter from guiding catheter
- 7) Review Data (refer to 5 Review Data)
- 8) Export Data (refer to 6 Export Data)

- **Prepare CPM system**

- 1) Preparations

Observe normal sterile procedure if CPM System is used within the sterile zone. This includes covering CPM System in sterile plastic and keeping the remote control in a sterile plastic bag.

- 2) Start system

Plug the attached power cable into the power outlet, and connect the power adapter with the power cable. Press the On/Off switch to turn on CPM System, and make it warm up for 5 minutes. CPM System will display the Patient Window automatically.

- 3) Select or Create new patient information

Create new patient information: Touch **New** button, and input patient information.

Create a new exam for existing patients: Choose an existing patient record and touch **New Exam** button.

CPM System will display the Live Window automatically after the patient information is created.

- **Prepare AO transducer**

- 1) Connect AO transducer with the AO IN signal cable. The  $P_a$  area on the screen will display “Connected”. Flush the tubing with solution to clean all air in the tubing system to ensure that no bubbles exist in the sensor or valve.

**Warning:** If bubbles exist in the tubing system, please flush the tubing with saline solution again. Bubbles may cause errors in pressure readings.

- 2) Place the AO transducer at heart level, then open the AO transducer manifold to air, and press the **Zero  $P_a$**  button. The  $P_a$  information area on the screen will display “0”.

**Figure 4. Zero  $P_a$**

#### ● **Prepare Pressure Microcatheter**

- 1) Unpack Pressure Microcatheter using regular sterile routine but leave Pressure Microcatheter in the tray. Place the carrier tubing flat on the table.
- 2) Fill the carrier tubing with saline solution (approximately 20 cc) through the flush port to make sure the sensor element is immersed in saline.
- 3) Take out the cable assembly of Pressure Microcatheter and connect the cable with the PressureCath IN port of CPM System. The  $P_d$  area on the screen will display “Not Zeroed”.

- 4) Press Zero  $P_d$  button to zero Pressure Microcatheter, then the  $P_d$  area on the screen will display “0”.

**Note:**

*If the  $P_d$  area displays “Zero Failed” instead of “0”, please disconnect and then connect the cable assembly with the PressureCath IN port of CPM System again. Otherwise, replace the Pressure Microcatheter with a new one.*

- 5) Pull out the Rapid Exchange Pressure Microcatheter carefully from carrier tubing.

**Figure 5. Zero Pd**

- **Equalize  $P_a$  and  $P_d$** 
  - 1) Advance Pressure Microcatheter along with guidewire carefully until the sensor element of Pressure Microcatheter is just outside the tip of the guiding catheter. Observe the marker position of Pressure Microcatheter through fluoroscopy.
  - 2) Press **Equalize** to make  $P_a$  equal to  $P_d$ .


Figure 6. Equalize

- **Measurement Mode**

**Measure FFR**

- 1) Advance Pressure Microcatheter to the distal part of lesion.
- 2) Inject adenosin or ATP from coronary or vein to induce maximum hyperemia.
- 3) Touch **Record** button to start recording.
- 4) To place a marker in the recording at any time, touch the Marker button.
- 5) Touch **Stop** button to stop recording until a steady state of maximum hyperemic condition is reached, then the measurement is completed. The recording is automatically saved and the Review window is opened. CPM system will display the minimum FFR value automatically.
- 6) Please click on the triangle module  if it is necessary to zero  $P_a$  or  $P_d$  during the procedure.

## Measure cRR

- 1) Advance Pressure Microcatheter to the distal part of lesion.
- 2) Touch **Record** button to start recording.
- 3) Record the region of interest.
- 4) To place a marker in the recording at any time, touch the Marker button.

*Note: When the cRR waveform curve appears a dotted line, it indicates that the current  $P_a$  and  $P_d$  waveform is abnormal and cannot be used to calculate cRR.*

- 5) Touch the **Stop** button. The recording is automatically saved and the Review window is opened. CPM system will display the minimum cRR value automatically.
  - 6) Please click on the triangle module  if it is necessary to zero  $P_a$  or  $P_d$  during the procedure.
- **Withdraw Pressure Microcatheter from guiding catheter.**


Figure 7. FFR/cRR Measurement

## 4. Create Patient Data

### 4.1 Patient Window

**Figure 8. Patient Window**

No.	Explanation
1	New Patient Button: Create Name, Patient ID, Birth, Gender and Doctor information.
2	Edit Button: Edit selected patient information.
3	Delete Button: Delete selected patient record.
4	New Exam Button: Return to Live interface.
5	Last Page Button: Move to the last page by pressing this button.
6	Next Page Button: Move to the next page by pressing this button
7	Page 1 of 1 (2) : Page information and Total patient quantity

## 4.2 Create, Edit, Delete Patient Information

New Patient: By pressing the **New Patient** button, you can input or select Patient Name, Birthday ID, Gender information, and doctor information to create new patient information. The system will automatically generate a time-stamped patient ID by default, which is convenient for users to quickly create new patients.

After the CPM System is turned on, the **New Exam** button is displayed as **Quick Exam**, and a time-stamped patient ID can be automatically generated by clicking it, which is convenient for users to quickly enter the real-time page to perform measurement procedures. One quick exam can be created only each time the CPM system is turned on.

### *Note:*

*The maximum length of Patient Name and ID is 64 characters, and doctor information can only be edited and deleted in the setup bar.*

Edit Patient: Select the target patient record and press the **Edit** button. The Patient Name, ID, Birthday and Gender information will be available for editing.

Delete Patient: Select the target record and press the **Delete** button, then the interface will pop up a warning window. The target record will be deleted by pressing the **OK** button.

## 4.3 Create a new exam for existing patient

New Exam: Select the target patient record, the system will switch to Live Window by pressing the **New Exam** button. For the newly created patient, click the information bar of Patient Bar to anonymize the patient name and patient ID; click again to cancel the anonymity of the patient information.

## 4.4 Rank Function

Rank function: By touching the column name of Last Name, First Name, Patient ID, Date of Birth, Gender or Doctor, you can put them in alphabetical or numerical order.

## 4.5 Button Display

Button display: Highlighted button is available. Dimmed button is not available.

## 5. Review Recordings

### 5.1 Review Window


Figure 9. Review Window

No.	Explanation
1	FFR cursor: Indicates FFR value
2	$P_d/P_a$ box: Indicates $P_d/P_a$ value at purple cursor position
3	FFR value box: Indicates FFR value at yellow cursor position. Reposition FFR value by long touching the FFR value box.
4	Exam summary box: Summary list of recordings in Exam.
5	Delete/Export button: Delete the recordings from Exam summary box; Export the recordings to a USB Memory Stick
6	Live button: Return to live window
7	Annotation button: Open annotation window
8	Mark button: Mark one in recordings. Delete all marks by long pressing the Mark button
9	Screenshot: Intercept the current screen and save it as a picture, and the button will show numbers if screenshot is performed successfully. Screenshots can only be viewed after exporting.
10	Zoom in
11	Zoom out
12	Page turning left button: turn the graph to the left pages
13	Page turning right button: turn the graph to the right pages
14	X-axis: measurement time ( shown in second unit)
15	Y-axis: $P_a$ value, $P_d$ value, $P_a-P_d$ value, FFR value
16	cRR box: Indicates cRR value at purple cursor position
17	cRR value box: Indicates cRR value at yellow cursor position. Reposition cRR value by long touching the cRR value box.

## 5.2 Review Recordings

The Review window displays the recording with a yellow FFR or cRR cursor line indicating the FFR or cRR value, which was defined as the point of the lowest  $P_d/P_a$  or cRR value during the recording cycle.

### Caution:

**The system may place the cursor at the wrong location (value) due to artifacts in  $P_a$  or  $P_d$  pressure. The responsible physician should confirm that the point selected by the system is a valid point for the selected index.**

Adjust FFR: Review the recording and place the index cursor by moving the cursor with a finger or the arrow buttons. If necessary, reposition the cursor at a new position to calculate a new index value. Then long press value box, and press the **OK** button upon the popup window. The recording will be updated with the adjusted index value.

Add an additional secondary cursor to view an additional index at another point in the recording.

Annotations: Annotations are entered via the **Annotation** button. Select the annotation for Vessel, Step and Drug from the respective lists. A free text comment can also be added. Use the **OK** button to save entries.

**Note:**

*The Vessel, Step, and Drug can only be added or deleted in the setup bar.*

Review and delete other recordings in Exam

All recordings in exam are listed on the exam summary box with procedure time and corresponding value. Select the recordings of interest for review. Touch the **Delete** button to delete the recordings. Start a new measurement for the same patient by pressing **Live** button.

Review, add, and delete marks in Exam

Select the recording of interest for mark. The purple line in graph is the mark line. Touch the interesting position in graph, and press the **Mark** button for marking. Delete all marks in the selected recording by long press the **Mark** button.

### 5.3 Export Recordings

1. To USB Memory

Insert a USB Memory Stick into the USB interface. Select the recordings of interest for review, then touch the **Export** button to export recordings to the USB Memory Stick.

**Note:**

If recording export fails, please check as the following steps:

- Is an external storage device connected?
- Do you use an available external storage device?
- Is the capacity of an external storage enough?

2. To DICOM Archive

Connect with PACS System by inserting a compatible hub into the USB interface, then set the local network IP and the IP of the connected PACS system (export IP). Select the recordings of interest for review, touch the **Export** button to export recordings to DICOM Archive, then the recordings will be synchronized to the PACS system.

**Note:**

If recording export fails, please check as the following steps:

- Is the external network port connected securely?
- Is the network setting correct?
- Is the network interrupted?

Please contact Insight Lifetech Co., Ltd. or your local distributor for the type of hub.

**Caution:**

The USB interface only supports regular USB Memory Sticks. It is not allowed to connect other devices.

## 6. Review Archive Studies

### 6.1 Archive Window


Figure 10. Archive Window

No.	Explanation
1	Edit button: Edit patient name, birth, gender and doctor
2	Delete button: Delete selected Exam
3	Export button: Export the selected Exam to a USB Memory Stick
4	Last Page button: Move to the last page by pressing this button
5	Next Page button: Move to the next page by pressing this button
6	Page 1 of 1 (2) : Page information and Total patient quantity

## 6.2 Review Files in Archive

Review files: Press the **Archive** button, select a patient information, and then press the **Review** button.

### **Note:**

*The patient files cannot be reviewed when Patient button is selected. Patient interface can only be used to manage the patient information but not for review.*

## 6.3 Edit and Delete Files in Archive

Edit recording: Press the **Edit** button to edit the recordings of interest. The patient name, birth, gender, and doctor can be edited.

Delete recording: Press the **Delete** button to delete the recordings of interest.

### **Note:**

*The recordings in archive are listed per exam date and patient information, so one patient may have several exams in one day. The recordings in review are listed per each exam.*

## 6.4 Export Files from Archive

### 1. To USB Memory

Insert a USB Memory Stick into the USB interface. Select the recordings of interest for review, then touch the **Export** button to export recordings to the USB Memory Stick.

### **Note:**

If recording export fails, please check as the following steps:

- Is an external storage device connected?
- Do you use an available external storage device?
- Is the capacity of an external storage enough?

### 2. To DICOM Archive

Connect with PACS System by inserting a compatible hub into the USB interface, then set the local network IP and the IP of the connected PACS system (export IP). Select the recordings of interest for review, touch the **Export** button to export recordings to DICOM Archive, then the recordings will be synchronized to the PACS system.

### **Note:**

If recording export fails, please check as the following steps:

- Is the external network port connected securely?
- Is the network setting correct?
- Is the network interrupted?

Please contact Insight Lifetech Co., Ltd. or your local distributor for the type of hub.

### **Caution:**

The USB interface only supports regular USB Memory Sticks. It is not allowed to connect other devices.

### **6.5 Rank Function**

By touching the column name of Last Name, First Name, Patient ID, Gender, Doctor, Exam Date, Times, FFR, cRR or Export, you can put them in alphabetical or numerical order.

### **6.6 Single Selection, Multi Selection, Select All Functions**

Touch select all/deselect all to select/deselect all recordings. Touch the square box to select more than one recordings. Under any condition, if you would like to select only one recording, just touch any position of this recording other than the square box.

## 7. Settings

### 7.1 User settings

The screenshot shows the 'User Setup Window' of the Insight Lifetech software. The window has a dark background with white text. At the top, there is a header bar with the text 'test 001 Room 09-01-2022 09:55:29' and the Insight Lifetech logo. Below the header, the settings are organized into two columns. The left column contains: Date (09-01-2022), Time (09 : 55), Language (English), Screen Light (100%), Speed(mm/s) (25), Scale(mmHg) (0 ~ 150), Average(beats) (3), and Max Length(min) (60). The right column contains: Test Pa (Off), Test Pd (Off), Display (Grid checked, Pa-Pd Trace checked), Available/Capacity: 1.52GB/1.53GB, and Software Release Version: V2. Below these settings are two buttons: 'Network Settings' and 'Restore Factory Settings'. At the bottom of the window, there is a 'User' tab and an 'Institution' tab, both with text input fields. To the right of these tabs is a 'Save' button. At the very bottom of the window, there is a navigation bar with icons for 'Patient', 'Live', 'Review', and 'Archive', and a large blue 'Live' button on the right.

**Figure 11. User Setup Window**

In this interface, you can:

- Set date and time: Input current date and time.
- Set user interface language: Chinese, English, Polish, Italian, Spanish and German are available. The default language is English.
- Set screen light: 20%, 40%, 60%, 80% and 100% are available. The default screen light is 100%.
- Set speed: 5 mm/s, 10 mm/s, 25 mm/s, and 50 mm/s are available. The default speed is 25 mm/s.
- Set scale: -20 ~ 80, 0 ~ 100, 0 ~ 150, 0 ~ 200, 0 ~ 250, -80 ~ 320 mmHg are available. The default scale is 0 ~ 150 mmHg.
- Set average (beats): Number of heart beats included in pressure averaging filter: 1, 2, 3, 4, 8 and 12 beats are available. The default setting: 3 beats.

**Note:**

*Choosing a high number of heartbeats will make the pressure averaging slower and less sensitive to artifacts, but it may also result in insensitive averaging which is noticeable when*

*there is a short hyperemic plateau. Choosing a low number of heartbeats makes the pressure averaging faster and more sensitive to pressure changes, which is desirable when using a short hyperemic plateau, but it may also result in an averaging overly sensitive to arrhythmia and pressure disturbances.*

- Set the maximum recording time for one recording: 5 minutes, 10 minutes, 30 minutes, and 60 minutes are available. The default setting is 60 minutes.

**Note:**

*If a recording passes the maximum recording time, the recording will automatically stop, and a new one will automatically start.*

- Set test  $P_a$ ,  $P_d$ : -30mmHg, 0mmHg, 50mmHg, 100mmHg, 200mmHg, and 300mmHg are available. The default setting is Off.

**Note:**

*This function is used for testing the consistency in displaying value and selecting value, and the output accuracy of  $P_a$  and  $P_d$  signal.*

- Set display mode: Grid mode and  $P_a$ - $P_d$  trace
- Restore default settings: By pressing this button, the system will reserve patient data, and set other parameters to the factory default settings. If the user chooses to delete all patient data, all existing patient data will be deleted after restoration.
- CPM system can store more than 1000 data.

**Note:**

*After the settings is completed, the system will move to the Live Window by pressing the **Save** button.*

## 7.2 Institution Setup

test 001 Room 09-01-2022 09:55:45 Insight Lifetech

Room

Add Doctor

All

Add/Del Doctor

Comment  
Vessel ▼

Add Comment

Current  
None  
RCA Prox  
RCA Mid  
RCA Distal  
PDA  
Left Main  
LAD Prox  
LAD Mid  
LAD Distal  
Other

Option

>  
<

Add Comment

Delete Comment

User Institution Save

⚙️ Patient Live Review Archive Live

Figure 12. Institution Setup Window

In this interface, you can:

- Room: Input Room information.
- Add doctor information: Input a doctor name, and press the **Add Doctor** button, then the doctor can be added in All box.
- Delete doctor information from the list: Select the doctor name of interest and press the **Delete Doctor** button, and the selected doctor name will be deleted from the list.
- Annotations: Vessel, Step, and Drug are available. Select **Vessel/Step/Drug**, then input information and press the **Add Comment** button, the input information will be added in the Current box. Select the interested information from the Current box and press the **Delete** button, and the selected information will be deleted from the Current box. The information in the Current box can be moved to the Option box.

**Note:**

*The existing vessel, step, and drug information in original system are not able to be deleted. After you complete the settings above, the system will enter the Live Window by pressing the **Save** button.*

## 8. Troubleshooting

### 8.1 On screen messages

When the error codes in the following table appear on the screen, please take actions according to the following instructions.

Error Code	Explanation	Possible Cause	Action
Error 01	Electrical ratings of +12 V DC are exceeded.	Electronic components damaged	Contact Insight Lifetech
Error 02	Electrical ratings of +6 V DC are exceeded.	Electric components damaged	Contact Insight Lifetech
Error 03	Electrical ratings of +5 VA DC are exceeded.	Electric components damaged	Contact Insight Lifetech
Error 04	Electrical ratings of +3.3VDD DC are exceeded.	Electric components damaged	Contact Insight Lifetech
Error 05	Electrical ratings of +5VF DC are exceeded.	Electric components damaged	Contact Insight Lifetech

### 8.2 Tracing EM interference

The CPM System is susceptible to interference generated by other RF energy sources such as medical devices, information technology products, or radio/television transmission towers.

Tracing the source of radiated interference can be difficult. In accordance with the standards identified in these instructions for use, no interference was observed. However, the trained user must determine if an artifact caused by radiated interference will negatively impact signal quality and the subsequent exam results.

To help identify the source of electromagnetic interference, ask the following questions:

- Is the power cable or grounding wire connected with protective earth?
- Is the interference intermittent or constant?
- Is the interference present if the CPM System is moved to a different location in the facility?

Moving cables or other medical equipment away from the system can reduce electromagnetic interference.

Please answer these questions before contacting your service representative. The answers will help a service representative determine if the problem is in the system or in the imaging environment.

## 9. Cleaning and Maintenance

### 9.1 Cleaning

#### Warning:

- Before cleaning CPM System, turn off the power supply, and disconnect the power cable from the power plug.
- Make sure that every part of the unit is dry before reconnecting it to the power supply.
- Do not use CPM System if liquid is suspected to have penetrated the enclosure or the power supply unit. This can cause the user or patient to be exposed to electric shock. Contact Insight Lifetech Co., Ltd. or your local distributor for further instructions.
- Clean CPM System once a week.

#### Caution:

Clean CPM System with cleaning agents. Using a clean soft cloth or tissue to wipe the residual cleaning agent.

The following cleaning agents are verified to be used for cleaning CPM System:

- Alcohol (75%)
- Isopropanol (70%)

Clean the surface of CPM System as the following steps:

1. Turn off the power supply and disconnect the power cable.
2. Use a soft cloth to adsorb some cleaning solution, and wipe all the surfaces of CPM System including the Screen.
3. Use a new cloth or tissue to adsorb some water to wipe the residual cleaning solution.
4. Dry CPM System in a cool place.
5. Reconnect the power supply, and press the on/off switch to start CPM System.

### 9.2 Maintenance

#### Warning:

The maintenance schedule should be performed by the hospital or medical agents who use this CPM System, or it may cause CPM System failure or unexpected consequence, or a threat to personnel safety.

#### 9.2.1 Review

To make sure that CPM System can run normally, it is recommended that one comprehensive review should be performed by the maintenance personnel in the following conditions:

- Before using CPM System
- After using CPM System for 24 months
- After completing maintenance and upgrade

The review items are listed as follows:

- Environment and the power supply comply with requirements.
- The power cable is in good conditions.
- The equipment has no mechanical damage.
- Use indicated accessories.
- The functions of CPM System are OK.

- Leakage current complies with requirements
- Batteries of the remote control are available, and take out the batteries if they are not used for a long time. Do not use the remote control if battery leakage or a nasty smell happen.

### 9.2.2 Maintenance Schedule

The following schedule should always be carried out by the maintenance personnel. Contact the maintenance personnel for helping if the following maintenance is necessary. Before testing or maintaining, CPM System should be cleaned.

Maintenance project	Frequency
Leakage current measurement according to IEC 60601-1	Once every two years. It should also be performed if CPM System falls off.

## **10. Warranty and After-sales Service**

### **10.1 Warranty**

Insight Lifetech Co., Ltd. warrants that CPM system meets the specifications of the product and will be free from defects in materials and workmanship during the warranty period. If CPM System covered by this warranty is determined to be defective because of defective materials, components, or workmanship, and the warranty claim is made within the warranty period, Insight Lifetech Co., Ltd. will, at its discretion, repair or replace the defective part(s) free of charge. Insight Lifetech Co., Ltd. will not provide a substitute product for use when the defective product is being repaired.

CPM System has no components that can be repaired by users themselves. All repairs should be performed by personnel authorized by Insight Lifetech Co., Ltd. The warranty is void in the case of:

- Components are disassembled, extended, or readjusted.
- Repair or modify CPM System by unauthorized personnel by Insight Lifetech Co., Ltd.
- The damage is caused by improper use or maintenance.
- Replace or remove original label or manufacturing mark.

### **10.2 After-sales Service**

Contact Insight Lifetech Co., Ltd. or your local distributor if you have any question during use of CPM System.

After-sales service agent (Manufacturer): Insight Lifetech Co., Ltd.

Contact Address:

Zone E, 3rd Floor, Building 3, Tinwe Business Park, No.6 Liufang Road, Baoan District, 518000  
Shenzhen, PEOPLE'S REPUBLIC OF CHINA

Tel: +86 (755) 23074845

<http://www.insight-med.com>

Email: [service@insight-med.com](mailto:service@insight-med.com)

## 11. Technical Specifications

### Electrical Specifications

Type of protection against electric shock	Class I
Degree of protection against electrical shock.	Defibrillation-proof type CF applied part
Operating voltage	100-240V, 50/60Hz
Rated Power Consumption	35VA

### Environmental Requirements

Requirement	Range-CPM System	Range-Pressure Microcatheter
Operating temperature	15 °C ~ 35 °C	15 °C ~ 35 °C
Operating relative humidity	30 %RH ~ 75 %RH	30 %RH ~ 75 %RH
Operating, ambient pressure	525 mmHg ~ 795 mmHg (700 hPa ~ 1060 hPa)	525 mmHg ~ 795 mmHg (700 hPa ~ 1060 hPa)
Transport temperature	-20 °C ~ +70 °C	-29 °C to +60 °C
Transport relative humidity	10 %RH ~ 95 %RH (non-condensing)	5 % RH ~ 85 % RH (non-condensing)
Transport ambient pressure	375 mmHg ~ 795 mmHg (500 hPa ~ 1060 hPa)	577.5 mmHg ~ 795 mmHg (770 hPa ~ 1060 hPa)
Storage temperature	-20 °C ~ +70 °C	Room Temperature
Storage relative humidity	10 %RH ~ 95 %RH (non-condensing)	5 % RH ~ 85 % RH (non-condensing)
Storage ambient pressure	375 mmHg ~ 795 mmHg (500 hPa ~ 1060 hPa)	577.5 mmHg ~ 795 mmHg (770 hPa ~ 1060 hPa)

### Dimensions

Devices	Dimensions
CPM System	337mm*245mm*108 mm (Length*Width*Height)

### Weight

Devices	Weight
CPM System	4.8 Kg

### Pressure Measurement

<b>Pressure Range</b>	-30 mmHg ~ +300 mmHg
<b>Pressure Accuracy</b>	Pressure Microcatheter Pressure Accuracy: $\pm 1$ mmHg+( $\pm 1\%$ of reading) (-30 mmHg~50 mmHg); $\pm 3\%$ of reading (50 mmHg~300 mmHg);  AO Transducer Pressure Accuracy: $\pm 1$ mmHg+( $\pm 1\%$ of reading) (-30 mmHg~50 mmHg); $\pm 3\%$ of reading (50 mmHg~300 mmHg);  FFR/cRR: $\pm 1$ mmHg+( $\pm 1\%$ of reading) (-30 mmHg~50 mmHg); $\pm 3\%$ of reading (50 mmHg~300 mmHg)
<b>Zero Drift</b>	<5 mmHg/h
<b>Frequency response</b>	0-25Hz

### Interfaces Specifications

<b>PressureCath IN</b>	Proprietary Pressure Microcatheter Interface  Isolation: IEC 60601-1 Type CF Defibrillation protected
<b>AO IN</b>	Different Analog input.  Isolation: IEC 60601-1 Type CF Defibrillation protected  Excitation voltage: 5VDC  Pressure sensor sensitivity: 5uV/V/mmHg
<b>PressureCath OUT</b>	Analog output.  Isolation: IEC 60601-1 Type CF Defibrillation protected  Excitation voltage: 4-8VDC or AC 0-5kHz  Pressure sensor sensitivity: 5uV/V/mmHg
<b>AO OUT</b>	Analog output.  Isolation: IEC 60601-1 Type CF Defibrillation protected  Excitation voltage: 4-8VDC or AC 0-5kHz  Pressure sensor sensitivity: 5uV/V/mmHg
<b>VGA</b>	Connect with a standard VGA cable to VGA adapter
<b>USB</b>	Version: 2.0; Connector: Type A

<b>SD Card</b>	Supports the 8GB SD card
----------------	--------------------------

#### Accessories from Other Suppliers

Accessories	Name	Model
Aortic pressure (AO) transducer	SHENZHEN ANT HI-TECH INDUSTRIAL CO., LTD	PT-01

#### **Note:**

*Both CPM System and Pressure Microcatheter do not include the AO transducer. Only model PT-01 from SHENZHEN ANT HI-TECH INDUSTRIAL CO., LTD was verified to be used with the CPM System.*

#### Applied Parts

Aortic pressure (AO) transducer
Pressure Microcatheter

#### **Lifetime**

The lifetime of CPM System is 5 years, and Pressure Microcatheter should be connected with CPM System during use. Pressure Microcatheter is designed as a sterile and single use device, and its shelf life is 2 years.

## 12. EMC Information

#### **Guidance and manufacturer's declaration - electromagnetic emissions**

CPM System is intended for use in the electromagnetic environment specified below. The customer or user of CPM System should assure that it is used in such an environment.

Emission Test	Compliance	Electromagnetic Environment
Radiated Emissions CISPR 11	Group 1	CPM System uses RF energy only for its internal functioning. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Conducted Emission CISPR 11	Class A	
Harmonic emissions IEC 61000-3-2	Class A	CPM System is suitable for use in all establishments other than domestic and those directly connected with the public low-voltage
Voltage	Complies	

fluctuations/flicker emissions IEC 61000-3-3		power supply network that supplies buildings for domestic purposes.
---	--	---

#### **Guidance and manufacturer's declaration - electromagnetic immunity**

CPM System is intended for use in the electromagnetic environment specified below. The customer or user of CPM System should assure that it is used in such an environment.

<b>Immunity Test</b>	<b>IEC 60601 Test Level</b>	<b>Compliance Level</b>	<b>Electromagnetic Environment - Guidance</b>
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/ output lines	±2 kV for power supply lines ±1 kV for input/ output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV for line to line ±2 kV for line to ground	±1 kV for line to line ±2 kV for line to ground	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips and Interruptions IEC 61000-4-11	0% for 0.5 cycle, 0% for 1 cycle; 50Hz: 70% for 25 cycle, 0% for 250 cycle; 60Hz:70% for 25 cycle, 0% for 250 cycle;	0% for 0.5 cycle, 0% for 1 cycle; 50Hz: 70% for 25 cycle, 0% for 250 cycle; 60Hz:70% for 25 cycle, 0% for 250 cycle;	Mains power quality should be that of a typical commercial or hospital environment. If the user of CPM System requires continued operation during power mains interruptions, it is recommended that CPM System be powered from an uninterruptible power supply or battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Portable and mobile RF communication equipment should be used no closer to any part of CPM System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.			
Conducted RF IEC 61000-4-6	3V 6V (in ISM bands)	3V 6V (in ISM	$d = 1,17 \sqrt{P}$

	80% Am at 1kHz	bands) 80% Am at 1kHz	
Radiated RF IEC 61000-4-3	3V/m, 80% Am at 1kHz 27V/m PM at 18Hz 28V/m FM $\pm$ 5 kHz deviation at 1kHz sine 9V/m PM at 217 Hz 28V/m PM at 18Hz 28V/m PM at 217 Hz 28V/m PM at 217 Hz 9V/m PM at 217 Hz	3V/m, 80% Am at 1kHz 27V/m PM at 18Hz 28V/m FM $\pm$ 5 kHz deviation at 1kHz sine 9V/m PM at 217 Hz 28V/m PM at 18Hz 28V/m PM at 217 Hz 9V/m PM at 217 Hz	$d = 1,17 \sqrt{P}$ 80 MHz to 800 MHz $d = 2,33 \sqrt{P}$ 800 MHz to 2,5 GHz Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey a), should be less than the compliance level in each frequency range b). Interference may occur in the vicinity of equipment marked with the following symbol: 

**Note:** At 80 MHz and 800MHz, the higher frequency range applies.

**Note:** These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

**Note:** UT is the a.c. mains voltage prior to application of the test level.

- Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which equipment is used exceeds the applicable RF compliance level above, equipment should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorientation or relocating equipment.
- Over the frequency range 150 KHz to 80 MHz, field strengths should be less than 3 V/m.

### Recommended separation distances between portable and mobile RF communications equipment and CPM System

CPM System is intended for use in the electromagnetic environment where radiated RF disturbances are controlled. The customer or the user of CPM System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and CPM System as recommended below, according to the maximum

output power of the communications equipment.


Rated maximum output power of transmitter (W)	Separation distance to frequency of transmitter (m)		
	150 kHz to 80 MHz $d = 1,17\sqrt{P}$	80 MHz to 800 MHz $d = 1,17\sqrt{P}$	800 MHz to 2,5 GHz $d = 233\sqrt{P}$
0.01	0.12	0.12	0.23
0.10	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.33
<p>For transmitters rated at a maximum output power not listed above, the recommended separation distance (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.</p> <p><b>Note:</b> At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.</p> <p><b>Note:</b> These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.</p>			

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
**Gambar 1. Cardiovascular Pressure Measurement (CPM) System**

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**6.6 Fungsi Single Selection, Multi Selection, Select All**

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## 7. Pengaturan

### 7.1 Pengaturan Pengguna

The screenshot shows the 'Pengaturan Pengguna' (User Settings) window. At the top, there is a status bar with 'test', '001', 'Room', and '09-01-2022 09:55:29'. The main area contains several settings:

- Date: 09-01-2022
- Time: 09 : 55
- Language: English
- Screen Light: 100%
- Speed(mm/s): 25
- Scale(mmHg): 0 ~ 150
- Average(beats): 3
- Max Length(min): 60
- Test Pa: Off
- Test Pd: Off
- Display: ☒ Grid, ☒ Pa-Pd Trace
- Available/Capacity: 1.52GB/1.53GB
- Software Release Version: V2

At the bottom, there are buttons for 'Network Settings' and 'Restore Factory Settings'. The 'User' and 'Institution' fields are visible, along with a 'Save' button. The bottom navigation bar includes icons for 'Patient', 'Live', 'Review', and 'Archive', with the 'Live' button highlighted in blue.

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**MedNet EC-REP GmbH**
**MedNet SWISS GmbH**