## Thrombuster II


CE 0123

### [Contraindications]


- For single use only. Do not reuse, reprocess or resterilize.
- These components are contraindicated in the following lesion
  - left main trunk lesion unprotected by bypass graft or collateral circulation
  - The venous system
- These components are contraindicated in the following patients
  - patients who are responsive to surgical treatment
  - patients with lateral distal stenoses of the stent installed at the bifurcated lesion
  - patients with serious blood coagulation abnormality

### [Shape and structure]

#### 1. Aspiration Catheter

#### 2) Flushing needle


#### NOTE

- Flushing needle is used for priming guidewire lumen of Thrombuster II.

#### 2. Syringe Kit


##### 1) Syringe Kit

#### NOTE

- Please do not use other syringe and/or tubes for suction except this Syringe Kit

#### 2) Filter

#### NOTE

Filter is used for observing thrombus aspirated by Thrombuster II. Aspirated blood including thrombus can be filtrated by the filter and thrombus can be observed.

### [Indication for use]

#### 1. Indications

Thrombuster II is indicated for the myocardial infarction, angina pectoris and peripheral vascular thrombosis by percutaneous suction.

#### 2. Intended use

Thrombuster II is intended for use to remove thrombus and debris in the coronary or peripheral arteries by percutaneous suction.

### [Operation method and instructions for use]

#### (A) Operation in the presence of core wire

##### 1. preparations

- Aseptically remove this catheter from the package container and carrier tube.
- Prime the suction lumen and guidewire lumen of this catheter (using heparinized saline to replace the air)

#### Procedures for priming

- Draw 10 ml of heparinized saline into the suction syringe and expel any air inside the suction syringe.
- Attach the Syringe kit onto the connector of the core wire of this catheter, as shown in Figure 1.


Figure 1. Installation of Thrombuster II and Syringe kit (with the core wire)

- Air inside this catheter and syringe kit is replaced with the full volume of heparinized saline inside the suction syringe. In the case of incomplete air replacement, procedures a) – c) should be repeated, and after flushing, the stopcock should be closed.
- Attach the Flushing needle onto another prepared syringe containing heparinized saline and fill the guidewire lumen of the distal tip of this catheter with the same heparinized saline.

##### 2. Procedures for insertion into the body and suction

- Prior to use, check if this catheter is compatible to the guide catheter (not supplied in this set), the Y-connector (not supplied in this set) and the guidewire (not supplied in this set) on the following table.

catalogue No.	Compatible guide catheter	Min. guide catheter I.D.	Max. guidewire diameter
T2R7SA4W	7Fr	0.080" / 2.03mm	0.014" / 0.36mm
T2R6SA4W	6Fr	0.070" / 1.78mm	0.014" / 0.36mm

- 2) Insert the guide catheter, and after the installation of the Y-connector, advance the guidewire to the distal of the lesion. In case of usage of distal protection device, advance it to the end of the lesion.
- 3) Insert the proximal end of the guidewire into the distal part of the guidewire lumen of this catheter.
- 4) Through the Y-connector, insert this catheter under fluoroscopy to make the radiopaque marker at the distal tip of the catheter reach the target site.
- 5) Moderately lock the Y-connector to prevent blood loss and resistance with catheter during the procedure.
- 6) Disconnect the Syringe kit from the connector, pull out the core wire from this catheter.
- 7) Reattach the Syringe kit to the hub of this catheter.
- 8) After verifying the location of the catheter under fluoroscopy, pull the plunger of the suction syringe to the desired volume and reduce the pressure inside the suction syringe after locking the plunger.
- 9) Open the stopcock and start suction  
<in case of usage without protection balloon>  
Open the stopcock and slowly advance this catheter towards the thrombus for thrombus suction. Slowly advance the catheter through the whole area of the thrombus. Then draw back the catheter to the proximal side of the lesion.  
<in case of usage with protection balloon>  
In accordance with Instruction for Use of protection balloon, temporarily occlude the target lesion. Then open the stopcock and slowly advance the catheter through the whole area of the thrombus for suction.
- 10) Close the stopcock and remove the suction syringe. In case of usage with protection balloon, if urgent recanalization is necessary, deflate the protection balloon immediately according to its Instruction for Use.
- 11) If another suction is needed, after removing the sucked materials from the suction syringe, reattach the suction syringe onto the stopcock and repeat procedures 8) - 10).

### 3. Removal

- 1) Ensure the suction of the thrombus is completed
- 2) Ensure the stopcock is closed.
- 3) Slowly remove the catheter.
- 4) Ensure no thrombus remains inside the Y-connector.  
Any thrombus inside the Y-connector should be removed.
- 5) If reinsertion of this catheter is necessary, the insertion should be performed again from the priming procedure after the insertion of the core wire.

## (B) Operation in the absence of core wire

### 1. preparations

- 1) Aseptically remove the catheter from the package container and carrier tube.
- 2) Prime the suction lumen and guidewire lumen of this catheter (using heparinized saline to replace the air).

#### Procedures for priming

- a) Loose the handheld connector of the catheter and disconnect it from the hub and then remove the core wire from the catheter.
- b) Draw 10 ml of heparinized saline into the suction syringe and expel any air inside the suction syringe.
- c) Attach the Syringe kit to the hub of this catheter, as shown in Figure 2.


Figure 2: Installation of Thrombuster II and Syringe kit (without the core wire).

- d) Air inside this catheter and syringe kit is replaced with the full volume of heparinized saline inside the suction syringe. In the case of incomplete air replacement, procedures b) – d) should be repeated, and after flushing, the stopcock should be closed.
- e) Attach the Flushing needle onto another prepared syringe containing heparinized saline and fill the guidewire lumen of the distal tip of this catheter with the same heparinized saline.

### 2. Procedures for insertion into body and suction

- 1) Prior to use, please check if this catheter is compatible to the guide catheter (not supplied in this kit), the Y-connector (not supplied in this kit) and the guidewire (not supplied in this kit) on the following table.

catalogue No.	Compatible guide catheter	Min. guide catheter I.D.	Max. guidewire diameter
T2R7SA4W	7Fr	0.080" / 2.03mm	0.014" / 0.36mm
T2R6SA4W	6Fr	0.070" / 1.78mm	0.014" / 0.36mm

- 2) Insert the guide catheter, and after the installation of the Y-connector, advance the guidewire to the distal of the lesion. In case of usage of distal protection device, advance it to the end of the lesion.
- 3) Insert the proximal end of the guidewire into the distal part of the guidewire lumen of this catheter.

- 4) Through the Y-connector, insert this catheter under fluoroscopy to make the radiopaque marker at the distal tip of the catheter reach the target site.
- 5) Moderately lock the Y-connector to prevent blood loss and resistance with catheter during the procedure.
- 6) After verifying the location of the catheter under fluoroscopy, pull the plunger of the suction syringe to the desired volume and reduce the pressure inside the suction syringe after locking the plunger.
- 7) Open the stopcock and start suction  
<in case of usage without protection balloon>  
Open the stopcock and slowly advance this catheter towards the thrombus for thrombus suction. Slowly advance the catheter through the whole area of the thrombus. Then draw back the catheter to the proximal side of the lesion.  
<in case of usage with protection balloon>  
In accordance with Instruction for Use of protection balloon, temporarily occlude the target lesion. Then open the stopcock and slowly advance the catheter through the whole area of the thrombus for suction.
- 8) Close the stopcock and remove the suction syringe. In case of usage with protection balloon, if urgent recanalization is necessary, deflate the protection balloon immediately according to its Instruction for Use.
- 9) If another suction is needed, after removing the sucked materials from the suction syringe, reattach the suction syringe onto the stopcock and repeat procedures 6) - 8).

### 3. Removal

- 1) Ensure the suction of the thrombus is completed.
- 2) Ensure the stopcock is closed.
- 3) Slowly remove the catheter.
- 4) Ensure the thrombus remains inside the Y-connector.  
Any thrombus inside the Y-connector should be removed.
- 5) If reinsertion of this catheter is necessary, insertion should be performed again.

### [Precautions related to procedures]

1. The insertion and withdrawal of the catheter should be performed slowly with caution. (If this procedure is performed too fast, it may lead to kinking of the catheter shaft, damage to the guidewire lumen, or dispersing of the thrombus.)
2. Wide-angle separation between the catheter shaft and guidewire should be avoided. (It may lead to kinking of the catheter shaft or damage to the guidewire lumen.)
3. Prior to thrombus suction, the core wire should be removed from the catheter. (Thrombus suction cannot be performed with the core wire set.)
4. If the core wire cannot be withdrawn from the catheter, it should not be forced out. Pull the catheter out of the body to verify the cause and then perform the appropriate procedures. (The catheter may be damaged in the case of forced withdrawal)
5. When sucking thrombus without protection balloon, the catheter should be slowly advanced while gradually sucking from the proximal end of the thrombus. Precautions should be taken during operation, especially for thrombi at bifurcated lesion and the distal ends of bifurcated lesion. (If it is done too quickly, the thrombus may break up and it may lead to dispersion of thrombi at distal ends or side branches)
6. When sucking thrombus without protection balloon, if blood flow into the suction syringe has been stopped, the catheter should be slowly pulled back till the restoration of blood flow. Once blood flow inside the catheter has been restarted, the distal end of the catheter should be advanced to the distal end of the thrombus. Without the restart of blood flow, pressure inside the catheter should be kept at a reduced level since the distal end of the catheter may be blocked by thrombus. After being withdrawn, the inside of the catheter must be rinsed. (If the inner pressure of the catheter becomes positive, the thrombus may be scattered to the distal tip.)
7. If the thrombus remains stuck at the distal end of catheter, the withdrawal of the catheter inside the guide catheter must be performed carefully and slowly to prevent any falling off of the attached thrombus. (The thrombus falling off at the distal end of the guide catheter may lead to a new vascular occlusion).
8. Moderately lock the Y-connector to prevent blood loss and resistance with catheter during the procedure.
9. During the withdrawal or pulling back of this catheter, the guidewire should be fixed at the handheld side. (Drawing or pulling the catheter back without fixing the guidewire may lead to the guidewire being pulled out together with it. Furthermore, advancing the catheter without fixing the guidewire may lead to release of the guidewire from the catheter.)
10. Unless the distal end of the catheter has been blocked, the withdrawal of the catheter must be performed with the stopcock closed. (Otherwise, thrombus may be dispersed inside the body).
11. Under the circumstance of the thrombus remaining inside the Y-connector, infusion of contrast media or reinsertion of the catheter is forbidden. (It may reintroduce the thrombus into the blood vessels).
12. When reducing pressure by using the suction syringe, precaution should be taken to prevent hands from being clamped between the plunger and tube.
13. This catheter should only be used after being moistened by heparinized saline. (Without being totally moistened, it may lead to residual contrast media or blood coagulation inside the guide catheter, and therefore lubricity may be damaged to make the sliding unsmooth sometimes.)
14. For concomitant medical devices, the Instruction for Use should be followed.

## [Precautions during usage]

### [Important basic precautions]

1. This catheter may be used only by physicians skilled in percutaneous transluminal angioplasty.
2. For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or crossinfection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
3. Do not use if packaging is damaged.
4. Do not use agents containing organic solvents or oleaginous contrast media. Contact with these agents may lead to damage of Thrombuster II.
5. Since serious complications might arise when using this catheter, operation should be done in the medical institution where emergency procedure can be executed.
6. When aspirating thrombus, adjust the aspiration amount considering the diameter of the vessel from which the thrombus is sucked and the flow of blood into the target lesion. [If the diameter of the vessel from which the thrombus is sucked and the flow of blood into the target site is too small compared to the aspiration amount of this catheter, negative pressure may arise in the distal tip of the catheter.]
7. If the guide catheter is wedged in the entrance of the coronary artery or there is stenosis in the upstream of the target lesion for suction, perform the procedure cautiously as the vessel or the implanted stent may be deformed by the negative pressure arising from the aspirating manipulation. [If the flow of blood into the target lesion from which the thrombus is sucked is too small compared to the aspiration amount of this catheter, negative pressure may arise in the distal tip of the catheter.]
8. Especially in the case of usage for carotid artery, protection device usage is mandatory to prevent distal embolism. Temporarily occlude distal side of the target lesion by such protection devices.
9. During insertion of this catheter into vessels, the guide catheter should be infused with heparinized saline for anticoagulation.
10. This catheter must be primed before use to remove any residual air inside the catheter and syringe. [Any residual air inside the catheter and syringe may lead to incomplete suction.]
11. This catheter can only be inserted with the use of a guidewire. [Insertion of this catheter alone may lead to damage to the vascular wall or perforation of vessels.]
12. The guidewire must be completely advanced to reach the end of the lesion. [If the guidewire is not completely advanced to the end, it may loosen from the guidewire lumen and it may even lead to damage to the vascular wall or perforation of vessels.]
13. If the guidewire loosens from the guidewire lumen during operation, remove the catheter and then reinsert the guidewire.
14. The locations of the catheter and thrombus must be verified under fluoroscopy during operation. [During the insertion of this catheter into the area of the thrombus, this catheter may move the thrombus and cause it to disperse. In addition, during thrombus suction, attachment of the distal end of the catheter onto the inner wall of the vessel may lead to damage to the inner wall of the vessel.]
15. If abnormal or strong resistance is experienced during the operation, the cause for such abnormality or resistance should be verified and appropriate measures should be performed before proceeding. [If such abnormality or resistance is ignored and excessive force is applied, it may lead to damage of the vessels or to the catheter shaft breaking and remaining inside the body].
16. In the case that this catheter is advanced to the area where the stent is placed, insertion or pulling out should be performed slowly to avoid the distal end of the catheter becoming stuck on the stent strut. The gap between this catheter and the guidewire may catch the stent strut, which could damage the stent and this catheter. If any resistance is encountered, do not insert or remove the catheter forcefully. Check how the catheter is caught and try to remove it from the strut by adjusting the position of the catheter and the guidewire. If it is still not possible to detach the catheter, first, remove the guidewire, and then, move the catheter back and forth to take it away from the strut.
17. Do not advance this catheter to the distal part of the site where venous filter is implanted. [The strut of the venous filter may damage the catheter or break the catheter shaft leaving debris inside the body.]
18. During usage, the catheter shaft should be replaced for any bend, break or kink. [If the catheter continues to be used and such occurrence is ignored, suction may fail or the catheter shaft may be damaged and remain inside body].
19. If a great resistance is encountered during insertion, movement, or pulling out of this catheter, it should be verified that the guidewire is not tangled. If so, the tangling of the guidewire should be removed. [Since the guidewire lumen of this catheter is short, the guidewire may wind around the catheter shaft. In addition, while drawing this catheter back into the guide catheter inside vessels, the wide-angle separation between the catheter shaft and guidewire may occur. Under this circumstance, a forced withdrawal may lead to damage to the guidewire or catheter.]

20. Due to the lack of conductivity of twist forces, the catheter shaft should not be twisted [If twisted, the catheter shaft may be damaged and then remain inside the body].
21. The volume of thrombus suction should be decided by physicians for each lesion and adjusted via the suction syringe.
22. Air and sucked blood or thrombus in the suction syringe should never be injected into the body via the catheter.
23. This catheter should not be unreasonably inserted into or swiftly pulled out from lesions with highly tortuous vessels, bifurcated lesions, or calcified lesions. [The shaft at the distal end may then kink or be damaged, leading to vascular damage.]
24. Precaution should be taken to prevent any damage to the catheter by the surgical knife or scissors. Do not use a damaged catheter.
25. A catheter with any sign indicating damage should not be used.
26. During the usage of this catheter, the temperature, blood pressure, pulse, and respiration of patients should be monitored. In case of any abnormality, the procedure should be stopped or appropriate measures taken based on the physician's judgment.
27. After use, dispose of product and packaging in accordance with hospital, administrative and/or relevant national regulations.

#### [Adverse events]

Adverse events related to percutaneous removal of thrombus include, but are not limited to, infarction caused by occlusion of distal vessels or side branch, vasospasm, stripping of vascular endothelium, dissection of vascular intima, re-occlusion, vascular perforation, blood pressure fluctuation, stroke, shock, reaction to drugs, reaction to contrast media, renal insufficiency, transient ischemia, air embolism, thromboembolism, internal bleeding, hematoma, infection, etc. These adverse events may cause emergent coronary bypass surgery, myocardial infarction, re-stenosis, cardiac tamponade, hemorrhage, emergent brain surgery cerebral infarction, formation of vessel fistula, aneurysm, arrhythmia, and even death.

#### [Additional precautions]

1. Do not use if the product or package is believed to be damaged.
2. Use immediately after the sterile package is opened.  
Dispose as medical waste after use.

#### [Storage, care and expiration date]

1. Store in a cool, dark, and dry place avoiding exposure to water and direct sunlight, extreme temperature, or high humidity.
2. The expiration date is indicated on the box. Do not use after the expiration date.

#### [Package]

1 set (two packagings) / box

- 1) Aspiration Catheter packaging contains, an aspiration catheter and a flushing needle, and
- 2) Syringe Kit packaging contains, a Syringe Kit and a filter

#### [Name and address of manufacturer or representative]

##### Manufacturer:

##### KANEKA CORPORATION

Address: 3-18, 2-Chome, Nakanoshima, Kita-ku,  
Osaka-city, OSAKA, 530-8288 JAPAN

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FAX No.: (+32)-(0) 14-256-298

#### [Description of symbols for use]

Symbol	Description
	Do not re-use
	Do not resterilize
	Sterilized using ethylene oxide
	Caution, consult accompanying documents
	Consult instructions for use
	Do not use if package is damaged
	Keep away from sunlight
	Keep away from rain

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